

## Medaval Accreditation Assessment

Volume 2016

Report 1614

05 August 2016

### **Accreditation assessment of the blood pressure measurement technology used in the Avita BPM63S upper arm monitor, as validated according to the European Society of Hypertension International Protocol revision 2010**

#### **Approved by the Medaval Advisory Board**

Eoin O'Brien (Chair)	Yutaka Imai	Andrew Shennan
George S. Stergiou (Deputy Chair)	Martin Myers	Jan Staessen
Roland Asmar	Gbenga Ogedegbe	Martin J. Turner
Alejandro de la Sierra	Takayoshi Ohkubo	Paolo Verdecchia
Peter W. de Leeuw	Paolo Palatini	Bernard Waeber
Eamon Dolan	Gianfranco Parati	J-Guang Wang
Geoffrey A. Head		

#### **Reference**

Medaval Ltd. Accreditation assessment of the blood pressure measurement technology used in the Avita BPM63S upper arm monitor, as validated according to the European Society of Hypertension International Protocol revision 2010. *Medical Device Assessment*. 2016 Aug 5;2016(1614). 5 p. Epub: 2019 Jan 31. Available from: <https://www.medaval.ie/MDA/2016/MDA1614.pdf>.

*Medical Device Assessment* is published by

Medaval Ltd., Unit 107, SBC, Serpentine Ave., Ballsbridge, Dublin D04 H522, IRELAND.

© 2016-2019 Medaval Ltd. All rights reserved.

**Permissions:** Requests for permissions to reproduce figures, tables, or portions of reports or articles originally published in *Medical Device Assessment* can be obtained by email request to [info@medaval.ie](mailto:info@medaval.ie).

# Accreditation assessment of the blood pressure measurement technology used in the Avita BPM63S upper arm monitor, as validated according to the European Society of Hypertension International Protocol revision 2010

*Medaval Accreditation-Assessment Report – 5<sup>th</sup> August 2016*

### Test Device Details

		Assessment
<b>Full Name</b>	Avita BPM63S	Requirement satisfactory
<b>Model</b>	BPM63S	Requirement satisfactory
<b>Measurement Site</b>	Upper Arm	Requirement satisfactory
<b>Client Use</b>	Suitable for self-measurement.	Requirement satisfactory
<b>Operation Method</b>	Oscillometry, automatic during deflation	Requirement satisfactory
<b>Measurement Occurrence</b>	Single Measurements Only	Requirement satisfactory
<b>Device Photograph</b>		Modification: No photograph in paper. Standard image shown in report.
<b>Manufacturer(s)</b>	AVITA Corporation, 9F, No.78, Sec.1, Kwang-Fu Road, Sanchong Dist., New Taipei City 24158, TAIWAN.	Requirement satisfactory
<b>Cuffs</b>	Medium 23 cm to 33 cm Large 33 cm to 43 cm	Cuffs Listed: Requirement satisfactory Arm Circumferences: Requirement satisfactory

### Study Details

<b>Original Publication</b>	Kang YY, Zeng WF, Liu M, Li Y, Wang JG. Validation of the AVITA BPM63S upper arm blood pressure monitor for home blood pressure monitoring according to the European Society of Hypertension International Protocol revision 2010. <i>Blood Press Monit.</i> 2014 Feb; <b>19</b> (1):46-9. doi: 10.1097/MBP.0000000000000014. PMID: 24322871.	
<b>Protocol</b>	The European Society of Hypertension International Protocol revision 2010 for the validation of blood pressure measuring devices in adults <sup>1</sup>	

		Assessment
<b>Adherence</b>	Followed Precisely	Requirement satisfactory
<b>Adjustments</b>	None	Requirement satisfactory
<b>Study Meas. Method</b>	Oscillometric	Requirement satisfactory
<b>Study Measurement Site</b>	Upper Arm	Requirement satisfactory
<b>Observers</b>		
<b>Supervisor + 2 Observers</b>	Yes	Requirement satisfactory
<b>Observer Training</b>	BHS training video	Requirement satisfactory
<b>Observer Familiarisation</b>	12 measurements	Requirement satisfactory
<b>Observers Blinded</b>	From device and each other	Requirement satisfactory
<b>Sample</b>		
<b>Population</b>	A general population	Requirement satisfactory
<b>Circumstances</b>	None	Requirement satisfactory
<b>HBP Subjects Selection</b>	Hospital patients	Requirement satisfactory
<b>NBP Subjects Selection</b>	Hospital staff	Requirement satisfactory

<b>Test Device Details and Study Details Assessment</b>	<b>Checks</b>	22
	<b>Permitted Modifications</b>	1
	<b>Violations</b>	0

**Procedure**

**Table 1: Screening and Recruitment Details**

<b>Screening and Recruitment</b>				<b>Assessment</b>	
Total Screened			<b>47</b>	Value within requirements	
Total Excluded			<b>14</b>	Value within requirements	
	Ranges Complete		<b>9</b>	Value within requirements	
	Range Adjustment		<b>0</b>	Value within requirements	
	Arrhythmias		<b>2</b>	Value within requirements	
	Device Failure		<b>0</b>	Value within requirements	
	Poor Quality Sounds		<b>1</b>	Value within requirements	
	Cuff Size Unavailable		<b>0</b>	Value within requirements	
	Observer Disagreement		<b>0</b>	Value within requirements	
	Distribution		<b>0</b>	Value within requirements	
	Other Reasons*		<b>2</b>	Value within requirements	
Total Recruited			<b>33</b>	Value within requirements	
*Explanation Summary					
	Personal reasons			No details required	
<b>Recruitment Ranges</b>					
SBP	Total		<b>33</b>	Value within requirements	
	Low	< 90 mmHg	<b>2</b>	Value within requirements	
		90 – 129 mmHg	<b>9</b>	Value within requirements	
		130 – 160 mmHg	<b>10</b>	Value within requirements	
	Medium	161 – 180 mmHg	<b>10</b>	Value within requirements	
		> 180 mmHg	<b>2</b>	Value within requirements	
	High		<b>12</b>	Value within requirements	
DBP	Total		<b>33</b>	Value within requirements	
	Low	< 40 mmHg	<b>0</b>	Value within requirements	
		40 – 79 mmHg	<b>11</b>	Value within requirements	
		80 – 100 mmHg	<b>11</b>	Value within requirements	
	Medium	101 – 130 mmHg	<b>11</b>	Value within requirements	
		> 130 mmHg	<b>0</b>	Value within requirements	
	High		<b>11</b>	Value within requirements	
Total Extremes			<b>4</b>	Value within requirements	
<b>On Treatment Ranges</b>					
SBP	Low	< 130 mmHg	<b>2</b>	Value within requirements	
	Medium	130 – 160 mmHg	<b>8</b>	Value within requirements	
	High	> 160 mmHg	<b>9</b>	Value within requirements	
DBP	Low	< 80 mmHg	<b>4</b>	Value within requirements	
	Medium	80 – 100 mmHg	<b>8</b>	Value within requirements	
	High	> 100 mmHg	<b>7</b>	Value within requirements	
<b>Table 1 Assessment</b>				<b>Checks</b>	<b>36</b>
				<b>Permitted Modifications</b>	<b>0</b>
				<b>Violations</b>	<b>0</b>

### Study Results

**Table 2: Subject Details**

			Assessment	
Sex	Male:Female	19:14	Value within requirements	Value within requirements
Age (years)	Range (Low:High)	26:74	Value within requirements	Value within requirements
	Mean (SD)	47.1 (13.9)	Value within requirements	Value within requirements
Arm Circumference (cm)	Range (Low:High)	22:34	Value within requirements	Value within requirements
	Mean (SD)	28.0 (2.6)	Value within requirements	Value within requirements
Cuff for Test Device (cm)	Medium (23 – 33)	33		
	Large (33 – 43)			
	Total	33	Value within requirements	
Recruitment SBP (mmHg)	Range (Low:High)	?:?	Modification: Missing value accepted by paper review	Modification: Missing value accepted by paper review
	Mean (SD)	142.1 (32.3)	Value within requirements	Value within requirements
Recruitment DBP (mmHg)	Range (Low:High)	?:?	Modification: Missing value accepted by paper review	Modification: Missing value accepted by paper review
	Mean (SD)	88.2 (19.6)	Value within requirements	Value within requirements
<b>Table 2 Assessment</b>			<b>Checks</b>	19
			<b>Permitted Modifications</b>	4
			<b>Violations</b>	0

**Table 3: Observer Measurements in each Recruitment Range**

			Assessment	
SBP	Overall Range mmHg (Low:High)	85:209	Value within requirements	Value within requirements
	Low (< 130 mmHg)	24 to 38	Modification: Generality accepted by paper review	
	Medium (130 – 160 mmHg)	24 to 38	Modification: Generality accepted by paper review	
	High (> 160 mmHg)	24 to 38	Modification: Generality accepted by paper review	
	Maximum Difference	≤ 14	Modification: Generality accepted by paper review	
DBP	Overall Range mmHg (Low:High)	44:123	Value within requirements	Value within requirements
	Low (< 80 mmHg)	24 to 38	Modification: Generality accepted by paper review	
	Medium (80 – 100 mmHg)	24 to 38	Modification: Generality accepted by paper review	
	High (> 100 mmHg)	24 to 38	Modification: Generality accepted by paper review	
	Maximum Difference	≤ 14	Modification: Generality accepted by paper review	
<b>Table 3 Assessment</b>			<b>Checks</b>	12
			<b>Permitted Modifications</b>	8
			<b>Violations</b>	0

**Table 4: Observer Differences**



			Assessment	
Observer 2 – Observer 1				
SBP (mmHg)	Range (Low:High)	-4:+4	Value within requirements	Value within requirements
	Mean (SD)	+0.1 (1.7)	Value within requirements	Value within requirements
DBP (mmHg)	Range (Low:High)	-4:+4	Value within requirements	Value within requirements
	Mean (SD)	-0.5 (1.9)	Value within requirements	Value within requirements
Repeated Measurements		3	Value within requirements	
<b>Table 4 Assessment</b>			<b>Checks</b>	9
			<b>Permitted Modifications</b>	0
			<b>Violations</b>	0

**Table 5: Validation Results**

Part 1	Pass Req.		Achieved		Assessment	
	Two of	All of	SBP	DBP		
≤ 5 mmHg	73	65	<b>68</b>	<b>75</b>	Value within lower passing criteria	Value within passing criteria
≤ 10 mmHg	87	81	<b>89</b>	<b>95</b>	Value within passing criteria	Value within passing criteria
≤ 15 mmHg	96	93	<b>96</b>	<b>97</b>	Value within passing criteria	Value within passing criteria
Grade 1			<b>Pass</b>	<b>Pass</b>	Value within lower passing criteria	Value within passing criteria
Mean mmHg			<b>+2.6</b>	<b>-0.1</b>	Value within requirements	Value within requirements
SD mmHg			<b>6.2</b>	<b>5.2</b>	Value within requirements	Value within requirements
<b>Part 2</b>		Pass Req.	Achieved			
			SBP	DBP		
2/3 ≤ 5 mmHg		≥ 24	<b>24</b>	<b>25</b>	Value within passing criteria	Value within passing criteria
0/3 ≤ 5 mmHg		≤ 3	<b>1</b>	<b>2</b>	Value within passing criteria	Value within passing criteria
Grade 2			<b>Pass</b>	<b>Pass</b>	Value within passing criteria	Value within passing criteria
Grade 3			<b>Pass</b>	<b>Pass</b>	Value within lower passing criteria	Value within passing criteria
<b>Part 3</b>	Result		Pass			Value within lower passing criteria

<b>Table 5 Assessment</b>	<b>Checks</b>	21
	<b>Permitted Modifications</b>	0
	<b>Violations</b>	0

**Plots**



SBP Plot Provided **Yes**  
 DBP Plot Provided **Yes**

Requirement satisfactory  
 Requirement satisfactory

<b>Plots Assessment</b>	<b>Checks</b>	2
	<b>Permitted Modifications</b>	0
	<b>Violations</b>	0

**Recommendations**

**Overall Summary**

Number of checks 121  
 Number of permitted modifications 13  
 Number of violations 0

**Assessment Summary**

The validation has been checked and is verified as having been conducted in accordance with the protocol requirements. Therefore, the results are considered to be valid, the null hypothesis, that the device is inaccurate in measuring blood

pressure, is rejected and the conclusion, that the device is accurate for self-measurement in adults, is correct.

**Certification Decision**

The Avita BPM63S, with the Medium 23 cm to 33 cm cuff, is certified by Medaval Ltd., for blood pressure measurement in adults, as it fulfilled the conditions required for a pass in a validation study carried out in accordance with the requirements of the International Protocol of the European Society of Hypertension 2010 Revision.

Date of Advisory Board Approval: 29th July 2016.

**Reference**

- O'Brien E, Atkins N, Stergiou G, Karpettas N, Parati G, Asmar R, Imai Y, Wang J, Mengden T, Shennan A; Working Group on Blood Pressure Monitoring of the European Society of Hypertension. European Society of Hypertension International Protocol revision 2010 for the validation of blood pressure measuring devices in adults. *Blood Press Monit.* 2010;**15**:23-38. doi: 10.1097/MBP.0b013e3283360e98. PMID: 20110786. Erratum in *Blood Press Monit.* 2010;**15**(3):171-2.