(D)	Blutdruckmessgerät	
	Gebrauchsanweisung	2-13
(GB)	Blood pressure monitor	
	Instruction for Use.....	14-24
(F)	Tensiomètre	
	Mode d'emploi	25-36
(E)	Tensiómetro parlante	
	Manual de instrucciones.....	37-48
(RUS)	Прибор для измерения кровяного давления на запястье	
	Инструкция по применению.....	49-62
	Electromagnetic Compatibility Information	63-67

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Dear Customer,

Thank you for choosing one of our products. Our name stands for high-quality, thoroughly tested products for the applications in the areas of heat, weight, blood pressure, body temperature, pulse, gentle therapy, massage and air.

Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

Best regards,
Your Beurer Team

1. Getting to know your instrument

The wrist blood pressure monitor is used for non-invasive measurement and monitoring of adults' arterial blood pressure.

You can use it to measure your blood pressure quickly and easily, storing the results and displaying the progression of readings together with the average.

A warning is issued for anyone suffering from cardiac arrhythmia.

The values determined are classified and graphically evaluated according to WHO guidelines.

The device also has an acoustic voice output.

Keep these instructions carefully for further use and also let other users have access to them.



2. Important information

Signs and symbols

The following symbols are used in these instructions for use, on the packaging and on the type plate for the device and accessories:

	Caution
	Note Note on important information
	Follow instructions for use
	Type BF applied part
	Direct current

	Disposal in accordance with EC Directive – WEEE (Waste Electrical and Electronic Equipment).
	Manufacturer
Storage 	Permissible storage temperature and humidity
Operating 	Permissible operating temperature and humidity
	Keep dry
SN	Serial number
	The CE labelling certifies that the product complies with the essential requirements of Directive 93/42/EEC on medical products.

Advice on use

- In order to ensure comparable values, always measure your blood pressure at the same time of day.

- Before every measurement, relax for about five minutes.
- If you want to perform several measurements on the same person, wait five minutes between each measurement.
- Do not take a measurement within 30 minutes after eating, drinking, smoking or exercising.
- Repeat the measurement if you are unsure of the measured value.
- The measurements taken by you are for your information only – they are not a substitute for a medical examination! Discuss the measurements with your doctor, and never base any medical decisions on them (e.g. medicines and their administration)!
- Do not use the blood pressure monitor on newborns or patients with preeclampsia. We recommend consulting a doctor before using the blood pressure monitor during pregnancy.
- In the case of restricted circulation on the arm as a result of chronic or acute vascular diseases (including vascular constriction), the accuracy of the wrist measurement is limited. In this case you should avoid using an upper arm blood pressure monitor.
- Cardiovascular diseases may lead to incorrect measurements or have a detrimental effect on measurement accuracy. The same also applies to very low blood

pressure, diabetes, circulatory disorders and arrhythmias as well as chills or shaking.

- The blood pressure monitor must not be used in connection with a high-frequency surgical unit.
- Only use the unit on people who have the specified wrist measurement for the device.
- Please note that when inflating, the functions of the limb in question may be impaired.
- During the blood pressure measurement, blood circulation must not be stopped for an unnecessarily long time. If the device malfunctions, remove the cuff from the arm.
- Do not allow sustained pressure in the cuff or frequent measurements. The resulting restriction of the blood flow may cause injury.
- Ensure that the cuff is not placed on an arm in which the arteries or veins are undergoing medical treatment, e.g. intravascular access or therapy, or an arteriovenous (AV) shunt.
- Do not use the cuff on people who have undergone a mastectomy.
- Do not place the cuff over wounds as this may cause further injury.
- The blood pressure monitor can only be operated with batteries.
- To conserve the batteries, the monitor switches off automatically if no buttons are pressed for one minute.

- The device is only intended for the purpose described in these instructions for use. The manufacturer is not liable for damage resulting from improper or careless use.

Storage and Care

- The blood pressure monitor is made up of precision electronic components. Accuracy of readings and the instrument's service life depend on careful handling.
 - You should protect the device from impact, moisture, dirt, major temperature fluctuations and direct exposure to the sun's rays.
 - Never drop the device.
 - Do not use near strong electromagnetic fields, i.e. keep it away from any radio systems and mobile phones.
 - Only ever use the cuffs provided with the monitor or original replacement cuffs. Otherwise erroneous results will be recorded.
- Do not press any buttons until the cuff is in position.
- If the instrument is not used for any length of time, we recommend removing the batteries.

Advice on batteries

- Batteries can be fatal if swallowed. You should therefore store the batteries and products where they are

inaccessible to small children. If a battery has been swallowed, call a doctor immediately.

- Batteries should not be charged or reactivated with any other means, nor should they be taken apart, thrown in the fire or short-circuited.
- Remove the batteries from the instrument if they are worn out or if you are not going to use the instrument for any length of time. This prevents any damage as a result of leakage. Always replace all the batteries at the same time.
- Never use different types of battery, battery brands or batteries with different capacities. You should preferably use alkaline batteries.

Repair and disposal

- Batteries do not belong in domestic refuse. Used batteries should be disposed of at the collection points provided.
- Never open the instrument. If these instructions are not heeded, the warranty will be null and void.
- Never attempt to repair the instrument or adjust it yourself. We can no longer guarantee perfect functioning if you do.
- Repairs may only be performed by Beurer Customer Service or authorized dealers. However, always check

the batteries and replace them if necessary prior to making any complaint.

- The appliance should be disposed of according to Regulation-WEEE (Waste Electrical and Electronic Equipment). In case of queries, please contact the municipal authorities responsible for waste disposal in your area.


3. Unit description

1. Battery cover
2. Display
3. Loudspeaker
4. START/STOP button 
5. Memory button **M**
6. Wrist cuff
7. Volume button
louder **+**/quieter **-**
8. WHO - LED

Icons in the display:

1. Systolic pressure
2. Diastolic pressure
3. Pulse rate
4. Pulse icon ♥ / Cardiac arrhythmia icon (♥)
5. Number of memory space
6. Time and date
7. Inflate/deflate arrow
8. 'Change battery' icon
9. WHO classification
10. Memory Zones

4. Prepare measurement


Inserting battery

- Remove the battery cover (left side of the unit).
- Insert two 1.5 V micro batteries (alkaline, type LR 03). Make absolutely sure that you insert the batteries with the correct polarity as marked. Do not use rechargeable batteries.
- Replace the battery cover carefully.

If the 'change battery' icon  is displayed, measurement is no longer possible and you must replace all the batteries.

Once batteries have been removed from the device, the time must be reset. Any measurement results that have been stored are however retained.

Note: You will find these markings on batteries containing harmful substances: Pb = battery containing lead, Cd = battery containing cadmium, Hg = battery containing mercury.

Setting date and time

It is vital to set the date and time. This is the only way that measurements can be saved and recalled with the correct date and time.

The time is shown in 24 hour format. The device also has five languages. The device is set to German on delivery. Proceed as follows to set the date, time and language:

- Press the memory key **M** for over 5 seconds.
- The year starts flashing. Use the +/- key to set the year and confirm with the memory key **M**.
- Set the year, month, day and minute and confirm in each case with the memory key **M**.
- After setting the date and time you can set the language.
- Use the + and - volume keys to select the following languages: L1 = German, L2 = English, L3 = French, L4 = Spanish, L5 = Russian, L0 = no language. Confirm your selection with the memory key **M**.

- You can then set the volume of the selected language:
Vol 3 = high, Vol 2 = medium, Vol 1 = low. Confirm your selection with the memory key **M**.

Note: Language off **L0** means that voice output is not selected and the measurement value will not be announced.


Select memory

Switch on the device with the **⏻** button. When the device is switched off only the clock in the display is shown. After the announcement you can select the desired memory space by pressing the **+** or **-** volume button.

Two memories with 60 memory spaces each are available to save the measurement results of two different people separately or to save morning and evening measurement values separately.

5. Measuring blood pressure

Positioning cuff

- Bare your left wrist, making sure that the circulation in the arm is not restricted by any clothes etc. that are too tight. Position the cuff on the inside of your wrist.
- Fasten the cuff with the Velcro fastening so that the upper edge of the monitor is positioned approx. 1 cm below the ball of your thumb.
- The cuff has to be fitted tightly around the wrist but should not constrict it.

Important: The instrument should only be operated with the original cuff.

Correct posture

- Rest for approx. 5 minutes before each measurement. Otherwise there may be divergences.
- You can perform the measurement while sitting or lying. To carry out a blood pressure measurement, make sure you are sitting comfortably with your arms and back leaning on something. Do not cross your legs. Place your feet flat on the ground.

Make sure to rest your arm and move it. Always make sure that the cuff is at heart level. Otherwise significant deviations can occur. Relax your arm and the palm of your hand.

- In order not to distort the result, it is important to keep still during the measurement and not talk.

Measuring blood pressure

- Put on the cuff as described previously and assume the position in which you want to carry out the measurement.
- Start the measurement process by pressing START/STOP button . After checking the display with all digits lighting up, the monitor will automatically inflate. Following the self-test, during which all display elements briefly appear, the measurement begins. During inflation the device already determines measured values used to estimate the required inflation pressure. If this pressure is insufficient, the device automatically inflates by another 40 mmHg (Real Fuzzy Logic).
- Then the pressure in the cuff is slowly released and the pulse is measured.
- When the measurement has been completed, the remaining air pressure is released very quickly. The pulse, the systolic and the diastolic blood pressure are displayed.
- You can adjust the volume during the announcement by pressing + and -.
- You can interrupt measurement at any time by pressing the START/STOP button .
- E_r appears if it has not been possible to perform the measurement properly. Observe the section in these

instructions on error messages/troubleshooting and repeat the measurement.

- The test result is saved automatically.

Evaluating results

Cardiac arrhythmia:

This instrument can identify possible cardiac arrhythmia disorders during measurement and if necessary indicates the measurement with the flashing icon . This may be an indicator for arrhythmia. Arrhythmia is a condition where the heart rhythm is abnormal as a result of defects in the bioelectrical system controlling the heart beat. The symptoms (omitted or premature heart beats, slow or excessively fast heart rate) may be caused, among other things, by heart disease, age, physical predisposition, excessive use of stimulants, stress or lack of sleep. Arrhythmia can only be ascertained through examination by your doctor. Repeat the measurement if the flashing icon  is displayed after the measurement. Please note that you should rest for 5 minutes between measurements and not talk or move during the measurement. If the icon , appears often, please contact your doctor. Any self-diagnosis and treatment based on the test results may be dangerous. It is vital to follow your doctor's instructions.

WHO classification:

In accordance with the guidelines/definitions of the World Health Organization and the latest findings, the measurements can be classified and assessed according to the following table.

However, these standard values serve only as a general guideline, as the individual blood pressure varies in different people and different age groups etc.

It is important to consult your doctor regularly for advice. Your doctor will tell you your individual values for normal blood pressure as well as the value above which your blood pressure is classified as dangerous.

Range of blood pressure values	Systolic (in mmHg)	Diastolic (in mmHg)	Measure
Grade 3: Severe hypertension (Red)	≥180	≥110	Seek medical advice
Grade 2: Moderate hypertension (Orange)	160-179	100-109	Seek medical advice
Grade 1: Mild hypertension (Yellow)	140-159	90-99	Have it checked regularly by doctor

High-normal (Green)	130-139	85-89	Have it checked regularly by doctor
Normal (Green)	120-129	80-84	Check it yourself
Optimal (Green)	<120	<80	Check it yourself

Source: WHO, 1999

The bar graph in the display and the scale on the unit indicate the range of the blood pressure which has been recorded. An acoustic evaluation of your measurement also follows.

If the values for systolic and diastolic pressure are in two different WHO ranges (e.g. systolic in the high-normal range and diastolic pressure in the normal range) the graphic WHO classification on the unit indicates the higher range (high-normal in the example described).

6. Saving, retrieving and deleting results

- The results of each successful measurement are stored together with date and time. With more than 60 items of measured data, the earliest items of data measured are lost.

- In order to retrieve measurement results again, press the memory key **M**. Initially, the average value of the last three blood pressure readings is displayed (Display **R**). Press the memory key **M** again to display the latest measurements with date, time and graphic WHO classification and to hear the values announced. You can adjust the volume during the announcement by pressing the volume key **+** and **-**.
- You can clear the memory by pressing and holding **+** and **-** simultaneously for 5 seconds after pressing the memory key **M** and selecting the user.
- If you wish to change the user memory, see “Select memory”.

7. Error messages/troubleshooting

In case of faults, the E_r _ message appears in the display. Error messages may appear if:

- systolic or diastolic pressure could not be measured ($E_r 1$ or $E_r 2$ appears on the display)
- systolic or diastolic pressure was outside the measurement range (H_i or L_a appears on the display)
- the cuff is fastened too tightly or loosely ($E_r 3$ or $E_r 4$ appears on the display)
- the pump pressure is higher than 300 mmHg ($E_r 5$ appears on the display)
- pumping up takes longer than 160 seconds ($E_r 6$ appears on the display)

- there is a system or device error ($E_r R$, $E_r U$, $E_r 7$ or $E_r 8$ appears on the display)
- the batteries are almost empty .

In the above cases, you must repeat the measurement. Make sure that the cuff tube is properly inserted and that you do not move or talk.

Re-insert the batteries if necessary, or else replace them.

Technical alarm – description

Should the recorded blood pressure (systolic or diastolic) lie outside the limits specified in the section “Technical specifications”, the technical alarm will appear on the display indicating either “ H_i ” or “ L_a ”. In such cases, you should seek medical assistance and check the accuracy of your procedure.

The limit values for the technical alarm are factory set and cannot be adjusted or deactivated. These alarm limit values are accorded second priority under the standard IEC 60601-1-8.

The technical alarm is a non-locking alarm and must not be reset. The signal shown on the display will disappear automatically after about 8 seconds.

8. Cleaning and storing the instrument

- Clean your blood pressure monitor carefully using a slightly damp cloth only.
- Do not use any detergents or solvents.
- Never hold the instrument under water as otherwise liquid can penetrate and damage the instrument.
- Never place any heavy objects on the instrument.

9. Technical details

Model no.	BC 19
Measurement method	Oscillometric, non-invasive blood pressure measurement on the wrist
Measurement range	Cuff pressure 0–300 mmHg, systolic 60–260 mmHg, diastolic 40–199 mmHg, Pulse 40–180 beats/minute
Display accuracy	Systolic ± 3 mmHg, diastolic ± 3 mmHg, pulse $\pm 5\%$ of the value shown
Measurement inaccuracy	Max. permissible standard deviation according to clinical testing: systolic 8 mmHg/diastolic 8 mmHg
Memory	2 x 60 memory spaces

Dimensions	L 85 mm x W 70 mm x H 31 mm
Weight	Approx. 140 g (without batteries)
Cuff size	14.0 to 19.5 mm
Permissible operating conditions	+10 °C to +40 °C, $\leq 85\%$ relative air humidity (non-condensing)
Permissible storage conditions	-20 °C to +55 °C, $\leq 85\%$ relative air humidity, 800-1050 hPa ambient pressure
Power supply	2 x 1.5V  AAA batteries
Battery life	For approx. 250 measurements, depending on the blood pressure level and/or pump pressure
Accessories	Instruction for use, 2 x 1.5V AAA batteries, storage box
Classification	Internal supply, continuous operation, type BF applied part, IPX0, no AP or APG

Technical information is subject to change without notification to allow for updates.

- This unit is in line with European Standard EN 60601-1-2 and is subject to particular precautions with regard to electromagnetic compatibility (EMC). Please note that portable and mobile HF communication sys-

tems may interfere with this unit. More details can be requested from the stated Customer Service address or found at the end of the instructions for use.

- This device is in line with the EU Medical Devices Directive 93/42/EC, the “Medizinproduktegesetz” (German Medical Devices Act) and the standards EN 1060-1 (non-invasive sphygmomanometers, Part 1: General requirements), EN 1060-3 (non-invasive sphygmomanometers, Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems) and IEC 80601-2-30 (Medical electrical equipment – Part 2–30: Particular requirements for the safety and essential performance of automated non-invasive blood pressure monitors).
- The accuracy of this blood pressure monitor has been carefully checked and developed with regard to a long useful life. If using the device for commercial medical purposes, it must be regularly tested for accuracy by appropriate means. Precise instructions for checking accuracy may be requested from the service address.

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Table 1

For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacturer's declaration – electromagnetic emissions

The BC 19 is intended for use in the electromagnetic environment specified below. The customer or the user of the BC 19 should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The BC 19 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The BC 19 is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable	

Table 2**For all ME EQUIPMENT and ME SYSTEMS***Guidance and manufacturer's declaration – electromagnetic immunity*

The BC 19 is intended for use in the electromagnetic environment specified below. The customer or the user of the BC 19 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: U_T is the a.c. mains voltage prior to application of the test level.

Table 3

For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING

Guidance and manufacturer's declaration – electromagnetic immunity

The BC 19 is intended for use in the electromagnetic environment specified below. The customer or the user of the BC 19 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the BC 19, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2.3 \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the BC 19 is used exceeds the applicable RF compliance level above, the BC 19 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the BC 19.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.

Table 4**For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING***Recommended separation distances between portable and mobile RF communications equipment and the BC 19*

The BC 19 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the BC 19 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the BC 19 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter (m)	
	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.23
0.1	0.38	0.73
1	1.2	2.3
10	3.8	7.3
100	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.