BM 70

(D) Gebrauchsanweisung	
Blutdruckmessgerät	2-13
(GB) Instructions for use	
Blood pressure monitor	14-24
(F) Mode d'emploi	
Tensiomètre	25-36
(E) Manual de instrucciones	
Tensiómetro	37-48
(I) Istruzioni per l'uso	
Misuratore di pressione	49-60
(TR) Kullanım kılavuzu	
Bilgisayarlı tansiyon ölçer	61-71
(RUS) Инструкция по применению	
Прибор для измерения артериального давления на предплечье	72-84
(PL) Instrukcja obsługi	
Ciśnieniomierz	85-95
Electromagnetic Compatibility Information	96-99

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Dear Customer,

thank you for choosing one of our products. Our name stands for high-quality, thoroughly tested products for the applications in the areas of heat, weight, blood pressure, body temperature, pulse, gentle therapy, massage and air.

Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

Best regards,
Your Beurer Team

1. Getting to know your instrument

The upper arm blood pressure monitor is used for non-invasive measurement and monitoring of adults' arterial blood pressure. You can use it to measure your blood pressure quickly and easily, storing the results and displaying the progression of readings together with the average.

A warning is issued for anyone suffering from cardiac arrhythmia. The values determined are classified and graphically evaluated according to WHO guidelines.

Additionally, this blood pressure monitor includes a haemodynamic stability indicator, which will subsequently be referred to as a resting indicator. This shows whether the circulatory system is sufficiently at rest when taking a blood pressure measurement and is therefore a more precise indicator of your resting blood pressure. Read more about this on page 21–22.

Keep these instructions carefully for further use and also let other users have access to them.

2. Important information

Signs and symbols

The following symbols are used in these instructions for use, on the packaging and on the type plate for the device and accessories:

	Caution
	Note Note on important information
	Follow instructions for use
	Type BF applied part
	Direct current
	Disposal in accordance with EC Directive 2002/96/EC – WEEE (Waste Electrical and Electronic Equipment).
	Manufacturer

	Temperature limitation
	Permissible storage humidity
	Keep dry
SN	Serial number
	The CE labelling certifies that the product complies with the essential requirements of Directive 93/42/EEC on medical products.

Advice on use

- In order to ensure comparable values, always measure your blood pressure at the same time of day.
- Before every measurement, relax for about five minutes.
- If you want to perform several measurements on the same person, wait five minutes between each measurement.
- Do not take a measurement within 30 minutes after eating, drinking, smoking or exercising.
- Repeat the measurement if you are unsure of the measured value.
- The measurements taken by you are for your information only – they are not a substitute for a medical examination! Discuss the measurements with your doctor, and never base any medical decisions on them (e.g. medicines and their administration)!

- Do not use the blood pressure monitor on newborns, pregnant women or patients with preeclampsia.
- Cardiovascular diseases may lead to incorrect measurements or have a detrimental effect on measurement accuracy. The same also applies to very low blood pressure, diabetes, circulatory disorders and arrhythmias as well as chills or shaking.
- The blood pressure monitor must not be used in connection with a high-frequency surgical unit.
- Only use the device on people who have the specified upper arm measurement for the device.
- Please note that when inflating, the functions of the limb in question may be impaired.
- During the blood pressure measurement, blood circulation must not be stopped for an unnecessarily long time. If the device malfunctions, remove the cuff from the arm.
- Avoid any mechanical restriction, compression or bending of the cuff line.
- Do not allow sustained pressure in the cuff or frequent measurements. The resulting restriction of the blood flow may cause injury.
- Ensure that the cuff is not placed on an arm in which the arteries or veins are undergoing medical treatment, e.g. intravascular access or therapy, or an arteriovenous (AV) shunt.
- Do not use the cuff on people who have undergone a mastectomy.
- Do not place the cuff over wounds as this may cause further injury.
- You can either use the blood pressure monitor with batteries or with a mains part. Please note that data transfer and data storage is only possible when your blood pressure monitor is

supplied with power. As soon as the batteries are empty or the mains part is disconnected from the power supply, the blood pressure monitor loses the date and time.

- To conserve the batteries, the monitor switches off automatically if no buttons are pressed for one minute.
- The device is only intended for the purpose described in these instructions for use. The manufacturer is not liable for damage resulting from improper or careless use.

Storage and Care

- The blood pressure monitor is made up of precision electronic components. Accuracy of readings and the instrument's service life depend on careful handling.
 - You should protect the device from impact, moisture, dirt, major temperature fluctuations and direct exposure to the sun's rays.
 - Never drop the device.
 - Do not use near strong electromagnetic fields, i.e. keep it away from any radio systems and mobile phones.
 - Only ever use the cuffs provided with the monitor or original replacement cuffs. Otherwise erroneous results will be recorded.
- Do not press any buttons until the cuff is in position.
- If the instrument is not used for any length of time, we recommend removing the batteries.

Advice on batteries

- Batteries can be fatal if swallowed. You should therefore store the batteries and products where they are inaccessible to small children. If a battery has been swallowed, call a doctor immediately.

- Batteries should not be charged or reactivated with any other means, nor should they be taken apart, thrown in the fire or short-circuited.
- Remove the batteries from the instrument if they are worn out or if you are not going to use the instrument for any length of time. This prevents any damage as a result of leakage. Always replace all the batteries at the same time.
- Never use different types of battery, battery brands or batteries with different capacities. You should preferably use alkaline batteries.


Repair and disposal


- Batteries do not belong in domestic refuse. Used batteries should be disposed of at the collection points provided.
- Never open the instrument. If these instructions are not heeded, the warranty will be null and void.
- Never attempt to repair the instrument or adjust it yourself. We can no longer guarantee perfect functioning if you do.
- Repairs may only be performed by Beurer Customer Service or authorized dealers. However, always check the batteries and replace them if necessary prior to making any complaint.
- The appliance should be disposed of according to Regulation 2002/96/EC-WEEE (Waste Electrical and Electronic Equipment). In case of queries, please contact the municipal authorities responsible for waste disposal in your area.

3. Unit description


1. Cuff tube
2. Cuff
3. Cuff connector
4. Display
5. Memory button **M**
6. ON/OFF button
7. User selection button
8. Date/time button
9. Adjustment button
10. Connector, power supply
11. Data interface
12. Cuff connector port (left side)

Icons in the display:


1. Date/Time indicator
2. WHO classification
3. Memory zones
4. Resting indicator
5. Average value (AVG)
6. Memory sequence number
7. Weak battery mark
8. Cardiac arrhythmia icon
9. Pulse mark
10. Pulse rate
11. Diastolic pressure
12. Systolic pressure

4. Prepare measurement

Inserting battery

- Remove the battery cover from the back of the monitor.
- Insert four AA 1.5V alkaline batteries, making absolutely sure that you insert them with the correct polarity as marked. Never use rechargeable batteries.
- Replace the battery cover carefully.


If the battery change is continuously illuminated, measurement is no longer possible and you must replace all the bat-

teries. Once batteries have been removed from the device, the time must be reset.

Used batteries do not belong in the household waste. You are legally obliged to dispose of the batteries. Dispose of them via your specialist electrical supplier or local collecting point for recyclable waste.

Note: Batteries containing pollutant substances are marked as follows: Pb = Battery contains lead, Cd = Battery contains cadmium, Hg = Battery contains mercury.

Setting date and time

It is vital to set date and time. Only in this way is it possible to save and subsequently retrieve your measured values with the right date and time.

Press the key and to adjust the month. Then press the and in turn to adjust the date, hour, minute and confirm the entry .

The time is being displayed in 12-hour format, which means a time from 13:00 onwards is being displayed as 01:00 PM.

Mains operation

You can also operate this device with a mains adapter. No batteries should be in the battery compartment for this. The mains adapter is available from retailers or from the service address under order no. 071.29. The blood pressure monitor should only be used with the mains adapters specified here. The mains adapter should only be connected to the mains voltage indicated on the rating plate. Once you have disconnected the device from the mains, the date and time on the blood pressure monitor are lost. Any measurement results that have been stored are however retained.

PC interface

The Beurer blood pressure monitor also allows you to transfer your measured values to the PC.

To do this you require a data cable and the Beurer “Health Manager” PC software.

A data cable can be obtained from the specified Customer Service address with order no. 162.044. The software can be downloaded free of charge from www.beurer.com/service/download.

System requirements for the Beurer “Health Manager” computer software

1. Supported operating systems:

- Windows XP SP3
- Windows Vista SP1 or later
- Windows 7
- Windows 7 SP1

2. Supported architectures:

- x86 (32 bit)
- x64 (64 bit)


3. Hardware requirements:

- Recommended: at least Pentium 1 GHz or faster with at least 1 GB RAM.
- Free memory on the primary partition of at least:
 - x86 – 600 MB
 - x64 – 1.5 GB
- Graphic resolution from: 1024 x 768 pixels.
- USB port 1.0 or later.


5. Measuring blood pressure

5.1 Positioning cuff


Fit the cuff round your bare left upper arm. Blood circulation in the arm should not be restricted by tight clothing or other objects.


The cuff should be placed on the upper arm so that the lower edge is 2 to 3 cm above the bend of the elbow and above the artery. The tube should be in line with the centre of the palm.

Now place the free end of the cuff snugly, but not too tightly, around the arm, and fix it with the Velcro fastener. The cuff should be fitted tight enough to allow just two fingers to fit beneath the cuff.



Insert the cuff tubing into the socket for the cuff attachment.

Important: The instrument should only be operated with the original cuff. The cuff is suitable for an arm circumference of 24 to 36 cm. A larger cuff for upper arm circumferences from 34 to 46 cm is available from retailers or the service address under order no. 162.797.

5.2 Correct posture

- Rest for approx. 5 minutes before each measurement. Otherwise there may be divergences.
- You can perform the measurement either sitting or lying down. Always make sure that the cuff is on a level with your heart.
- To carry out a blood pressure measurement, make sure you are sitting comfortably with your arms and back leaning on something. Do not cross your legs. Place your feet flat on the ground.
- In order not to distort the result, it is important to keep still during the measurement and not talk.

5.3 Select memory

You have two memories (60 memory spaces each) in order to save the test results of 2 different people separately, or else save measurements in the morning and evening separately. Select the required memory space by pressing the user selection button . Confirm your selection with the On/Off button .

5.4 Measuring blood pressure

- Put on the cuff as described previously and assume the position in which you want to carry out the measurement.
- Press the user selection button  to select a memory zone, and press the  key 2x to start measurement in the chosen memory zone. After checking the display with all digits lighting up, the monitor will automatically inflate. Following the

self-test, during which all display elements briefly appear, the measurement begins. During inflation the device already determines measured values used to estimate the required inflation pressure. If this pressure is insufficient, the device automatically inflates by another 40 mmHg (Real Fuzzy Logic).

- Then the pressure in the cuff is slowly released and the pulse is measured.
- It displays pulse, systolic and diastolic blood pressure as well as the stop indicator (see Chapter 5.6).
- You can interrupt measurement at any time by pressing the ON/OFF button .
- To switch off and release the pressure, press the ON/OFF button  again. If you forget to switch off the device, it switches off automatically after approx. 1 minute.

Wait at least 5 minutes before taking another measurement!

5.5 Evaluating results

Cardiac arrhythmia:

This instrument can identify possible cardiac arrhythmia disorders during measurement and if necessary indicates the measurement with the flashing icon .

This may be an indicator for arrhythmia. Arrhythmia is a condition where the heart rhythm is abnormal as a result of defects in the bioelectrical system controlling the heart beat. The symptoms (omitted or premature heart beats, slow or excessively fast heart rate) may be caused, among other things, by heart disease, age, physical predisposition, excessive use of stimulants, stress or lack of sleep. Arrhythmia can only be ascertained through examination by your doctor.

Repeat the measurement if the flashing icon  is displayed after the measurement. Please note that you should rest for

5 minutes between measurements and not talk or move during the measurement. If the icon  appears often, please contact your doctor. Any self-diagnosis and treatment based on the test results may be dangerous. It is vital to follow your doctor's instructions.

WHO classification:

The World Health Organization (WHO) and National High Blood Pressure Education Program Coordinating Committee have developed a blood pressure standard, according to which areas of low- and high-risk blood pressure are identified. This standard, however, is a general guideline as individual's blood pressure varies among different people and different age groups ... etc. It is important that you consult with your physician regularly. Your physician will tell you your normal blood pressure range as well as the point at which you will be considered at risk.

Table for classification of blood pressure values (unit of measurement mmHg) for adults:

Range	Systolic	Diastolic	Action
Hypotonia (degraded blood pressure)	lower than 105	lower than 60	Check with doctor
Normal range	between 105 and 120	between 60 and 80	Self check
Pre-hypertensive blood pressure range* (P)	between 120 and 140	between 80 and 90	Check with doctor
High blood pressure - degree 1 (1)	between 140 and 160	between 90 and 100	Consult your doctor
High blood pressure - degree 2 (2)	higher than 160	higher than 100	Consult your doctor

(*) Blood pressure range that might transfer to high blood pressure.

Adapted from JNC 2003

The WHO classification in the display shows the range of the measured blood pressure.

If the values for systolic and diastolic pressure are in two different WHO ranges (e.g. systolic in the grade 1 hypertension range and diastolic pressure in the normal range), the WHO classification on the unit always indicates the higher range (grade 1 hypertension in the example described).

5.6 Resting indicator measurement (using HSD diagnostics)

The most frequent error made when measuring blood pressure is taking the measurement when not at rest (haemodynamic stability), which means that both the systolic and the diastolic blood pressures are incorrect in this case. During blood pressure measurement, the device automatically determines whether the circulatory system is sufficiently at rest or not.

If there is no indication that the circulatory system is not sufficiently at rest,  (haemodynamic stability) is displayed and the measurement can be recorded as a reliable resting blood pressure value.

: Haemodynamic stability

Measurement of the systolic and diastolic pressure is increased when the circulatory system is sufficiently at rest and is a very reliable indicator of resting blood pressure. However, if the circulatory system is not sufficiently at rest (haemodynamic instability), the symbol  is displayed.

In this case, the measurement should be repeated after a period of physical and mental rest. The blood pressure measurement must be taken when the patient is physically and mentally rested, as it will be the basis for a diagnosis and regulation of the patient's medical treatment.

: Lack of haemodynamic stability

It is very probable that the measurement of systolic and diastolic blood pressure has not been carried out at rest and has therefore distorted the measurement.

Repeat the measurement after a rest and relaxation period of at least 5 minutes. Go to a sufficiently calm and comfortable spot and remain there calmly; close your eyes, breathe deeply and evenly and try to relax.

If the following measurement shows a lack of stability again, repeat the measurement after another period of rest. In case of further unstable measurements, identify these measurements as they were taken when the circulatory system had not been sufficiently rested.

In this case, inner unrest of the nervous system could be the cause, which cannot be cured by brief periods of rest. Developing cardiac arrhythmia can also prevent stable blood pressure measurements.

A lack of resting blood pressure can have various causes, such as physical or mental strain or distraction and speaking or experiencing heart rhythmic disturbance during the measurement.

In an overwhelming number of cases, the HSD diagnosis will give a very good guide as to whether the circulatory system is rested when taking the measurement. Certain patients suffering from cardiac arrhythmia or chronic mental conditions can remain haemodynamically unstable in the long-term, something which persists even after repeated periods of rest. The precision of the results of the resting blood pressure is reduced in these users. Like any medical measurement method, the precision of the HSD diagnosis is limited and can lead to incorrect results in some cases. The blood pressure measurements taken when the circulatory system was at rest represent particularly reliable results.

6. Saving, retrieving and deleting results

The computer automatically saves the blood pressure readings of the last 60 measurements. When the 60 memory locations are exceeded, the current oldest value is deleted.

You can display the memory as follows:

- Press the user selection button  to select a memory zone, and press the **M** key to switch from memory to memory with sequence numbers.
- The average value **AVG** of the last three measurements is displayed first.
- Continue pressing the memory button to recall other individual measurement results - the first result to appear is that of the most recent value measured.
- Press the On/Off switch  to turn the unit off.
- If you forget to switch off the computer, it will switch off automatically after 1 minute.

Press the user selection button  to select a memory zone, then press and hold the **M** key for approximately 5-seconds to clear the data in the pre-designated memory zone.

7. Cleaning and storing the instrument

- Clean your device and cuff carefully only with a slightly moistened cloth.
- Do not use detergents or solvents.
- On no account must you immerse the computer in water, otherwise liquid can enter it and cause damage.
- When storing the device, make sure that no heavy objects are placed on top of it. Remove the batteries. The cuff tube should not have any sharp kinks.

8. Rectifying faults

Error messages can occur when

- measurement error (EE appears in the display),
- you move or talk during the measurement (EE appears in the display),

- the cuff tube is not properly inserted (E l appears in the display),
- inflation takes longer than 15 seconds (E l appears in the display),
- the inflation pressure is higher than 300 mmHg (E 2 appears in the display),
- an error occurs when storing the measured values (E 3 appears in the display),
- exceeding measurement range (E r appears in the display).

In the about cases, you must repeat the measurement. Make sure that the cuff tube is properly inserted and that you do not move or talk. Put the batteries back in if necessary, or else replace them.

9. Specifications

Model no.	BM 70
Measurement method	Oscillometric, non-invasive blood pressure measurement on the upper arm
Measurement range	Cuff pressure 0–300 mmHg, systolic 30–260 mmHg, diastolic 30–260 mmHg, Pulse 40–199 beats/minute
Display accuracy	Systolic ± 3 mmHg, diastolic ± 3 mmHg, pulse ± 5 % of the value shown
Measurement inaccuracy	Max. permissible standard deviation according to clinical testing: systolic 8 mmHg/diastolic 8 mmHg

Memory	2 x 60 memory spaces
Dimensions	L 156 mm x W 117 mm x H 80 mm
Weight	Approx. 535 g (without batteries)
Cuff size	24 to 36 cm
Permissible operating conditions	+10 °C to +40 °C, 40–85 % relative air humidity (non-condensing)
Permissible storage conditions	-10 °C to +60 °C, 10–90 % relative air humidity, 800–1050 hPa ambient pressure
Power supply	4 x 1,5V  AA batteries
Battery life	For approx. 300 measurements, depending on the blood pressure level and/or pump pressure
Accessories	Instruction for use, 4 x 1.5V AA batteries, storage pouch
Classification	Internal supply, IPX0, no AP or APG, continuous operation, type BF applied part

These specification are subject to change without notice for purpose of improvement.

10. Adapter

Model No.	FW 7333SM/12
Input	100–240V, 50–60 Hz
Output	12V DC, 700 mA, only in connection with beurer blood pressure monitor
Supplier	Friwo Gerätebau GmbH
Protection	This device is double insulated and protected against short circuit and overload by a primary thermal fuse. The green LED lights under normal conditions. Make sure to take the batteries out of the compartment before using the adapter.
	Double insulated / equipment class 2
Enclosures and Protective Covers	Equipment enclosed to protect against contact with live parts, and with parts which can become live (finger, pin, hook test). The operator shall not contact the patient and the output plug of AC adaptor simultaneously.

- This unit is in line with European Standard EN 60601-1-2 and is subject to particular precautions with regard to electromagnetic compatibility (EMC). Please note that portable and mobile HF communication systems may interfere with this unit. More details can be requested from the stated Customer Service address or found at the end of the instructions for use.

- This device is in line with the EU Medical Devices Directive 93/42/EC, the “Medizinproduktegesetz” (German Medical Devices Act) and the standards EN 1060-1 (non-invasive sphygmomanometers, Part 1: General requirements), EN 1060-3 (non-invasive sphygmomanometers, Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems) and IEC 80601-2-30 (medical electrical equipment – Part 2–30: Particular requirements for the safety and essential performance of automated non-invasive blood pressure monitors).
- The accuracy of this blood pressure monitor has been carefully checked and developed with regard to a long useful life. If using the device for commercial medical purposes, it must be regularly tested for accuracy by appropriate means. Precise instructions for checking accuracy may be requested from the service address.

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Table 1

For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacture's declaration – electromagnetic emissions

The BM 70 is intended for use in the electromagnetic environment specified below.

The customer or the user of the BM 70 should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The BM 70 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The BM 70 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Table 2*Guidance and manufacturer's declaration – electromagnetic immunity*

The BM 70 is intended for use in the electromagnetic environment specified below. The customer or the user of the BM 70 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines N/A	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV differential mode N/A	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 s	<5 % U_T (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the BM 70 requires continued operation during power mains interruptions, it is recommended that the BM 70 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: U_T is the a.c. mains voltage prior to application of the test level.

Table 3
Guidance and manufacturer's declaration – electromagnetic immunity

The BM 70 is intended for use in the electromagnetic environment specified below. The customer or the user of the BM 70 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3V _{rms} 150 kHz to 80 MHz	3V _{rms}	Portable and mobile RF communications equipment should be used no closer to any part of the BM 70, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1,2 \sqrt{P}$ $d = 1,2 \sqrt{P}$ 80 MHz to 800 MHz
Radiated RF IEC 61000-4-3	3V/m 80 MHz to 2,5 GHz	3V/m	$d = 2,3 \sqrt{P}$ 800 MHz to 2,5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the BM 70 is used exceeds the applicable RF compliance level above, the BM 70 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the BM 70. ^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.			

Table 4**Recommended separation distances between portable and mobile RF communications equipment and the BM 70**

The BM 70 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the BM 70 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the BM 70 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.