
ForaCare[®] Test N'GO[™]

Blood Glucose Monitoring System

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Owner's Manual
Manual del Usuario

Dear Owner of ForaCareTest N'GO™ Blood Glucose Monitoring System:

The system consists of three main products: the meter, test strips and control solutions. These products have been designed, tested, and proven to work together as a system to produce accurate blood glucose test results. Only use **ForaCare Test N'GO™** test strips and control solutions with the **ForaCare Test N'GO™** Blood Glucose Monitoring System.

TABLE OF CONTENTS

IMPORTANT SAFETY INSTRUCTIONS	04
Read Before Use	04
Important Information	06
Intended Use	07
Test Principle	07
Appearance And Key Functions of The Meter	08
Meter Display	09
Test Strip	10
Setting the Meter And Deleting the Memory	11
CHECKING THE SYSTEM WITH CONTROL SOLUTIONS	14
How to Perform A Control Solution Test	14
TESTING YOUR BLOOD	17
The Four Measuring Modes	17
How to Perform A Blood Glucose Test	18
Preparing the Puncture Site	19
Performing A Blood Glucose Test	16
VIEWING THE METER MEMORY	23
Viewing Results	23
Viewing Day Average Results	24
DOWNLOADING RESULTS ONTO A COMPUTER	26
CLEANING AND DISINFECTING YOUR METER	27
Caring for Your Meter	27
Disinfecting Procedures	28
BATTERY RECHARGE	30
DISPLAY MESSAGES AND PROBLEM-SOLVING GUIDE	31
Troubleshooting	33
DETAILED INFORMATION	34
Reference Values	34
Comparing Meter and Laboratory Results	34
SPECIFICATIONS	36
Federal Communications Commision (FCC) Statement	37

IMPORTANT SAFETY INSTRUCTIONS

Read Before Use

- The meter and lancing device are for single patient use. Do not share them with anyone including other family members! Do not use on multiple patients!
- All parts of the kit are considered biohazardous and can potentially transmit infectious diseases, even after you have performed cleaning and disinfection.

For more information, please visit

1. “FDA Public Health Notification: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens: Initial Communication” (2010)

<http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/ucm224025.htm>

2. “CDC Clinical Reminder: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens” (2010)

<http://www.cdc.gov/injectionsafety/Fingerstick-DevicesBGM.html>

1. Use this device **ONLY** for the intended use described in this manual.
2. Do **NOT** use accessories which are not specified by the manufacturer.
3. Do **NOT** use the device if it is not working properly or or if it is damaged.
4. Do **NOT** use the equipment in places where aerosol sprays are being used, or where oxygen is being administered.
5. This device does **NOT** serve as a cure for any symptoms or diseases.
The data measured is for reference only.
6. Before using this device to test blood glucose, read all instructions thoroughly and practice the test. Carry out all the quality control checks as directed.
7. Keep the device and testing equipment away from young children.
Small items such as the battery cover, batteries, test strips, lancets and vial caps are choking hazards.
8. Use of this device in a dry environment, especially if synthetic materials are present (synthetic clothing, carpets etc.) may cause damaging static discharges that may cause erroneous results.
9. Do **NOT** use this device in close proximity to sources of strong electromagnetic radiation, as these may interfere with the accurate operation.

KEEP THESE INSTRUCTIONS IN A SAFE PLACE

Important Information

- Severe dehydration and excessive water loss may cause readings which are lower than actual values. If you believe you are suffering from severe dehydration, consult a healthcare professional immediately.
- If your blood glucose results are lower or higher than usual, and you do not have symptoms of illness, first repeat the test. If you have symptoms or continue to get results which are higher or lower than usual, follow the treatment advice of your healthcare professional.
- Use only fresh whole blood samples to test your blood glucose. Using other substances will lead to inaccurate results.
- If you are experiencing symptoms that are inconsistent with your blood glucose test results and you have followed all the instructions given in this owner's manual, contact your healthcare professional.
- We do not recommend using this product on severely hypotensive individuals or patients who are in shock. Readings which are lower than actual values may occur for individuals in a hyperglycemic-hyperosmolar state, with or without ketosis. Please consult your healthcare professional before use.

Intended Use

ForaCare Test N'GO™ Blood Glucose Monitoring System is intended for self testing outside the body (*in vitro* diagnostic use) by people with diabetes at home as an aid to monitoring the effectiveness of diabetes control. It is intended to be used for the quantitative measurement of glucose (sugar) in fresh capillary whole blood samples drawn from the fingertips. This system is intended for single-patient use (lay-users at home) and should not be shared.

It should not be used for the diagnosis of or screening for diabetes mellitus, or testing on newborns.

This system provides you with plasma equivalent results.


Test Principle

Your system measures the amount of sugar (glucose) in whole blood. The glucose testing is based on the measurement of electrical current generated by the reaction of glucose with the reagent of the strip. The meter measures the current, calculates the blood glucose level, and displays the result. The strength of the current produced by the reaction depends on the amount of glucose in the blood sample.

NOTE:

Check your system to ensure that it is unopened prior to use and that it contains all parts listed above. If anything is missing or damaged, please return your system to the place of purchase.

Appearance And Key Functions of The Meter


- 1 Display Screen**
- 2 Test Slot**
- 3 Main Button**
- 4 Charging indicator/Bluetooth indicator**
- 5 SET Button**
- 6 Battery Cover**
- 7 USB Port (only for charge)**
For battery recharge.

NOTE:

The meter will turn off automatically after 180 seconds without any action or you can press and hold the main button for 3 seconds to turn off the meter.

Meter Display

1 Error Message Symbol

2 Measurement Unit

3 Testing Result Display

4 Measuring Mode

Gen, AC, PC, QC.

5 Memory Symbol

Indicates a test result stored in memory.

6 Blood Drop Symbol

Flashes when meter is ready for sample application.


7 Time Display

8 Alarm Symbol

9 Low Battery Symbol

Test Strip

Your system measures the amount of sugar (glucose) in whole blood. Blood is applied to the absorbent hole on the test strip and is automatically drawn into the reaction cell where the reaction takes place.

Setting the Meter And Deleting the Memory

Start with the meter off (no test strip inserted). Press **SET**.

1. Setting the date

With the year flashing, press Main button until the correct year appears. Press SET.

With the month flashing, press Main button until the correct month appears. Press SET.

With the day flashing, press Main button until the correct day appears. Press SET.

2. Setting the time format

Press Main button to select the desired time format — 12h or 24h. Press SET.

After you have finished selecting the 12h or 24h clock, press SET and then the hour will flash on the screen.

3. Set the hour and minute

With the hour flashing, press Main button until the correct hour appears. Press SET.

With the minute flashing, press Main button until the correct minute appears. Press SET.

4. Deleting the memory

With "dEL" and a flashing "M" on the display, press Main button and select "no" to keep the results in memory then press SET to skip.

To delete all the results, press Main button and "yes" and "M" are displayed on the meter, press SET to delete the memory.

If you do not want to delete the results, do not press the Main button but press the SET button again to turn off the meter.

5. Set the reminder alarm

You may set up any or all of the reminder alarms (1-4). The meter displays "On" or "OFF" and "M", press Main button to turn on or turn off to set the first reminder alarm.

Press Main button to select "On", then press SET to set the hour. When the hour is flashing, press Main button to add an hour. Press SET to confirm and go to minutes, press Main button to add one minute. Hold Main button longer to add

faster. Press SET to confirm and go to the next alarm setting.

If you do not want to set an alarm, press SET to skip this step.

If you want to turn off an alarm, find the alarm number by pressing SET in the setting mode, press Main button to change from “ON” to “OFF”.

6. Enter the Bluetooth pairing

With the “no” symbol flashing on the display, press SET to skip the pairing.

If you wish to enter the pairing mode, press Main button once and the meter will display “yes”. Then press SET to confirm and the meter will show “CLr”. The meter will then automatically turn off and the Bluetooth indicator will start blinking, which means the meter is in the process of pairing with a Bluetooth receiver.

NOTE:

This step is recommended when the user needs to pair this meter to a Bluetooth receiver for the first time, or when user needs to pair this meter to another new Bluetooth receiver.

CHECKING THE SYSTEM WITH CONTROL SOLUTIONS

Our control solutions contain a known amount of glucose that reacts with test strips. By comparing the result of your control solution test with the expected range printed on the test strip vial label, you can check that the meter and the test strips are working together as a system and that you are performing the test correctly. It is very important that you perform this simple check routinely to make sure you get accurate results.

Do a control solution test when:

- you first received the meter,
- at least once a week for routinely check the meter and test strips,
- you begin using a new vial of test strips,
- you suspect the meter or test strips are not working properly,
- your blood glucose test results are not consistent with how you feel, or if you think the results are not accurate,
- practicing the testing process, or
- you dropped or think you have damaged the meter.


How to Perform A Control Solution Test

Our control solutions contain a known amount of glucose that reacts with test strips and is used to ensure that the meter and test strips are working together correctly.

Test strips, control solutions, or sterile lancets may not be included in the kit (please check the contents on your product box). They can be purchased separately. Please make sure you have all items prior to testing your blood glucose. Please contact the Customer Care Line at 1-866-469-2632 or the place of purchase for availability.

1. Insert test strip


Insert a test strip into the test slot with the contact bars end first and facing up. (Contact bars must be inserted all the way into the meter or you may get an inaccurate test result.) The meter turns on automatically and displays the following in sequence:

When the “” appears on the display, press the Main button and “QC” will appear on the display. When the “QC” sign is displayed, the meter will not store your test result in the memory. If you decide not to perform a control solution test, press the Main button again and the “QC” sign will disappear.

CAUTION:

Every time you perform a control solution test you must enter “QC” test mode so that the test result will not be stored in the meter memory.

2. Apply control solution

Shake the control solution vial thoroughly before use. Squeeze out a drop and wipe it off, then squeeze out another drop and place it on the tip of the vial cap. Hold the meter to move the absorbent hole of the test strip to touch the drop. Once the confirmation window fills completely, the meter will begin counting down. To avoid contaminating the control solution, do not directly apply control solution onto a strip.

3. Read and compare the results

After the meter counts to 0, the result of the control solution test will appear. Compare the result with the range printed on the test strip vial. The result should fall within this range.

Out-of-range results

If you continue to have test results fall outside the range printed on the test strip vial, the meter and strips may not be working properly. Do **NOT** test your blood. Contact the local customer service or place of purchase for help.

TESTING YOUR BLOOD

Before testing, be sure to read this section and the test strip insert found in the test strip box carefully. Make sure you have all the items needed for the test.

The Four Measuring Modes

The meter provides you with four modes for measuring: General, AC, PC and QC.

MODES	USE WHEN
General tests	any time of day without regard to time since last meal
AC	no food intake for at least 8 hours
PC	2 hours after a meal
QC	testing with the control solution

You can switch between the two modes mode by:

1. Start with the meter switched off. Insert a test strip to turn on the meter. The screen will display a “”.

2. Press Main button to switch between General, AC, PC and QC mode.

How to Perform A Blood Glucose Test

Sharing the lancing device and lancets may increase the risk of contracting infectious diseases. Lancing device must not be used for more than one person.

To reduce the chance of infection:

- **Wash and dry your hands before starting.**
- Clean the puncture site using cotton moistened with 70% alcohol and **let it air dry.**
- Never share a lancet or the lancing device.
- Always use a new, sterile lancet. Lancets are for single use only.
- Wash and dry your hands thoroughly after handling the meter, lancing device and test strips to prevent infection. For more information, please refer to the “Cleaning and Disinfection” section.
- If the meter is being operated by a second person who is providing

testing assistance to the user, the meter and lancing device should be decontaminated prior to use by the second person.

Preparing the Puncture Site

Rubbing the puncture site before blood extraction has a significant influence on the accuracy of your test results.

Blood from an alternate a site that has not been rubbed exhibits a measurably different glucose concentration than blood from the finger. When the puncture site was rubbed prior to blood extraction, the difference was significantly reduced.

**Please follow the suggestions below before obtaining a drop of blood:
Wash and dry your hands before starting.**

- Select the puncture site on fingertips.
- Clean the puncture site using cotton moistened with 70% alcohol and **let it air dry.**
- Rub the puncture site for about 20 seconds before penetration.
- Use a clear cap while setting up the lancing device.

● **Fingertip testing**

Press the tip of the lancing device firmly against the lower side of your fingertip. Press the release button to puncture your finger, a click indicates that the puncture is complete.

NOTE:

- Choose a different spot each time you test. Repeated punctures at the same spot may cause soreness and calluses.
- It is recommended that you discard the first drop of blood as it might contain tissue fluid, which may affect the test result.

Performing A Blood Glucose Test



1. Insert the test strip to turn on the meter

Wait for the meter to display “”.

2. Select the appropriate measuring mode by pressing Main button

For selecting the measurement mode, please refer to the “FOUR MEASURING MODES”.

3. Apply blood sample

Obtain a drop of blood of at least 0.5 μL with a lancing device. Use the clear cap for sites other than fingers and refer to the strip package insert for more details. When the “” symbol appears on the display, touch the drop of blood with the absorbent hole of the test strip. Do not rub the test strip with your finger or try to apply a smeared sample.

Gently apply the drop of blood to the absorbent hole of the test strip at a titled angle. Confirmation window should be completely filled if enough blood sample has been applied. Do **NOT** remove your finger until you hear a beep sound.

If the confirmation window is not filled completely before the meter begins to count down, do not add more blood to the test strip. Discard the test strip and start again. If you have trouble filling the test strip, please call your local customer service number for assistance.

NOTE:

If you do not apply a blood sample within 3 minutes, the meter will automatically turn off. You must remove the test strip and insert it back into the meter to restart the test.

4. Get result

The result of your blood glucose test will appear after the meter counts down to 0. The blood glucose result will be stored in the memory automatically.

WARNING:

- Please do not change your treatment based on the result without first consulting your health care professional.
- Turn the meter off by removing the test strip. Discard the used test strip and lancet carefully according to your local regulations.
- Wash your hands thoroughly with soap and water after handling the meter, lancing device and test strips to avoid contamination.

VIEWING THE METER MEMORY

Your Meter stores the 450 most recent blood glucose test results with date and time in its memory. You can review the test results with these easy steps.

Viewing Results

1. Enter the memory mode

With the meter turned off, press the Main button. The first test result will appear, indicating that you are in the memory mode.


When using the meter for the first time or if the results have been deleted, “- - -” will appear, indicating that there are no test results in the memory.

2. Recall test results

After the last test result, the most recent test result appears with the date and time. Press the Main button to review your last 450 test results in order. When the memory is full, the oldest result will be deleted when the newest is added.

3. Exit the memory mode

Press the Main button for three seconds to turn off the meter or after the last result, the display will show “End” and it will turn off automatically.

NOTE:

The results of control solution tests are **NOT** stored in the memory. Stored results are blood glucose results only.

Viewing Day Average Results

1. Press and release Main button

Keep pressing Main button for 3 seconds until the flashing “AVG” appears. Release Main button and then your 7-day average result measured in general mode will appear on the display.

2. Press Main button to review

14-, 21-, 28-, 60- and 90- day average results stored in each measuring mode in the order of Gen, AC, then PC.

3. Exit the meter memory

Keep pressing the Main button and the meter will turn off after displaying the last test result.

DOWNLOADING RESULTS ONTO A COMPUTER

Data Transmission via Bluetooth

You can use your device with iOS (5.0.1 or higher) or Android system to download data from your meter via Bluetooth. Follow the steps below to transmit data from your meter. Please contact local customer services or the place of purchase for assistance.

Please note that it is required to do the pairing between meter and Bluetooth Receiver before transmitting data.

1. Make sure your meter is already paired with your device with iOS (5.0.1 or higher) or Android system by following instructions of the setting model in Step 6.

2. Install “iFORA Diabetes Manager App” e on your device with iOS (5.0.1 or higher) or Android system.

3. Data transmission

Turn on the Bluetooth on your device with iOS (5.0.1 or higher) or Android system and keep it on, and your meter will be able to connected to the device after measurement. You can view the data on your device.

WARNING:

- While the meter is in transmission mode, it will be unable to perform a blood glucose test.
- Make sure your your device with iOS (5.0.1 or higher) or Android system before transmitting the data and the meter is within the receiving range.
- The USB port of ForaCare Test N'GO™ is only for battery recharge. It does not have the function of data transmission.

CLEANING AND DISINFECTING YOUR METER

Caring for Your Meter

To avoid the meter and test strips attracting dirt, dust or other contaminants, please wash and dry your hands thoroughly before use.

Why the cleaning and disinfection should be performed

Cleaning and disinfection are different. Cleaning is the process of removing dirt (e.g. food debris, grease, dust), disinfection is the process of killing germs (e.g. bacteria and viruses).

When to clean and disinfect the meter

Clean the meter when you see any dirt on it. You should disinfect the meter at least once a week to prevent infection.

How to clean and disinfect the meter

The meter must be cleaned prior to the disinfection. Use one disinfecting wipe to clean exposed surfaces of the meter thoroughly and remove any visible dirt or blood or any other body fluid with the wipe. Use a second wipe to disinfect the meter. Do **NOT** use organic solvents to clean the meter.



We recommend for meter cleaning and disinfection you should use the disinfecting wipe/towelette from below. The following product with isopropyl alcohol as the active ingredient has been shown to be safe for use with the ForaCare Test N'GO™ glucose meter.

- ▶ Micro-Kill+™ (Micro-Kill Plus™) by Medline
(EPA Reg. No. 59894-10-37549)

To obtain disinfecting wipes and other information, please contact Medline at 1-800-MEDLINE (1-800-633-5463) or visit www.medline.com.

Disinfecting Procedures

1. Take out one disinfecting wipe from the package and squeeze out any excess liquid in order to prevent damage to the meter.

2. Wipe all meter's exterior surface display and buttons. Hold the meter with the test strip slot pointing down and wipe the area around the test slot but be careful not to allow excess liquid to get inside. Keep the meter surface wet with disinfection solution for a minimum of 2 minutes for Micro-Kill+™ wipes. Follow the instructions on the package label of disinfecting wipe.

3. Allow the meter surface to dry completely.

4. Discard the used wipes and never reuse them. Wash your hands thoroughly with soap and water after handling the meter, lancing device and test strips to avoid contamination.

Improper system cleaning and disinfection may result in meter malfunction. If you have a question, please contact customer service at 1-866-469-2632 for assistance.

This device has been validated to withstand 5,000 cleaning and disinfection cycles using the recommended disinfecting wipe/towelette. The tested number of cycles is estimated by 2 disinfection cycles per day for 5 years the expect life of device.

Stop using the meter if you see any signs of deterioration. For example:

- meter can not be turned on,
- LCD display cracks or becomes cloudy,
- buttons no longer function,
- meter outer casing cracks,
- data can not be transmitted to PC,
- color or paint/printing on housing is abnormal, or
- scratches or abrasions on meter are higher than acceptable.

Please contact the customer service for a replacement meter if any of the signs of deterioration are noticed.

NOTE :

- Do **NOT** clean and disinfect the meter while performing tests.
- If the meter is being operated by a second person, the meter and lancing device should be decontaminated prior to use by the second person.
- Do **NOT** allow cleaning and disinfecting solution to get in the test slot.
- If you do get moisture in the test strip slot, wipe it away with a corner of tissue.
- Always dry the meter thoroughly before using it.
- Do not spray the meter directly with cleaning solutions especially those containing water (i.e. soapy water), as this could cause the solution to enter the case inside and damage the electronic components or circuitry.

BATTERY RECHARGE

Your meter comes with a built-in rechargeable Li-polymer battery. If the low battery symbol “” appears on the screen indicating that the battery is low and it is time to recharge the battery.

1. Connect the USB cable to the data port of the meter.
2. Connect the other end of cable to a free USB port on your personal computer.
3. “USB” will appear on the display and the charging indicator will light up in red when battery is recharging. After the battery is fully charged, the red light will turn green.
4. Remove the USB cable and the meter will automatically switch off.

NOTE :

- It should take approximately 2 hours to be fully recharged. Recharging the battery does not affect the test results stored in the memory.
- The USB cable is for battery recharging and data transmission with a computer. If you wish to recharge through a wall outlet, please obtain a USB adapter.

DISPLAY MESSAGES AND PROBLEM-SOLVING GUIDE

The following is a summary of display messages. If your meter displays an error message, please follow the actions for the error message as described in the table below. If the problem persists, please contact your local customer service agent for help.

MESSAGE	WHAT IT MEANS	ACTION
	Appears when glucose < 20 mg/dL (1.1 mmol/L).	This indicates hypoglycemia (low blood glucose.) Please seek immediate medical assistance.
	Appears when glucose > 600 mg/dL (33.3mmol/L).	This indicates severe hyperglycemia (high blood glucose). Please seek immediate medical assistance.
	The “  ” symbol appears with this message. This means that the battery does not have enough power remaining for a test.	Replace the battery immediately.
	A used test strip is inserted, or there is an electrical problem.	Repeat test with a new test strip. If the error message re-appears, please contact your local customer service agent for help.
	You may have removed the strip after applying blood to the absorbent hole.	Keep the test strip in the meter after it has absorbed the blood until the meter displays the test result.

MESSAGE	WHAT IT MEANS	ACTION
	<p>The temperature of the environment, meter, or test strip is outside the system's operating range. You cannot perform a test until the meter and test strip are within the operating range of 50°F to 104°F (10°C to 40°C).</p>	<p>Repeat the test after the meter and test strip are within the operating temperatures.</p>
	<p>Problem with the meter.</p>	<p>Review the instructions and try again with a new test strip. If the problem persists, please contact your local customer service agent for help.</p>

For questionable results, please see test strip package insert for important information.

Troubleshooting

1. If the meter does not display a message after inserting a test strip:

POSSIBLE CAUSE	WHAT TO DO
Batteries exhausted.	Replace the batteries.
Test strip inserted upside down or incompletely.	Insert the test strip with contact bars end first and facing up.
Defective meter or test strips.	Please contact customer service.

2. If the test does not start counting down after applying the sample:

POSSIBLE CAUSE	WHAT TO DO
Insufficient blood sample.	Repeat the test using a new test strip with larger volume of blood sample.
Defective test strip.	Repeat the test with a new test strip.
Sample applied after automatic switch-off (3 minutes after last user action).	Repeat the test with a new test strip. Apply sample only when flashing “  ” appears on the display.
Defective meter.	Please contact customer service.

3. If the control solution testing result is out of range:

POSSIBLE CAUSE	WHAT TO DO
Error in performing the test.	Read instructions thoroughly and repeat the test again.
Control solution vial was poorly shaken.	Shake the control solution vigorously and repeat the test again.
Expired or contaminated control solution.	Check the expiry date of the control solution.
Control solution that is too warm or too cold.	Control solution, meter, and test strips should be at room temperature (20°C to 25°C/68°F to 77°F) before testing.
Defective test strip.	Repeat the test with a new test strip.
Meter malfunction.	Please contact customer service.

DETAILED INFORMATION

Reference Values

Blood Glucose

Blood glucose monitoring plays an important role in diabetes control. A long-term study showed that maintaining **blood glucose levels close to normal** can reduce the risk of diabetes complications by up to 60%*¹. The results provided by this system can help you and your healthcare professional monitor and adjust your treatment plan to gain better control of your diabetes.

Time of day	Normal plasma glucose range for people without diabetes (mg/dL)
Fasting and before meal	< 100 mg/dL (5.6 mmol/L)
2 hours after meals	< 140 mg/dL (7.8 mmol/L)

Source: American Diabetes Association. Standards of Medical Care in Diabetes – Diabetes Care 2012; 35 (Suppl 1):S1-S100.

Please work with your doctor to determine a target range that works best for you.

References:

*¹: American Diabetes Association position statement on the Diabetes Control and Complications Trial (1993).

Comparing Meter and Laboratory Results

The meter provides you with whole blood equivalent results. The result you obtain from your meter may differ somewhat from your laboratory result due to normal variation. Meter results can be affected by factors and conditions that do not affect laboratory results in the same way. To

make an accurate comparison between meter and laboratory results, follow the guidelines below.

Before going to the lab:

- Perform a control solution test to make sure that the meter is working properly.
- Fast for at least eight hours before doing comparison tests, if possible.
- Take your meter with you to the lab.

While staying at the lab:

Make sure that the samples for both tests are taken and tested within 15 minutes of each other.

- Wash your hands before obtaining a blood sample.
- Never use your meter with blood that has been collected in a gray-top test tube.
- Use fresh capillary blood only.

You may still have a variation from the result because blood glucose levels can change significantly over short periods of time, especially if you have recently eaten, exercised, taken medication, or experienced stress^{*2}. In addition, if you have eaten recently, the blood glucose level from a finger prick can be up to 70 mg/dL (3.9 mol/L) higher than blood drawn from a vein (venous sample) used for a lab test^{*3}. Therefore, it is best to fast for eight hours before doing comparison tests. Factors such as the amount of red blood cells in the blood (a high or low hematocrit) or the loss of body fluid (dehydration) may also cause a meter result to be different from a laboratory result.

References:

*2: Surwit, R.S., and Feinglos, M.N.: Diabetes Forecast (1988), April, 49-51.

*3: Sacks, D.B.: "Carbohydrates." Burtis, C.A., and Ashwood, E.R. (ed.), Tietz Textbook of Clinical Chemistry. Philadelphia: W.B. Saunders Company (1994), 959.

SPECIFICATIONS

Dimension: 93.0 (L) x 26.03 (W) x 5.57 (H) mm

Weight: 26.8 g (without battery)

Power Source: Li-Polymer battery

Display: LCD

Memory: 450 measurement results with date and time

Automatic detection of electrode insertion

Automatic reaction time count-down

Automatic blood detection

Temperature warning

Operating Condition:

50 °F to 104 °F (10 °C to 40 °C), below 85% R.H.

Storage/Transportation Conditions:

-4 °F to 140 °F (-20 °C to 60 °C), below 95% R.H.

Measurement Units: mg/dL

Measurement Range: 20 to 600 mg/dL (1.1 to 33.3 mmol/L)

This device has been tested to meet the electrical and safety requirements of: IEC/EN 61010-1, IEC/EN 61010-2-101, EN 61326-1, EN 61326-2-6.

Statement

15.21

You are cautioned that changes or modifications not expressly approved by the part responsible for compliance could void the user's authority to operate the equipment.

15.105(b)

Federal Communications Commission (FCC) Statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

- 1) This device may not cause harmful interference and
- 2) This device must accept any interference received, including interference that may cause undesired operation of the device.

FCC RF Radiation Exposure Statement:

- 1) This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.
- 2) This equipment complies with FCC RF radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with a minimum distance of 20 centimeters between the radiator radiation source and your body.

"Made for iPod", "Made for iPhone", and "Made for iPad" mean that an electronic accessory has been designed to connect specifically to iPod, iPhone or iPad, respectively, and has been certified by the developer to meet Apple performance standards. Apple is not responsible for the operation of this device or its compliance with safety and regulatory standards. Please note that the use of this accessory with iPod, iPhone or iPad may affect wireless performance.

iPod touch®, iPhone® and iPad® are trademarks of Apple Inc., registered in the U.S. and other countries.

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1. “FDA Public Health Notification: Use of Fingertick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens: Initial Communication” (2010)

<http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/ucm224025.htm>

2. “CDC Clinical Reminder: Use of Fingertick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens” (2010)

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Fora Care Inc.