
JC-604 (TMB-1581-S)

JOYCARE
PHARMA **JL**

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GB

Thank you very much for selecting JOYCARE Blood Pressure Monitor.

IMPORTANT INFORMATION! PLEASE KEEP!

Carefully read these user instructions and in particular the safety information prior to using this device and keep these user instructions for future use. When handing the device over to third parties, please also include these user instructions. Non-compliance with these instructions may cause injuries or damage the device.

Table of Contents

INTRODUCTION	42
General Description	
Indications for Use	
Contraindications	
Measurement Principle	
Safety Information	
LCD Display Signal	
Monitor Components	
List	

BEFORE YOU START	53
Installing and Replacing the Batteries	
Setting Date, Time and Measurement Unit	
MEASUREMENT	57
Tie the Cuff	
Start the Measurement	
DATA MANAGEMENT	60
Recall the Records	
Delete the Records	
INFORMATION FOR USER	63
Tips for measurement	
Maintenance	
ABOUT BLOOD PRESSURE	64
What are systolic pressure and diastolic pressure?	
What is the standard blood pressure classification?	
Irregular heartbeat detector	
Why does my blood pressure fluctuate throughout the day?	

Why do I get a different blood pressure at home compared to the hospital?
Is the result the same if measuring on the right wrist?

SPECIFICATIONS	69
TROUBLESHOOTING	72
COMPLIED STANDARDS LIST	74
EMC GUIDANCE	115

INTRODUCTION

General Description

Thank you for selecting JOYCARE wrist type blood pressure monitor.

Readings taken by this blood pressure monitor are equivalent to those obtained by a trained observer using the cuff and stethoscope auscultation method.

This manual contains important safety and care information, and provides step by step instructions for using the product.

Read the manual thoroughly before using the product.

Features:

Systolic blood pressure

Diastolic blood pressure

Pulse rate

60 records for one user

Indications for Use

The Blood Pressure Monitor is digital monitors intended for use in measuring blood pressure and heartbeat rate with a wrist circumference ranging from 13.5cm to 21.5 cm (about 5 $\frac{1}{3}$ " -8 $\frac{1}{2}$ ").It is intended for adult indoor use only.

Contraindications

- 1.The device should not be used by any person who may be suspected of, or is pregnant .
- 2.The device is not suitable for use on patients with implanted, electrical devices, such as cardiac pacemakers, defibrillators.

Measurement Principle

This product uses the Oscillometric Measuring method to detect blood pressure.

Before every measurement, the unit establishes a "zero pressure" equivalent to the atmospheric pressure. Then it starts inflating the wrist cuff, meanwhile, the unit detects pressure oscillations generated by beat-to-beat pulsatile, which is used to determine the systolic and diastolic pressure, and also pulse rate.

Safety Information

The signs below might be in the user manual, labelling or other component.

	Symbol for "THE OPERATION GUIDE MUST BE READ"		Symbol for "TYPE BF APPLIED PARTS"
CE 0123	Symbol for "COMPLIES WITH MDD 93/42/EEC REQUIREMENTS"		Symbol for "ENVIRONMENT PROTECTION - Electrical waste products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local authority or retailer for recycling advice"
	Symbol for "MANUFACTURER"		
	Symbol for "SERIAL NUMBER"		Symbol for "DIRECT CURRENT"
	Symbol for "MANUFACTURE DATE"		Symbol for "Authorised Representative in the European Community"
	Caution: These notes must be observed to prevent any damage to the device.		

CAUTION

This device is intended for adult use in homes only.

The device is not suitable for use on neonatal patients, pregnant women, patients with implanted, electronical devices, patients with pre-eclampsia, premature ventricular beats, atrial fibrillation, peripheral, arterial disease and patients undergoing intravascular therapy or arterio-venous shunt or people who received a mastectomy. Please consult your doctor prior to using the unit if you suffer from illnesses.

The device is not suitable for measuring the blood pressure of children. Ask your doctor before using it on older children.

The device is not intended for patient transport outside a healthcare facility.

The device is not intended for public use.

This device is intended for no-invasive measuring and monitoring of arterial blood pressure. It is not intended for use on extremities other than the wrist or for functions other than obtaining a blood pressure measurement.

Do not confuse self-monitoring with self-diagnosis. This unit allows you to monitor your blood pressure. Do not begin or end medical treatment without asking a physician for treatment advice.

If you are taking medication, consult your physician to determine the most appropriate time to measure your blood pressure. Never change a prescribed medication without consulting your physician.

Do not take any therapeutic measures on the basis of a self measurement. Never alter the dose of a medicine prescribed by a doctor. Consult your doctor if you have any question about your blood pressure.

When the device was used to measure patients who have common arrhythmias such as atrial or ventricular premature beats or atrial fibrillation, the best result may occur with deviation. Please consult your physician about the result.

Warning: Do not apply the cuff over a wound; otherwise it can cause further injury.

Do not inflate the cuff on the same limb which other monitoring ME equipment is applied

around simultaneously, because this could cause temporary loss of function of those simultaneously-used monitoring ME equipment.

On the rare occasion of a fault causing the cuff to remain fully inflated during measurement, open the cuff immediately. Prolonged high pressure (cuff pressure > 300mmHg or constant pressure > 15mmHg for more than 3 minutes) applied to the wrist may lead to an ecchymosis.

Please check that operation of the device does not result in prolonged impairment of patient blood circulation.

When measurement, please avoid compression or restriction of the connection tubing.

The device cannot be used with HF surgical equipment at the same time.

The ACCOMPANYING DOCUMENT shall disclose that the SPHYGMOMANOMETER was clinically investigated according to the requirements of ISO 81060-2:2013.

To verify the calibration of the AUTOMATED SPHYGMOMANOMETER, please contact the manufacturer.

This device is contraindicated for any female who may be suspected of, or is pregnant. Besides providing inaccurate readings, the effects of this device on the fetus are unknown. Too frequent and consecutive measurements could cause disturbances in blood circulation and injuries.

This unit is not suitable for continuous monitoring during medical emergencies or operations. Otherwise, the patient's wrist and fingers will become anaesthetic, swollen and even purple due to a lack of blood.

When not in use, store the device in a dry room and protect it against extreme moisture, heat, lint, dust and direct sunlight. Never place any heavy objects on the storage case. This device may be used only for the purpose described in this booklet. The manufacturer cannot be held liable for damage caused by incorrect application.

This device comprises sensitive components and must be treated with caution. Observe the storage and operating conditions described in this booklet.

The maximum temperature that the applied part can be achieved is 42.5°C while the environmental temperature is 40°C.

The equipment is not AP/APG equipment and not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.

Warning: No servicing/maintenance while the ME equipment is in use.

The patient is an intended operator.

The patient can measure data and change battery under normal circumstances and maintain the device and its accessories according to the user manual.

To avoid measurement errors, please avoid the condition of strong electromagnetic field radiated interference signal or electrical fast transient/burst signal.

The blood pressure monitor and the cuff are suitable for use within the patient environment. If you are allergic to polyester, nylon or plastic, please don't use this device. During use, the patient will be in contact with the cuff. The materials of the cuff have been tested and found to comply with requirements of ISO 10993-5:2009 and ISO 10993-10:2010. It will not cause any potential sensitization or irritation reaction.

If you experience discomfort during a measurement, such as pain in the wrist or other complaints, press the START/STOP button to release the air immediately from the cuff. Loosen the cuff and remove it from your wrist.

If the cuff pressure reaches 40 kPa (300 mmHg), the unit will automatically deflate. Should the cuff not deflate when pressures reaches 40 kPa (300 mmHg), detach the cuff from the wrist and press the START/STOP button to stop inflation.

Before use, make sure the device functions safely and is in proper working condition. Check the device, do not use the device if it is damaged in any way. The continuous use of a damaged unit may cause injury, improper results, or serious danger.

Do not wash the cuff in a washing machine or dishwasher!

The service life of the cuff may vary by the frequency of washing, skin condition, and storage state.

The typical service life is 10000 times.

It is recommended that the performance should be checked every 2 years and after maintenance and repair, by retesting at least the requirements in limits of the error of the cuff pressure indication and air leakage (testing at least at 50mmHg and 200mmHg). Please dispose of ACCESSORIES, detachable parts, and the ME EQUIPMENT according to the local guidelines.

Manufacturer will make available on request circuit diagrams, component part lists, descriptions, calibration instructions, etc., to assist to service personnel in parts repair.

The operator shall not touch output of batteries and the patient simultaneously.

Cleaning: dust environment may affect the performance of the unit. Please use the soft cloth to clean the whole unit before and after use. Don't use any abrasive or volatile cleaners.

The device doesn't need to be calibrated within two years of reliable service.

If you have any problems with this device, such as setting up, maintaining or using, please contact the SERVICE PERSONNEL of Joycare. Don't open or repair the device by yourself in the event of malfunctions. The device must only be serviced, repaired and opened by individuals at authorized sales/service centres.

Please report to Joycare if any unexpected operation or events occur.

Keep the unit out of reach of infants, young children or pets to avoid inhalation or swallowing of small parts. It is dangerous or even fatal.

Be careful to strangulation due to cables and hoses, particularly due to excessive length. At least 30 min required for ME equipment to warm from the minimum storage temperature between uses until it is ready for intended use. At least 30 min required for ME equipment to cool from the maximum storage temperature between uses until it is ready for intended use.

This equipment needs to be installed and put into service in accordance with the information provided in the ACCOMPANYING DOCUMENTS.

Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies can affect this equipment and should be kept at least a distance d away from the equipment. The

distance d is calculated by the MANUFACTURER from the 80 MHz to 5.8 GHz column of Table 4 and Table 9 of IEC 60601-1-2:2014, as appropriate.


Please use ACCESSORIES and detachable parts specified/ authorised by MANUFACTURE.

Otherwise, it may cause damage to the unit or danger to the user/patients.

There is no luer lock connectors are used in the construction of tubing, there is a possibility that they might be inadvertently connected to intravascular fluid systems, allowing air to be pumped into a blood vessel.


Please use the device under the environment which was provided in the user manual. Otherwise, the performance and lifetime of the device will be impacted and reduced.

LCD display signal

SYMBOL	DESCRIPTION	EXPLANATION
SYS	Systolic blood pressure	High blood pressure
DIA	Diastolic blood pressure	Low blood pressure
	Pulse display	Pulse in beats per minute
	Excessive Body Motion Detector	Motion may result in an inaccurate measurement.
	Low battery	Batteries are running low and need to be replaced
kPa	kPa	Measurement unit of the blood pressure (1 kPa=7.5mmHg)
mmHg	mmHg	Measurement unit the blood pressure (1 mmHg=0.133kPa)
	Irregular heartbeat	Blood pressure monitor is detecting an irregular heartbeat during measurement.
	Current time	Month/Day/Year, Hour : Minute
	Blood pressure level indicator	Indicate the blood pressure level
	Heartbeat	Blood pressure monitor is detecting a heartbeat during measurement.

Monitor Components

- 1. LCD DISPLAY
- 2. GRADE
- 3. MEM BUTTON SET BUTTON
- 4. SYSTOLIC
- 5. DIASTOLIC
- 6. PULSE RATE
- 7. START/STOP BUTTON
- 8. CUFF (Type BF applied part)
- 9. BATTERY COMPARTMENT

Component List pressure measuring system:

PCBA
Air Pipe
Pump
Valve
Cuff

List

Wrist Blood Pressure Monitor
2×AAA batteries
User manual


BEFORE YOU START

Installing and Replacing the Batteries
Slide off the battery cover.

Install the batteries by matching the correct polarity, as shown below. Always use the correct battery type (2 x AAA batteries).

Replace the battery cover.

Replace the batteries whenever the below happen:

The  shows
The display is dim.
The display does not light up

CAUTION


Information about handling batteries
Do not mix new and used batteries.
Do not use different types of batteries at the same time.
Do not throw batteries into a fire. Batteries may explode or leak.
Remove the batteries if the device is not used for a long time.

Battery disposal


Worn batteries are harmful to the environment. Do not dispose with daily garbage.
Remove the old batteries from the device following your local recycling guidelines

Setting Date, Time and Measurement Unit


It is important to set the clock before using your blood pressure monitor, so that a time stamp can be assigned to each record that is stored in the memory. (year: 2014 —2054,time format: 12H /24 H).



1. When the monitor is off, hold pressing “MEM” button for about 3 seconds to enter into setting mode. The blinking numeral represents [YEAR].



2. Press the “MEM” button to change the [YEAR]. Each press will increase the numeral by one in a cycling manner.



3. When you get the right year, press “START/STOP” button to confirm your selection and it will turn to the next step.



4. Repeat steps 2 and 3 to confirm [MONTH] and [DAY]..

5. Repeat steps 2 and 3 to confirm the time format [12H] and [24H].

6. Repeat steps 2 and 3 to confirm [HOUR] and [MINUTE].

7. Repeat steps 2 and 3 to confirm the [MEASUREMENT UNIT].

After confirming the measurement unit, the LCD will display all the settings you have done once again and then turn off.

MEASUREMENT

Tie the Cuff

Remove all accessories (watch, bracelet, etc) from your wrist. If your physician has diagnosed you with poor circulation in your wrist, use the other one.

Roll or push up your sleeve to expose the skin.

Apply the cuff to your wrist with your palm facing up.

Position the edge of the cuff about 1cm~1.5cm from wrist joints.

Fasten the wrist cuff around your wrist, leaving no extra room between the cuff and your skin. If the cuff is too loose, the measurement will not be accurate.

Sit comfortably with your tested wrist resting on a flat surface. Place your elbow on a table so that the cuff is at the same level as your heart. Turn your palm upwards. Sit upright in a chair, and take 5-6 deep breaths.

Patients with Hypertension:

The middle of the cuff should be at the level of the right atrium of the heart; Before starting measurement, please sit comfortably with legs uncrossed, feet flat on the floor, back and wrist supported.

Rest for 5 minutes before measuring.

Wait at least 3 minutes between measurements. This allows your blood circulation to recover.

Take the measurement in a silent room.

The patient must relax as much as possible and do not move and talk during the measurement procedure

The cuff should maintain at the same level as the right atrium of the heart.

Do not cross your legs and keep your feet on the ground.

Keep your back against the backrest of the chair.

For a meaningful comparison, try to measure under similar conditions.


For example, take daily measurements at approximately the same time, on the same wrist, or as directed by a physician.

Start the Measurement


1. When the monitor is off, press "START/STOP" button to turn on the monitor, and it will finish the whole measurement.

LCD Display

Adjust the zero

Inflating and measuring


Display and save the result


2. Press "START/STOP" button to power off, otherwise it will turn off within 1 minute.

DATA MANAGEMENT

Recall the Recor

1. When the monitor is off, press "MEM" button to show the average value of the latest three measurement records. If the records are less than 3 groups, it will display the latest record instead.

2. Press "MEM" button again, it will display the latest measurement result, date and time. Press "MEM" button again, it will display the next record, and so on. During the process of displaying the results, if there is no operation, the blood pressure monitor will turn off in one minute. Or you can press "START/STOP" button to turn it off.

The record number, date and time will display alternately.

1/5

It means the total records is 5, the current No. is No. 1.

5/14

The corresponding date is May 14th

8:16

The corresponding time is 8:16.


CAUTION

The most recent record (1) is shown first. Each new measurement is assigned to the first (1) record. All other records are pushed back one digit (e.g., 2 becomes 3, and so on), and the last record (60) dropped from the list.

Delete the Records

If you did not get the correct measurement, you can delete all results by

following steps below.


1. Hold pressing "MEM" button for 3 seconds when the monitor is in the memory recall mode, the "dEL All" will flash on the display.

2. Press "MEM" to confirm deleting, the LCD displays "dEL dOnE" and the monitor will turn off.

Note: To exit out of delete mode without deleting any records, press START/STOP button before pressing "MEM" to confirm any delete commands.

3. If there is no record, the right display will show.

MAINTENANCE AND CARE OF THE DEVICE

Storage information

If the device is not being used, store it in a dry location away from heat, dust, fibres, lint, and out of direct sunlight. Do not place heavy objects on the storage box.

Clean the device using a cloth dampened with water or a delicate detergent, and dry using a supple cloth. Never use aggressive detergents or sprays to clean the device. Never immerse the wrist cuff or the device in water.

INFORMATION FOR USER

Tips for Measurement

Measurements may be inaccurate if taken in the following circumstances.

Within 1 hour after dinner or drinking

Immediate measurement after tea, coffee, smoking

Within 20 minutes after taking a bath

When talking or moving your fingers

In a very cold environment

When you want to discharge urine

Maintenance

In order to get the best performance, please follow the instructions below.

Put in a dry place and avoid the sunshine

Avoid touching water, clean it with a dry cloth in case.

Avoid intense shaking and collisions

Avoid dusty and unstable temperature environment


Using wet cloths to remove dirt

Do not attempt to clean the reusable cuff with water and never immerse the cuff in water.

ABOUT BLOOD PRESSURE

What are systolic pressure and diastolic pressure?


When ventricles contract and pump blood out of the Diastolic heart, the blood pressure reaches its maximum value in the cycle, which is called systolic pressure. When the ventricles relax, the blood pressure reaches its minimum value in the cycle, which is called diastolic pressure.

What is the standard blood pressure classification?

The blood pressure classification published by World Health Organization (WHO) and

International Society of Hypertension (ISH) in 1999 is as follows:

CAUTION

Only a physician can tell your normal BP range. Please contact a physician if your measuring result falls out of the range. Please note that only a physician can tell whether your blood pressure value has reached a dangerous point.

Level Blood Pressure (mm Hg)	Optimal	Normal	High-normal	Mild	Moderate	Severe
SYS	<120	120-129	130-139	140-159	160-179	180
DIA	<80	80-84	85-89	90-99	100-109	110

Irregular Heartbeat Detector

An irregular heartbeat is detected when a heartbeat rhythm varies while the device is measuring systolic pressure and diastolic pressure. During each measurement, blood

pressure monitor will keep a record of all the pulse intervals and calculate the average value of them. If there are two or more pulse intervals , the difference between each interval and the average is more than the average value of $\pm 25\%$, or there are four or more pulse intervals ,the difference between each interval and the average is more than the average value of $\pm 15\%$, then the irregular heartbeat symbol will appear on the display with the measurement result.

CAUTION

The appearance of the IHB icon indicates that a pulse irregularity consistent with an irregular heart-beat was detected during measurement. Usually this is NOT a cause for concern. However, if the symbol appears often, we recommend you seek medical advice. Please note that the device does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

Why does my blood pressure fluctuate throughout the day?

1. Individual blood pressure varies multiple times every day. It is also affected by the way you tie your cuff and your measurement position, so please take the measurement under the same conditions.
- 2.If the person takes medicine, the pressure will vary more.
- 3.Wait at least 3 minutes for another measurement.

Why do I get a different blood pressure at home compared to the hospital?

The blood pressure is different even throughout the day due to weather, emotion, exercise etc, Also, there is the “white coat” effect, which means blood pressure usually increases in clinical settings.

What you need to pay attention to when you measure your blood pressure at home:

If the cuff is tied properly.

If the cuff is too tight or too loose.

If the cuff is tied on the wrist.

If you feel anxious.

Taking 2-3 deep breaths before beginning will be better for measuring.

Advice: Relax yourself for 4-5 minutes until you calm down.

Is the result the same if measuring on the right wrist?

It is ok for both wrists, but there will be some different results for different people. We suggest you measure the same wrist every time.

MAINTENANCE AND CARE OF THE DEVICE

Storage information

If the device is not being used, store it in a dry location away from heat, dust, fibres, lint, and out of direct sunlight. Do not place heavy objects on the storage box.

Clean the device using a cloth dampened with water or a delicate detergent, and dry using a supple cloth. Never use aggressive detergents or sprays to clean the device. Never immerse the cuff or the device in water.

DISPOSAL

The device (including its removable parts and accessories) must not be disposed of together with municipal waste at the end of its life, but in compliance with European Directive. Since it must be handled separately from household waste, either carry it to a separately-collected waste disposal centre for electrical and electronic appliances or give it back to the retailer on purchasing a new device with the same purpose. Any infringement will be severely prosecuted.

The batteries used in this device must be disposed of in the special bins at the end of their life.

All the instructions and drawings are based on the latest information available at the time of going to print and may be subject to change.

CE0123 This appliance has been manufactured in compliance with the European directives 93/42/ECC

Guangdong Transtek Medical Electronics Co., Ltd. - Zone A, No.105 ,Dongli Road, Torch Development District, Zhongshan, 528437, Guangdong, China

MDSS - Medical Device Safety Service GmbH - Schiffgraben 41, 30175 Hannover, Germany

Distributed under the brand JOYCARE, by JOYCARE S.R.L. sole partner, Via Massimo D'Antona, 28 - 60033 Chiaravalle (AN) - Italy

SPECIFICATIONS

Model	JC-604 (TMB-1581-S)
Description	Wrist blood pressure monitor
Power supply	Battery powered mode: 2*AAA batteries

Display mode	Digital LCD V.A.32mmx45mm
Measurement mode	Oscillographic testing mode
Measurement range	Rated cuff pressure: 0mmHg~299mmHg(0kPa ~ 39.9kPa) Measurement pressure: SYS: 60mmHg~230mmHg (8.0kPa~30.7kPa) DIA: 40mmHg~130mmHg (5.3kPa~17.3kPa) Pulse value: (40-199)beat/minute
Accuracy	Pressure:5°C-40°C within ± 3 mmHg(0.4kPa)
Working condition	A temperature range of :+5°C to +40°C A relative humidity range of 15% to 90%, non-condensing, but not requiring a water vapor partial pressure greater than 50 hPa An atmospheric pressure range of : 700 hPa to 1060 hPa
Storage & transportation condition	Temperature:-20°C to +60°C A relative humidity range of $\leq 93\%$, non-condensing, at a water vapour pressure up to 50hPa
Measurement perimeter of the wrist	About 13.5cm-21.5cm

Weight	Approx.104g(Excluding the batteries)
External dimensions	Approx.85mm×67mm×23mm(Excluding the cuff)
Attachment	2*AAA batteries, user manual
Mode of operation	Continuous operation
Degree of protection	Type BF applied part
Device Classification	Internally Powered ME Equipment
IP Classification	IP22: The first number 2: Protected against solid foreign objects of 12,5mm Φ and greater. The second number: Protected against vertically falling water drops when enclosure titled up to 15°. Vertically falling drops shall have no harmful effects when the enclosure is titled at any angle up to 15° on either side of the vertical.
Software Version	A01

WARNING: No modification of this equipment is allowed

TROUBLE SHOOTING

This section includes a list of error messages and frequently asked questions for problems you may encounter with your blood pressure monitor. If the product is not operating as you think it should, check here before arranging for servicing.

PROBLEM	SYMPTOM	CHECK THIS	REMEDY
No power	Display will not light up.	Batteries are exhausted	Replace with new batteries
		Batteries are inserted incorrectly	Insert the batteries correctly
Low batteries	Display is dim or Show 	Batteries are low.	Replace with new batteries
Error message	E 1 shows	The inflation is too low.	Refasten the cuff and then measure again.
	E 2 shows	The cuff is not secure.	Refasten the cuff and then measure again.
	E 3 shows	Can't inflate the air during the process of measurement.	Refasten the cuff and then measure again.

Error message	E 4 shows	The cuff is too tight	Refasten the cuff and then measure again.
	E 5 shows	The pressure of the cuff is excess.	Relax for a moment and then measure again.
	E 6 shows	The voltage of the batteries is too low during the inflation.	Measure again . If E6 appears again, please replace with new batteries.
	E 10 or E11 shows	The monitor detected motion, talking while measuring.	Movement can affect the measurement. Relax for a moment and then measure again.
	E 20 shows	The measurement process does not detect the pulse signal	Loosen the clothing on the arm and then measure again.
	E 21 shows	Measure incorrectly.	Relax for a moment and then measure again.

Error message	EExx, shows on the display.	A calibration error occurred.	Retake the measurement. If the problem persists, contact the retailer or our customer service department for further assistance. Refer to the warranty for contact information and return instructions.
Warning message	"out " shows	Out of measurement range	Relax for a moment. Refasten the cuff and then measure again. If the problem persists, contact your physician.

COMPLIED STANDARDS LIST

Risk management	EN ISO 14971:2012 / ISO 14971:2007 Medical devices - Application of risk management to medical devices
Labeling	EN ISO 15223-1:2016 / ISO 15223-1:2016 Medical devices. Symbols to be used with medical device labels, labeling and information to be supplied. Part 1 : General requirements

User manual	EN 1041:2008 Information supplied by the manufacturer of medical devices
General Requirements for Safety	EN 60601-1:2006+A1:2013/ IEC 60601-1:2005+A1:2012 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance EN 60601-1-11:2015/ IEC 60601-1-11:2015 Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment
Electromagnetic compatibility	EN 60601-1-2:2015/ IEC 60601-1-2:2014 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic disturbances - Requirements and tests
Performance requirements	EN ISO 81060-1:2012 Non-invasive sphygmomanometers - Part 1: Requirements and test methods for non-automated measurement type EN 1060-3:1997+A2:2009 Non-invasive sphygmomanometers - Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems IEC 80601-2-30:2009+A1:2013 Medical electrical equipment- Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers

Clinical investigation	EN 1060-4:2004 Non-invasive sphygmomanometers - Part 4: Test procedures to determine the overall system accuracy of automated non-invasive sphygmomanometers ISO 81060-2:2013 Non-invasive sphygmomanometers - Part 2: Clinical validation of automated measurement type
Usability	EN 60601-1-6:2010+A1:2015/IEC 60601-1-6:2010+A1:2013 Medical electrical equipment - Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability IEC 62366-1:2015 Medical devices - Part 1: Application of usability engineering to medical devices
Software life-cycle processes	EN 62304:2006/AC: 2008 / IEC 62304: 2006+A1:2015 Medical device software - Software life-cycle processes
Bio-compatibility	ISO 10993-1:2009 Biological evaluation of medical devices- Part 1: Evaluation and testing within a risk management process ISO 10993-5:2009 Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity ISO 10993-10:2010 Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization

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EMC Guidance

- 1) This product needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided, and this unit can be affected by portable and mobile RF communications equipment.
- 2) Do not use a mobile phone or other devices that emit electromagnetic fields, near the unit. This may result in incorrect operation of the unit.
- 3) Caution: This unit has been thoroughly tested and inspected to assure proper performance and operation!
- 4) Caution: This machine should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this machine should be observed to verify normal operation in the configuration in which it will be used.

Table 1

Guidance and manufacturer's declaration – electromagnetic emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance

RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Table 2

Guidance and manufacturer's declaration – electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance

Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	power supply lines: ± 2 kV input/output lines: ± 1 kV	power supply lines: ± 2 kV	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC61000-4-5	line(s) to line(s): ± 1 kV line(s) to earth: ± 2 kV 100 kHz repetition frequency	line(s) to line(s): ± 1 kV 100 kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.

Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT; 1 cycle and 70% UT; 25/30 cycles Single phase: at 0° 0% UT; 300 cycle	0% UT; 0.5 cycle At 0°, 45°, 90°, 35°, 180°, 225°, 270° and 315° 0% UT; 1 cycle and 70% UT; 25/30 cycles Single phase: at 0° 0% UT; 300 cycle	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE UT is the a.c. mains voltage prior to application of the test level.			

Table 3

Guidance and manufacturer's declaration – electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity test	EC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	150 kHz to 80 MHz: 3 Vrms 6Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	150 kHz to 80 MHz: 3 Vrms 6Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter. Recommended separation distances: $d=0.35 \sqrt{P}$; $d=1.2 \sqrt{P}$

Radiated RF IEC 61000- 4-3	10V/m, 80% Am at 1kHz	10V/m, 80% Am at 1kHz	80 MHz to 800 MHz: $d=1.2$ $\sqrt{P}$ 800 MHz to 2.7 GHz: $d=2.3$ $\sqrt{P}$	where, “P” is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer, <i>d</i> is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, “a” should be less than the compliance level in each frequency range. “b” Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.				

"a" - Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

"b" Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

Table 4

Recommended separation distances between portable and mobile RF communications equipment and the device.

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 3.5 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.7 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

Table 5

Guidance and manufacturer's declaration - electromagnetic immunity							
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device, should assure that it is used in such an environment.							
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	Test Frequency (MHz)	Band a) (MHz)	Service a)	Modulation b)	Modulation b) (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
	385	380-390	TETRA 400	Pulse modulation b) 18Hz	1.8	0.3	27
	450	380-390	GMRS 460 FRS 460	FM c) ± 5kHz deviation 1kHz sine	2	0.3	28
	710	704-787	LTE Band 13, 17	Pulse modulation b) 217Hz	0.2	0.3	9
	745						
	780						

	810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation b) 18Hz	2	0.3	28
	870						
	930						
	1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation b) 217Hz	2	0.3	28
	1845						
	1970						
	2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation b) 217Hz	2	0.3	28

	5240	5100-5800	WLAN 802.11 a/n	Pulse modulation b) 217Hz	0.2	0.3	9
	5240						
	5785						

NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

a) For some services, only the uplink frequencies are included.
b) The carrier shall be modulated using a 50% duty cycle square wave signal.
c) As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

The MANUFACTURER should consider reducing the minimum separation distance, based on RISK MANAGEMENT, and using higher IMMUNITY TEST LEVELS that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher IMMUNITY TEST LEVELS shall be calculated using the following equation:

$$E = \frac{6}{d} \sqrt{P}$$

Where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVEL in V/m.

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GB: WARRANTY

Joycare S.r.l., single-member company, having its registered office in Via Massimo D'Antona 28, 60033 Chiaravalle (AN), guarantees this product for a period of 24 months from the date of delivery.

The consumer shall have the rights, provided by the applicable national laws, ruling the sales of consumer goods. This warranty shall leave such rights unaffected.

To require intervention during warranty period, the purchaser, as soon as detecting the flaw, must contact the retailer/local importer by written means or over the telephone.

In order to obtain assistance during the warranty period, it is necessary to return the faulty product to the retailer/local importer, together with the warranty slip, duly and thoroughly filled out, the purchase invoice and delivery note (purchase receipt, invoice or receipt, bill of lading).

However, the warranty slip alone does not replace in any way the document, proving the delivery date of the product.

Note: any attempt of repairing, recovering or changing the product, beyond the provisions of the user's manual, by an individual, not authorized by the seller, shall invalidate the warranty.

We recommend to carefully and thoroughly read the instructions manual of the product.

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Santion , Importator: S.C. PERFECT MEDICAL S.R.L.

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