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
pressure logic

professional KS 551

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
**GB Important**

1. Read these directions for use carefully in order to obtain accurate measuring results.
2. This blood pressure monitor is only suitable for measuring the blood pressure of adults.
3. This blood pressure monitor measures blood pressure on the basis of pulse rate. If you have a very weak or irregular pulse, the blood pressure monitor may not be able to determine the right pressure value.
4. Eating, smoking, drinking (especially alcoholic beverages) and the use of medicines will influence the measurement.
5. If in doubt, consult a doctor.
6. Only a doctor is capable of interpreting your blood pressure values. Never change your medication without consulting a doctor.
7. Your blood pressure may fluctuate over time and during day. It may even vary more than 10mmHg within minutes. This means that successive measurements may show different blood pressure values. (See Shape 1)


Shape 1: Alteration diagram of blood pressure towards 12 hours of a 35 years old

**GB Conditions for accurate blood pressure measurement**

1. Sit down and relax for 10 minutes before you use the blood pressure monitor. Breath deeply in and out five times before you start measuring.
2. Wait at least 3 minutes between two successive measurements to avoid incorrect readings due to the temporary congestion of the arteries caused by the cuff.

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3. Perform the measurement every day at the same time, preferably early in the morning on empty stomach, since your blood pressure fluctuates in the course of the day.
4. Do not eat, smoke or drink before measuring your blood pressure. Alcoholic beverages in particular will influence the measurement.
5. Do not talk, laugh, cough or move your arm during the measurement, as this may influence the measurement results.
6. Do not cross your legs during measurement.
7. Do not attach too much importance to the results of a single measurement. Performing measurements regularly over a longer period of time and recording the results will give you a better insight into your blood pressure than a single measurement.

#### **GB** General information about blood pressure

Blood pressure is measured in mmHg (mm mercury pressure). Blood pressure consists of two values. Systolic pressure (SYS) is created by the contractions of the heart muscle, which pumps the blood into the arteries. Diastolic pressure (DIA) is created when the heart relaxes and fills with blood again. The World Health Organization (WHO) has laid down the following standard values with regards to the blood pressure of adults.

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## GB Introduction

Blood pressure measurements determined with KS 551 are equivalent to those obtained by a trained observer using cuff/stethoscope auscultation method, within the limits prescribed by the American National Standard, Electronic or Automated Sphygmomanometers. This unit is to be used by adult consumers in a home environment. Do not use this device on infants or neonates. KS 551 is protected against manufacturing defects by an established International Warranty Program. For warranty information, you can contact the manufacturer, Rossmax International Ltd. or your local distributors.

**Attention:** Consult the accompanying documents.

Please read this manual carefully before use. For specific information on your own blood pressure, contact your physician. Please be sure to keep this manual.

## GB Preliminary Remarks

This Blood Pressure Monitor complies with the European regulations and bears the CE mark "CE 0120". The quality of the device has been verified and conforms to the provisions of the EC council directive 93/42/EEC dated 14 June 1993 concerning medical devices:

EN 1060-1: 1995/A1: 2002 Non-invasive sphygmomanometers - Part 1 General requirements, EN 1060-3: 1997/A1: 2005 Non-invasive sphygmomanometers - Part 3 Supplementary requirements for electromechanical blood pressure measuring systems, and EN 1060-4: 2004 Non-invasive sphygmomanometers - Part 4: Test Procedures to determine the overall system accuracy of automated non-invasive sphygmomanometers.

This blood pressure monitor was designed for long service time. To ensure accurate measurements, it is recommended that calibration be completed every two years.

## GB How " KS 551 " Works

The **KS** Series uses oscillometric method to detect your blood pressure. Before the cuff starts inflating, the device will establish a baseline cuff pressure equivalent to the air pressure. Your **KS** unit will determine the appropriate inflation level based on one's pressure oscillations, followed by cuff deflation. During the deflation, the device will detect the amplitude and slope of the pressure oscillations and thereby determine the systolic pressure, diastolic pressure, and pulse for you.

Systolic (SYS.)	< 120	120 – 139	140 – 159	≥ 160
Diastolic (DIA.)	< 80	80 – 89	90 – 99	≥ 100
<b>GB</b>	Normal	Prehypertension	Stage1 Hypertension	Stage2 Hypertension
<b>D</b>	Normal	Prähypertonie	Hypertension Bluthochdruck, Stufe 1	Hypertension Bluthochdruck, Stufe 2
<b>F</b>	Normale	Préhypertension	Hypertension stade 1	Hypertension stade 2
<b>E</b>	Normal	Prehipertensió	Fase 1 Hipertensió	Fase2 Hipertensió
<b>P</b>	Normal	Pré-hipertensão	Hipertensão de estágio 1	Hipertensão de estágio 2
<b>PL</b>	Ciśnienie normalne	Ciśnienie podwyższone	Nadciśnienie pierwszego stopnia	Nadciśnienie drugiego stopnia
<b>I</b>	Normal	Pre-iptensione	Stadio1 Iptensione	Stadio2 Iptensione
<b>NL</b>	Normaal	Prehypertensie	Fase 1 Hypertensie	Fase 2 Hypertensie
<b>TR</b>	Normal	Prehipertansiyon	1. Evre Hipertansiyon	2. Evre Hipertansiyon
<b>AR</b>	المرحلة الطبيعية	ما قبل زيادة ضغط الدم	المرحلة الثانية من زيادة ضغط الدم	المرحلة الأولى من زيادة ضغط الدم
<b>R</b>	Норма	Предгипертоническое состояние	1стадия гипертонии	2стадия гипертонии

**GB** For reliable monitoring and reference of blood pressure, keeping long- term records is recommended.

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### **GB** About Blood Pressure

The National High Blood Pressure Education Program Coordinating Committee has developed a blood pressure standard, according to which areas of low- and high-risk blood pressure are identified. This standard, however, is a general guideline as individual's blood pressure varies among different people and different age groups...etc.

It is important that you consult with your physician regularly. Your physician will tell you your normal blood pressure range as well as the point at which you will be considered at risk.

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### **GB** Hypertension Classification Indicator

The National High Blood Pressure Education Program Coordinating Committee has developed a blood pressure standard according to which areas of low- and high-risk blood pressure are identified, classifying blood pressure ranges into 4 stages. This unit is equipped with innovative blood pressure classification indicator, which visually indicates the assumed risk level (prehypertension / stage 1 hypertension / stage 2 hypertension) of the result after each measurement. (Please note above.)

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
### GB Applying the cuff

1. Plug the connector into the unit.
2. Put the "Index" end of the cuff through the D-ring, your arm through the cuff loop, and wrap the cuff with tube pointing your palm, the artery mark over your main artery.
3. The edge of the cuff should be at approximately 2 to 3 cm above the inner side of the elbow joint. If the Index line falls within the OK range, the cuff circumference is suitable, otherwise you may need cuff with other circumferences.

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## GB Error Messages

Error Message	Root Causes	Corrections
EE	Measurement Error	Exceeding measurement range. Make sure the L-plug is securely connected to the air socket and measure again quietly. Wrap the cuff correctly and keep arm steady during measurement. If error keeps occurring, return the device to the local distributor or service center.
E1	Air Circuit Abnormality	Make sure the L-plug is securely connected to the air socket and measure again quietly. If error keeps occurring, return the device to the local distributor or service center.
E2	Pressure Exceeding 300 mmHg	Switch the unit off and measure again quietly. If error keeps occurring, return the device to the local distributor or service center.
E3	Data Error	Remove and reload the batteries. If error keeps occurring, return the device to the local distributor or service center.
Er	Exceeding Measurement Range	Measure again quietly. If error keeps occurring, return the device to the local distributor or service center.
CA	Recalibration	The device has completed a measurement number after which it is recommended to proceed on recalibration in order to secure device correct function. Please contact local Kessler service to arrange recalibration.

## GB Cautionary Notes

- Wait about 5 minutes before taking the next measurement.
- The unit contains high-precision assemblies. Therefore, avoid extreme temperatures, humidity, and direct sunlight. Avoid dropping or strongly shocking the main unit, and protect it from dust.
- Leaky batteries can damage the unit. Remove the batteries when the unit is not used for a long time.
- The unit should not be operated by children so to avoid hazardous situations.
- The unit is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please your Kessler local service.
- For users diagnosed with common arrhythmia (atrial or ventricular premature beats or atrial fibrillation), diabetes, poor circulation of blood, kidney problems, or for users suffered from stroke, or for unconscious users, the device may have difficulty in determining the proper blood pressure.
- Clean the blood pressure monitor body and the cuff carefully with a slightly damp, soft cloth. Do not press. Do not wash the cuff or use chemical cleaner on it. Never use thinner, alcohol or petrol (gasoline) as cleaner.
- To stop inflation at any time, press the ON/OFF/START key, and the air in the cuff will be rapidly exhausted.
- Once the inflation reaches 300 mmHg, the unit will start deflating rapidly for safety reason.
- Please note that this is a home healthcare product only and it is not intended to serve as a substitute for the advice of a physician or medical professional.
- Do not use this device for diagnosis or treatment of any health problem or disease. Measurement results are for reference only. Consult a healthcare professional for interpretation of pressure measurements. Contact your physician if you have or suspect any medical problem. Do not change your medications without the advice of your physician or healthcare professional.
- Electromagnetic interference: The device contains sensitive electronic components. Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave ovens). These may lead to temporary impairment of measurement accuracy.
- Dispose of device, batteries, components and accessories should according to local regulations.
- This monitor may not meet its performance specification if stored or used outside temperature and humidity ranges specified in Specifications.

### Notes:

Batteries are hazardous waste.

Do not dispose them together with the household garbage.

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**GB Specifications**

Measurement Method	: Oscillometric
Measurement Range	: Pressure: 40~250mmHg; Pulse: 40~199 beats/ minute
Pressure Sensor	: Semi conductor
Accuracy	: Pressure: $\pm 3$ mmHg; Pulse: $\pm 5\%$ of reading
Inflation	: Pump Driven
Deflation	: Automatic Air Release Valve
Memory capacity	: 99 memories
Auto-shut-off	: 1 minute after last button operation
Operation	: Temperature: 10°C ~40°C
Environment	(50°F ~104°F); Humidity: 40% ~ 85% RH max
Storage Environment	: Temperature: -10°C ~ 60°C ( 14°F ~140°F ); Humidity: 10% ~ 95% RH max
Power Source	: 6V DC four R06 (AA) Batteries
AC Power Source	: 12V DC, 800mA,(Plug size:Outer(-) isΦ3.8, inner(+))isΦ1.3)
Dimensions	: 140 (L) X 137 (W) X 68 (H) mm
Weight	: 410g (G.W.) (w/o Batteries)
Arm circumference	: Adult: 24 ~36cm (9.4"~14.2")
<input type="checkbox"/>	: Type BF Device and cuff are designed to provide special protection against electrical shocks.
IPX0	: Ordinary equipment (enclosed equipment without protection against ingress of water)
NO AP/APG	: Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.

\*Specifications are subject to change without notice.

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