## CS 410

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# AUTOMATIC ARM SPHYGMOMANOMETER

## Features:

- Fuzzy logic technology
- 90 memories
- One Touch easy measurement
- Easy-to-read display

## Introduction


Blood pressure measurements taken with the CS 410 are equivalent to those obtained by a trained person using the armband/stethoscope auscultation method, within the limits prescribed by the American National Standard for Electronic or Automated Sphygmomanometers. The CS 410 is designed for home use by adults.

## Caution:


Please read this manual carefully before use. For specific information on your own blood pressure, please contact your doctor. Keep this manual in a safe place for future reference.

## Name/Function of each part:

4 x "AA" size  
1.5 V batteries

**AC adapter** producing: 6V DC 1°  
(5.5Ø/2.1 Ø) Mod. A30610G

## Fuzzy Logic measuring technology

The CS 410 uses the oscillometric method to measure arterial blood pressure. Before the arm-band is inflated, the unit will establish a baseline armband pressure equivalent to the air pressure. At the press of a button, the armband is automatically inflated. As it inflates, the CS 410 will measure your arterial blood pressure, analysing its oscillations. Once the amplitude and gradient of the pressure oscillations have been determined, the unit will display the systolic and diastolic pressure and pulse rate values.

## Preliminary remarks

This sphygmomanometer complies with the European regulations and bears the EC mark "CE 0068". The quality of the unit has been verified and conforms with the provisions of European Council directive 93/42/EEC dated 14th June 1993 regarding medical equipment, as well as EMC directive 89/336/EEC:

|                       |                                                                                                                                    |
|-----------------------|------------------------------------------------------------------------------------------------------------------------------------|
| EN 1060-1 (Nov. 1995) | Non-invasive blood pressure measuring equipment. General requirements                                                              |
| EN 1060-3 (Feb. 1997) | Non-invasive blood pressure measuring equipment. Supplementary requirements for electromechanical blood pressure measuring systems |

To ensure accurate measurements, the calibration of the unit should be checked every two years.

## Arterial blood pressure

*What is arterial blood pressure?*

Arterial blood pressure is the pressure exerted on the artery as blood flows through it. The pressure measured when the heart contracts and pumps out blood is the systolic pressure (highest). The pressure measured when the heart dilates and blood flows back is called diastolic (lowest) blood pressure.


*Why measure your blood pressure?*

Amongst the various health problems afflicting modern man, problems linked to high blood pressure are by far the most common. The dangerously close relationship between high blood pressure and heart disease, characterised by high morbidity, has made it necessary to monitor blood pressure in order to identify those at risk.

## Blood pressure standard

The World Health Organization (WHO) has established a blood pressure standard, identifying risk areas for hypotension (low blood pressure) and hypertension (high blood pressure). This standard, however, is intended as a general guideline because individual blood pressure measurements vary according to population, age group, etc.


It is important to consult your doctor regularly. Your doctor will be able to tell you your normal blood pressure range as well as the point at which you should consider yourself at risk.

## Blood Pressure Fluctuation

*Blood pressure fluctuates all the time!*

Do not worry if two or three measurements are high. Blood pressure fluctuates during the month and even throughout the day. It is also influenced by the season and temperature

### Key to display Display:

Systolic pressure

Diastolic pressure

Pulse rate

— N° of memory sequence

### Symbols

**mem.**

Showing when the measurement value is memorised or recalled from the memory.

Showing when the batteries need replacing


Showing the pulse rate per minute

Showing if a mistake is made during measurement

## Dual power supply


### Inserting batteries

1. Insert or replace 4 “AA” size batteries in the battery compartment, according to the diagram shown on the inside.
2. Replace the cover of the battery compartment by placing the end tabs of the battery cover in the slots on the top of the battery compartment and pushing down the other end to secure the battery cover.

### Replace the batteries when:

1. The flat battery symbol appears on the display
2. Nothing appears on the display when the ON/OFF/START button is pressed

### AC adapter


Mod. A30610G

This appliance can also be powered by an AC adapter. The AC adapter should be connected as shown in the diagram (left).



Note: Batteries are hazardous waste. Do not dispose of them together with domestic refuse.

## Applying the armband


1. Connect the armband to the unit.
2. Open the armband, leaving the “Index” end in ring D.

3. Apply armband to left arm. The “OK” interval mark should be turned towards you, whilst the rubber tube should be positioned in the same direction as your arm. Place the artery mark  $\Phi$  on the main arteries (on the inside of the arm) on the upper part of the arm. Tighten the armband by pulling the “Index” end.

- Turn the palm of your left hand upwards and position the edge of the armband about 2-3 cm above the inside of your elbow.

- If the Index line falls within the OK interval, shown on the upper edge of the armband, then this armband is suitable for you. If, on the other hand, the "index" line is situated outside the OK interval, you may require a special armband.

## Taking a measurement

- Press the ON/OFF/START button. The entire display will be shown for about a second before returning to "0".

- The unit will automatically start inflating the armband to measure your blood pressure. During inflation the display will show a variable reading of the pressure inside the armband.

**It is very important to remain still and quiet during measurement. The slightest movement could result in an inaccurate measurement.**

- When the measurement is complete, the systolic and diastolic pressure and pulse rate will be shown simultaneously and automatically saved. The memory can store up to 90 measurements.
- Press the ON/OFF button to switch the unit off. If no button is pressed, the unit will turn itself off automatically after 1 minute.

## Recalling values from the memory

- Press the Memory button to display the average values of the last three measurements.
- Press the Memory button a second time to recall the previous series of measurements.
- Each measurement stored in the memory is displayed with the relative sequence number.

## Helpful tips

Here are a few helpful tips to help you obtain more accurate readings:

- Do not measure your blood pressure immediately after eating a large meal. To obtain more accurate readings, wait at least one hour before measuring your blood pressure.
- Do not smoke or consume alcohol before measuring your blood pressure.
- Do not measure your blood pressure when you are physically tired.
- It is important to relax during measurement. Try to rest for 15 minutes before measuring your blood pressure.
- Do not measure your blood pressure if you are tense or under stress.
- Measure your blood pressure at normal body temperature. If you are feeling hot or cold, wait a while before taking a measurement.
- If the monitor is stored at very low temperature (near freezing), put it in a warm place for at least an hour before use.
- Leave at least 5 minutes between each measurement.

## Troubleshooting

If any faults should occur during use, please check the following points:

| Symptom                                                                                     | Check point                                              | Correction                                      |
|---------------------------------------------------------------------------------------------|----------------------------------------------------------|-------------------------------------------------|
| Nothing appears on the display when the ON/OFF/START button is pressed                      | Have the batteries run down?                             | Replace them with 4 new alkaline batteries.     |
|                                                                                             | Are the batteries inserted correctly?                    | Reinsert the batteries in the correct position. |
| The EE symbol appears on the display or the blood pressure value is excessively high or low | Is the armband correctly applied in the proper position? | Take the measurement again without moving.      |
|                                                                                             | Did you talk or move during measurement?                 | Take the measurement again without moving.      |
|                                                                                             | Did you move your wrist whilst wearing the armband?      |                                                 |

**Note:** If the unit still does not work, take it to an authorised Artsana service centre. Under no circumstances should you disassemble and attempt to repair the unit yourself.

## Warnings:

1. The appliance contains high-precision components. Therefore, avoid extreme temperatures, humidity and direct sunlight. Avoid dropping or vigorously shaking the main unit and protect it from dust.
2. Clean the monitor body and armband carefully with a soft and slightly damp cloth. Do not press. Do not wash the armband or use chemical cleaners on it. Never use thinner, alcohol or petrol as cleaning fluids
3. Batteries that leak can damage the unit. Remove the batteries if the unit is not used for long periods.
4. For safety reasons do not let children use the device.
5. If the unit is stored at temperatures near freezing point, allow it to acclimatize at room temperature before use.
6. The CS 410 is not field serviceable. Do not attempt to open the device using tools or try to adjust anything inside it. If you have any problems, please contact an authorised Artsana service centre.

7. The unit is not suitable for patients with common arrhythmia (atrial or ventricular premature beats or atrial fibrillation), circulatory disorders or kidney problems, those who have suffered a stroke or are in a state of unconsciousness.
8. Do not use the appliance in the vicinity of devices which use radio frequencies, such as mobile telephones.

## Specifications\*:

|                        |                                                                                               |
|------------------------|-----------------------------------------------------------------------------------------------|
| Measurement method:    | oscillometric                                                                                 |
| Measurement range:     | pressure: 20-280 mmHg; pulse 40-180 beats per minute                                          |
| Pressure sensor:       | Semiconductor                                                                                 |
| Accuracy:              | blood pressure: $\pm 3$ mmHg; pulse: $\pm 5$ % of reading                                     |
| Inflation:             | pump driven                                                                                   |
| Deflation:             | automatic pressure release valve                                                              |
| Memory capacity:       | 90                                                                                            |
| Auto shut-off:         | 1 minute after last button was pressed                                                        |
| Operating environment: | temperature: 10 - 40 °C (50 - 104 °F); relative humidity: 40 - 85 %                           |
| Storage environment:   | temperature: -5 °C - 60 °C (23°F - 140°F); relative humidity: 10 - 95%                        |
| DC power source:       | 4 x LR6 (AA) batteries 1.5V DC 6V 1A DC                                                       |
| AC power source:       | MOD - A30610G<br>230V 50 Hz 0.2 A<br>(plug dimensions: internal (-) Ø 5.5 external (+) Ø 2.1) |
| Dimensions:            | 137(L) x 140.5(W) x 63.5(H) mm                                                                |
| Weight:                | 430 g (gross weight) (without batteries)                                                      |
| Arm circumference:     | adults: 22-32 cm (8.7" - 12.6") (standard)<br>large: 32-42 cm (12.6" - 16.5") (optional)      |
| Max. usage time:       | max. 2 min.                                                                                   |
|                        | BF type classification                                                                        |

\*Specifications are subject to change without prior notice

## Blood Pressure log:

| Name:          | Age: |  |  |  |  | Weight: |  |  |  |  |  |  |  |  |  |
|----------------|------|--|--|--|--|---------|--|--|--|--|--|--|--|--|--|
| Date           |      |  |  |  |  |         |  |  |  |  |  |  |  |  |  |
| Time           |      |  |  |  |  |         |  |  |  |  |  |  |  |  |  |
| mmHg           |      |  |  |  |  |         |  |  |  |  |  |  |  |  |  |
| 220            |      |  |  |  |  |         |  |  |  |  |  |  |  |  |  |
| 200            |      |  |  |  |  |         |  |  |  |  |  |  |  |  |  |
| 180            |      |  |  |  |  |         |  |  |  |  |  |  |  |  |  |
| 160            |      |  |  |  |  |         |  |  |  |  |  |  |  |  |  |
| 140            |      |  |  |  |  |         |  |  |  |  |  |  |  |  |  |
| 120            |      |  |  |  |  |         |  |  |  |  |  |  |  |  |  |
| 100            |      |  |  |  |  |         |  |  |  |  |  |  |  |  |  |
| 80             |      |  |  |  |  |         |  |  |  |  |  |  |  |  |  |
| 60             |      |  |  |  |  |         |  |  |  |  |  |  |  |  |  |
| Pulse          |      |  |  |  |  |         |  |  |  |  |  |  |  |  |  |
| Body Condition |      |  |  |  |  |         |  |  |  |  |  |  |  |  |  |

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