Model: **pM-B51**

Blood Pressure Monitor - Automatic Wrist
Otomatik Tansiyon Aleti - Otomatik, Bilek Tipi

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EN Instruction Manual

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EN Medical Disclaimer

This manual and product are not meant as a substitute for advice provided by your doctor.

You are not to use the information contained herein, or this product for diagnosing or treating a health problem or prescribing any medication. If you have or suspect that you have a medical problem, promptly consult your healthcare provider.

EN Intended Use

This device uses the oscillometric method to automatically measure systolic and diastolic blood pressure as well as heart rate.

The measurement position is at human being's wrist.

All values can be read out in one LCD panel.

The device is designed for home use and recommended for use by adults aged 18 years and older with wrist circumference ranging 135 ~ 195 mm (approx. 5.3 ~ 7.7 inch).

EN About Blood Pressure

What is blood pressure?

Blood pressure is the measurement of the force of blood pushing against the walls of the arteries. Arterial blood pressure is constantly fluctuating during the course of the cardiac cycle. The highest pressure in the cycle is called the systolic blood pressure, and represents the pressure in the artery when the heart is beating. The lowest pressure is the diastolic blood pressure, and represents the pressure in the artery when the heart is at rest. Both the systolic and the diastolic pressure are necessary for a physician to evaluate the status of a patient's blood pressure.

Many factors such as physical activity, anxiety or the time of day, can influence your blood pressure. Blood pressure is typically low in the mornings and increases from the afternoon to the evening. It is on average lower in the summer and higher in the winter.


Why is it useful to measure blood pressure at home?

Having one's blood pressure measured by a doctor in a hospital or a clinic, is often associated with a phenomenon called "White Coat Hypertension" where the patient becomes nervous or anxious, thus raising his blood pressure. There are also numerous other factors that might cause your blood pressure to be raised at a specific time of day. This is why medical practitioners recommend home monitoring as it is important to get readings of blood pressure during different times of the day to really get an idea of your real blood pressure.

Medical practitioners generally recommend the "Rule of 3", where you are encouraged to take your blood pressure three times in a row (at 3 ~ 5 minute interval), three times a day for three days. After three days you can average all the results and this will give you an accurate idea of what your blood pressure really is.


WHO blood pressure classifications

Standards for assessment of high or low blood pressure without regard to age, have been established by the World Health Organization (WHO), as shown in the chart.

However this chart is not exact for classification of blood pressure and it's intended to be used as a guide in understanding non-invasive blood pressure measurements. Please consult with your physician for proper diagnosis.

Variations in blood pressure

Individual blood pressures vary greatly both on a daily and a seasonal basis. These variations are even more pronounced in hyper tense patients. Normally the blood pressure rises while at work and is at its lowest during sleeping period.

(Hyper tense: means a person who has high blood pressure symptom.)

The graph above illustrated the variations in blood pressure over a whole day with measurement taken every five minutes. The thick line represents sleep. The rise in blood pressure at 4 PM (A in the graph) and 12 PM (B in the graph) correspond to an attack of pain.

EN Precautions

- Do not use this manual and product as a substitute for advice, diagnosing or treating a health problem or prescribing any medication by your doctor. If you have a medical problem, promptly consult your healthcare provider.
- Read the Instruction Manual thoroughly before measuring and keep it at hand for your reference at any time.
- This device uses the oscillometric method to measure systolic and diastolic blood pressure as well as your heart rate. It's recommended for use by people over the age of 18 and not to be used on infant or children.
- The device is designed for home use and not suitable for clinical use.
- Do not take a measurement in a low (less than 41 °F/5 °C) and high (more than 104 °F/40 °C) temperature, nor in a place outside humidity ranges (15% ~ 93% R.H.), or you may get inaccurate readings.
- Wait 30 ~ 45 minutes before measurement if you've just consumed caffeinated beverages or smoked cigarettes.
- Rest at least 5 ~ 10 minutes before taking a measurement.
- To allow your blood vessels to return to the condition prior to taking the measurement, please wait at least 3 ~ 5 minutes in between measurements. You may need to adjust the wait time according to your personal physiological situation.
- We recommend you using the same wrist (preferably the left wrist) and measuring around the same time each day.

- Sit down comfortably and place your elbow on the table with your feet flat on the floor. Please do not cross your legs during measurements.
- Keep the device at heart level. Relax your hand with the palm facing up.
- Perform measurements in a quiet and relaxed environment at room temperature.
- Do not move or shake the device during a measurement. Please keep quiet and do not talk during measurements.
- Keep in mind that blood pressure naturally varies from time to time through out the day and is affected by lots of different factors such as stress, eating, smoking, alcohol consumption, medication, and physical activity, etc.
- Normally the blood pressure rises while at work and is at its lowest during sleeping period.
- Blood pressure measurements should be interpreted by a physician or a trained health professional who is familiar with your medical history. Using the unit and recording the results regularly for your physician to interpret, you will keep your physician informed of the continuing changes in your blood pressure.
- If you have one of the circulatory problems as arteriosclerosis, diabetes, liver disease, kidney disease, severe hypertension, peripheral circulation, please consult your healthcare professional before using the device.
- This product is not suitable for people with arrhythmias and pregnant women.
- Blood pressure measurements taken with this device are equivalent to those obtained by a trained observer using the cuff / stethoscope auscultation method and are within the accuracy limits prescribed by the Standard of EN 1060-4.

Attention!

1. Do not use the device on infants, children, or those who cannot express their own intention.
2. The device is equipped with sensitive electronic components. While measuring, avoid strong electrical or electromagnetic fields, e.g. mobile phones, microwave ovens, etc; or it may lead to temporary reading error or inaccuracy.
3. Over high frequency measurements may result in blood flow interference, which is likely to cause uncomfortable sensations, such as partial subcutaneous hemorrhage, or temporary numbness to your wrist. In general, these symptoms should not last long. However, if you do not recover in time, please seek your medical practitioners for help.

EN Device Overview

Product components

Caution!

Substitution of a component different from that supplied might result in measurement error.

Unit display

1. Hour : Minute (Also represents YEAR in Setting Mode)
2. Month / Date
3. Blood Pressure Unit
4. Error Symbol
5. Memory Record Number
6. Low Battery Symbol
7. Pulse Symbol
8. Irregular Heartbeat Detector
9. Heart Rate
10. Diastolic Rate
11. Risk Category Indicator
12. Systolic Rate

EN Symbol Definitions

Symbols	Definitions
 Low Battery Symbol	This symbol appears when the battery power is excessively low or the polarity reverses. We suggest you replace all batteries with new ones, and make sure the +/- polarities are properly positioned.
 Pulse Symbol	Once pulse is detected, the symbol flashes with each pulse beat. Our suggestion: Please do not talk or move during measurements.
 Irregular Heartbeat Detector	This symbol appears for 1 minute when the user was talking, moving, shaking, or an irregular heart beat was detected during measurements. Our suggestion: Please do not talk or move during measurements. Repeat the measurement after resting for at least 5 minutes, and restart your measurement while sitting down comfortably and quietly.
 Risk Category Indicator Bar	The arrowhead points out the specific Risk Category that your measurement reading fits in.


Risk Category Indicator

This device is equipped with Risk Category Indicator which classifies your blood pressure measurements into six stages (Optimal to Severe hypertension) as shown in below chart:

Stages of Blood Pressure Levels		Systolic (mmHg)	Diastolic (mmHg)	Color	Recommendation
Grade 3	Severe Hypertension	≥ 180	≥ 110	Red	Confirm immediately and repeat BP in one day and again within one week depending on clinical situation.
Grade 2	Moderate Hypertension	160 ~ 179	100 ~ 109	Red	Serial blood pressures repeated within one month.
Grade 1	Mild Hypertension	140 ~ 159	90 ~ 99	Red	Provide advice about lifestyle modification and confirm within two months.
High-Normal		130 ~ 139	85 ~ 89	Orange	Provide advice about lifestyle modification and recheck in one year.
Normal		120 ~ 129	80 ~ 84	Yellow	Recheck in 2 - 5 years. (patients aged > 75 years offered annual health check)
Optimal		< 120	< 80	Green	

*Source: WHO, 2003

After each measurement is completed, LCD display will show your position automatically on the six segments of the bar indicator which corresponds to Risk Category Indicator.

Note!

When a person's systolic and diastolic pressures fall into different categories, the higher category should apply.

e.g. systolic rate 181 & diastolic rate 99 > Red category (Severe Hypertension)

e.g. systolic rate 110 & diastolic rate 95 > Red category (Mild Hypertension)

Note!

The above table is not exact for classification of blood pressure and it's intended to be used as a guide in understanding non-invasive blood pressure measurements. Usually this is not a cause for concern; however we recommend you consult with your physician for proper diagnosis or seek medical advice according to our recommendation mentioned above. Please note that the device does not appropriate to diagnose hypertension, and it is only for user reference on blood pressure monitoring.

Irregular Heartbeat Detector

The symbol  will appear on screen indicating a certain heartbeat irregularity was detected during measurement. The heartbeat rhythm that is more than or less than 25% from the average rhythm is usually defined as an irregular heartbeat rhythm. Talking, moving, shaking or an irregular pulse during the measurement can result in the appearance of this symbol. Usually this is not a cause for concern, however if the symbol appears often, we recommend you seek medical advice. And please note that the device does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.


Note!

- The pulse display is not suitable for checking the frequency of heart pacemakers. If a certain pulse irregularity is detected during measurement often, we recommend you seek medical advice
- As a safeguard, we recommend that if you have arrhythmias such as atrial or ventricular premature beats and atrial fibrillation or any other special conditions you should check with your physician before using your device.
- The IHB function is not designed for use by people with arrhythmias nor for diagnosing or treating an arrhythmic problem. In order to filter the unstable status of user and avoid affecting the detection of heart rate from any movement, shaking or talking in the beginning of measurement, the method of averaging heart beat intervals of subject device is calculated with the three proper heart beat pulses detected in the beginning of measurement and that is different from a strict mathematical averaging of all recorded intervals.
- At least 3 beats with at least 25 % difference from the average heart beat interval will generate the IHB icon on the screen.

EN Installing Batteries

When LOW BATTERY SYMBOL  appears on the display, or no reaction toward operation, please change batteries.

Replace all worn-out batteries with new ones and do not mix new and used batteries. Do not mix alkaline, standard (carbon-zinc) or rechargeable (cadmium) batteries either. Such action may shorten the battery life or cause the device to malfunction. Slide the battery cover and insert 2 AAA alkaline batteries into the battery compartment as shown on the figure below. Make sure the polarities "+" and "-" ends are properly positioned.

Attention!

- Batteries are hazardous waste. Do not dispose of them together with the household garbage. Please discard worn-out batteries to the recycling site according to local regulations.
- Keep the battery away from children in case they choke on it.
- If the device is not to be used for over 2 months, please remove the batteries from its compartment for power-saving.
- After replacing the batteries, reset date and time.

EN Applying the Cuff

Do not place the pressure cuff over a jacket or sweater sleeve. Wrap the pressure cuff around the bare wrist with the monitor facing you.


Wrap the cuff snugly. Do not make it too tight.

Fold the remaining part of the cuff back out of the way.

Leave approximately 0.4 inch (10 mm) between the cuff and the bottom of your hand palm.

0.4 inch (10 mm)

5.3 ~ 7.7 inch (135 ~ 195 mm)

Note!

- Do not use this device if your wrist has any wound or injury.
- Do not wrap the cuff around any body part other than your wrist.

EN Positioning Guide

It is extremely important that the cuff be at the same height as the heart. Having the cuff higher or lower may cause inaccurate results.


1. Sit down comfortably with your feet flat on the floor.
2. Position the blood pressure monitor on your wrist.
3. Place your elbow on the table and rest the back of your hand on the device storage case or other object
4. Rest your wrist on the armrest until it's at the same height as your heart.
5. Relax your hand and turn your palm upwards.

EN Measurement Procedure

Switch on the monitor

- Press **⏻** button to switch on the monitor. (**⏻** button displays blue backlight for 10 seconds.)
- All segments appear on the screen.

Setting year, date and time

- Press **⌚** button („YEAR“ flashes). Press **▲** or **▼** button to adjust YEAR value.
- Press **⌚** button („MONTH“ flashes). Use **▲** or **▼** button to adjust MONTH (1, 2, 3,....., 12).
- Adjust DATE (1, 2, 3,..., 31), HOUR (1, 2, 3,.....12PM, 1 PM,..., 12) and MINUTE (00,01,02,03,.....59) as described in Step A above. When settings are done, press **⌚** button to confirm the entries. The device is ready to measure.

Taking a Measurement

A. Before measurement, press **▲** or **▼** button to select User 1, 2, or 3.

- B. With the cuff wrapped around your wrist, press **⏻** button to start measurement. All display units appear on the screen.

Note!

Do not inflate the cuff until it is wrapped around your wrist.

- After all symbols disappear, the display will show "00". The monitor is "Ready to Measure" and will automatically inflate to the level that is right for you.

- C. After inflation of the cuff, the pressure will slowly decrease. When pulse is detected, PULSE SYMBOL ♥ flashes.

**Note!**

If the cuff does not stop inflating, remove the cuff at once. To stop measurement, press ⓘ button.

- D. LCD screen displays your systolic rate, diastolic rate, pulse, Risk Category Indicator Bar, and Irregular Heartbeat Detector symbol (if any) with date and time for 1 minute.

- E. Without any operation for 1 minute, device automatically shuts off. Please do not try to close it manually.

EN Memory Function

Storing data

After each measurement, the systolic and diastolic pressure, heart rate with date and time will be automatically stored. The monitor can store up to 120 memory sets for 3 users, and automatically replace the oldest data with new one.


Recalling data

- Press **A** or **V** button to select User 1, 2, or 3.
- Press M button to enter Memory Mode. LCD displays average of last 3 measuring results first.
- Press M button again, LCD displays the latest measuring result. Use **A** or **V** button to scroll through all stored measuring results.
- To stop reading memories, press **I** button, and switch to Standby Mode.

Erasing data

- Press **A** or **V** button to select User 1, 2, or 3.
- Press M button to enter Memory Mode.
- Press and hold **A** and **I** buttons at the same time, the data will be erased automatically.
- To confirm the data in the selected user has been erased, press M button and no data should appear.

Note!

Once deleted, your data can NOT be restored.

EN Storage and Maintenance

General Use

- Do not in any way twist the cuff.
- Do not press  button if the cuff is not wrapped around the wrist.
- Do not drop the product and avoid any strong impacts.

Maintenance

- Use a piece of cloth with water or mild cleansing agent to wipe the device and dry it immediately with a dry cloth.
- Do not use detergent or any strong chemicals to clean the device.
- Use only a dry cloth to wipe the cuff.
- Do not attempt to disassemble or change any parts of the monitor, including wrist cuff, due to substitution of a component different from that supplied might result in measurement error.
- If any suggestion or service is requested, please consult your service station.

Storage

- If the device is not to be used for a long time, please remove the batteries from the device (leaking of battery acid can cause the device to malfunction).
- Always store the unit in the storage case after use.
- Do not place the device directly under sunlight, in high temperature, or in humid or dusty places.
- Do not store the device in extremely low (less than $-13^{\circ}\text{F}/-25^{\circ}\text{C}$) and high (more than $158^{\circ}\text{F}/70^{\circ}\text{C}$) temperature, nor in a place its humidity exceeds 93% R.H.

EN Troubleshooting

SYMBOLS & SYMPTOMS	CONDITIONS & CAUSES	INDICATION & CORRECTION
Unit does not turn on when ① button is pushed.	Worn-out batteries. Battery polarities have been positioned incorrectly.	Replace them with 2 new AAA alkaline batteries. Re-insert the batteries in the correct positions.
 Measuring Error Symbol appears when blood pressure value displayed is excessively low or high.	Cuff has been placed incorrectly. Did you talk or move during measurement? Shaking of the arm with the cuff on.	Wrap the cuff properly so that it is positioned correctly. Measure again. Keep arm steady during measurement. Measure again. Keep arm steady during measurement.
 Measuring Error Symbol	Air circuit abnormality. Cuff tube may not be plugged into monitor correctly.	Check cuff connection. Measure again.
 Measuring Error Symbol	Inflation pressure exceeding 300 mmHg.	Switch the unit off, then measure again.
 Measuring Error Symbol	Error determining measurement data.	Measure again.
Note! If „EP“ appears on the display, just return the device to your local distributor.		

EN Limited Warranty

To ensure continued measurement precision, all digital blood pressure monitors require recalibration regularly.

After 2 years from the manufacturing date, we recommend you have your monitor recalibrated at the local distributor or importer.

Please contact your distributor/importer for the details about the recalibration service and the charge of shipping and handling.

Please also note that this service does not cover damage caused by misuse or abuse; accident; the attachment of any unauthorized accessory; alteration to the product; improper installation; unauthorized repairs or modifications; improper use of electrical/power supply; loss of power; dropped product; malfunction or damage of an operating part from failure to provide manufacturer's recommended maintenance; transportation damage; theft; neglect; vandalism; or environmental conditions; loss of use during the period the product is at a repair facility or otherwise awaiting parts or repair; or any other conditions whatsoever that are beyond the control of importers or distributors.

EN Specifications

Model Number	pM-B51
Measurement Method	Oscillometric
Measurement Range	Pressure: 0 ~ 300 mmHg Pulse: 40 ~ 199 Beats / Minute
Accuracy	Pressure: ± 3 mmHg Pulse: ± 5 % Max.
Inflation	Automatic Inflation (Air Pump)
Deflation	Automatic Air Release Control Valve
Display	Liquid Crystal Display
Memory	120 Memory Total for 3 Users
Unit Dimensions	2.85 x 78 x 32.8 mm (L X W X H) 2.87 x 3.07 x 1.29 inch (L X W X H)
Unit Weight (Cuff & Batteries Excluded)	142.5 g $\pm$ 5 g (5.03 oz $\pm$ 0.18 oz)
Cuff Size	135 ~ 195 mm (5.3 ~ 7.7 inch)
Storage/Transportation Environment	Temperature: -25 °C ~ 70 °C (-13 °F ~ 158 °F) Humidity: ≤ 93 % R.H.
Operation Environment	Temperature: 5 °C ~ 40 °C (41 °F ~ 104 °F) Humidity: 15 % ~ 93 % R.H.
Power Supply	DC 3 V, AAA (LR03) (1.5V) Alkaline Battery x 2
Battery Life	Approx. 250 Measurements
Power-saving Mode	Without any operation for 1 minute, device automatically shuts off.
Accessories	Instruction Manual, 2 AAA Alkaline Batteries, Storage Box

**The contents of this manual and the specifications of the device covered by this manual are subject to change for improvement without notice.*

EN Note

This blood pressure monitor complies with the EC Directive (93/42/EEC) and bears the CE mark This blood pressure monitor also complies with mainly following standards (included but not limited):

Safety standard:

EN 60601-1 Medical electrical equipment part 1: General requirements for safety

EMC standard:

EN 60601-1-2 Medical electrical equipment part 1-2: General requirements for safety- Collateral standard: Electromagnetic compatibility- Requirements and tests

Performance standards:

EN 1060-1 Non-invasive sphygmomanometers - General requirements

EN 1060-3 Non-invasive sphygmomanometers - Supplementary requirements for electromechanical blood pressure measuring systems.

EN 1060-4 Non-invasive sphygmomanometers - Test procedures to determine the overall system accuracy of automated non-invasive sphygmomanometers.

Follow instructions for use.

BF Classification:

- Internally powered equipment
- BF type applied part
- IPX0
- Not suitable for use in presence of flammable anesthetic mixture with air or with Oxygen or nitrous oxide
- Continuous operation with short-time loading

Keep dry.

To avoid inaccurate results caused by electromagnetic interference between electrical and electronic equipments, do not use the device near a mobile phone or microwave oven. At least keep a maximum output power of 2 W yields and a distance 3.3m away from this equipment.

Discard the used product to the recycling collection point according to local regulations.

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