EN English Introduction

Blood pressure measurements determined with CG155F are equivalent to those obtained by a trained observer using cuff/stethoscope auscultation method, within the limits prescribed by the American National Standard, Electronic or Automated Sphygmomanometers. This unit is to be used by adult consumers in a home environment. Do not use this device on infants or neonates, CG155F is protected against manufacturing defects by an established International Warranty Program. For warranty information, you can contact your local distributors.

Attention: Consult the accompanying documents. Please read this manual carefully before use. For specific information on your own blood pressure, contact your physician. Please be sure to keep this manual.

Real Fuzzy Measuring Technology

This unit uses the oscilometric method to detect your blood pressure. Before the cuff starts inflating, the device will establish a baseline cuff pressure equivalent to the air pressure. This unit will determine the appropriate inflation level based on pressure oscillations, followed by cuff deflation.

During the deflation, the device will detect the amplitude and slope of the pressure oscillations and thereby determine for you the systolic blood pressure, diastolic blood pressure, and pulse.

Preliminary Remark

This Blood Pressure Monitor complies with the European regulations and bears the CE mark "CE 0120". The quality of the device has been verified and conforms to the provisions of the EC council directive 93/42/EEC (Medical Device Directive), Annex I essential requirements and applied harmonized standards. EN 1060-1: 1995/A2: 2009 Non-invasive sphygmomanometers - Part 1 - General requirements EN 1060-3: 1997/A2: 2009 Non-invasive sphygmomanometers - Part 3 - Supplementary requirements for electro-mechanical blood pressure measuring systems EN 1060-4: 2004 Non-invasive sphygmomanometers - Part 4: Test Procedures to determine the overall system accuracy of automated non-invasive sphygmomanometers.

This blood pressure monitor was designed for long service time. To ensure accurate measurements, this monitor is recommended to be re-calibrated every two years.

Blood Pressure Standard

The National High Blood Pressure Education Program Coordinating Committee has developed a blood pressure standard, classifying blood pressure ranges into 4 stages. (Ref. The Seventh Report of the Joint National Committee on Prevention, Detection, Evaluation, and Treatment of High Blood Pressure-Complete Report, JNC-7, 2003). This blood pressure classification are based on historical data, and may not be directly applicable to any particular patient. It is important that you consult with your physician regularly. Your physician will tell you your normal blood pressure range as well as the point at which you will be considered at risk. For reliable monitoring and reference of blood pressure, keeping long-term records is recommended. Please download the blood pressure log at our website www.accumed.ch.

Display Explanations

EE / Measurement Error: Make sure the L-Plug is securely connected to the air socket and measure again quietly. Wrap the cuff correctly and keep arm steady during measurement. If the error keeps occurring, return the device to your local distributor or service center.

E1 / Air Circuit Abnormally: Make sure the L-Plug is securely connected to the air socket on the side of the unit and measure again quietly. If the errors still occur, return the device to your local distributor or service center.

E2 / Pressure Exceeding 300 mmHg: Switch the unit off and measure again quietly. If the error keeps occurring, return the device to your local distributor or service center.

E3 / Data Error: Remove the batteries, wait for 60 seconds, and re-load. If the error keeps occurring, return the device to your local distributor or service center.

Er / Exceeding Measurement Range: Measure again quietly. If the error keeps occurring, return the device to your local distributor or service center.

Movement Detection

The "Movement Detection" helps reminding the user to remain still and is indicating any body movement during measurement. The specified icon appears once a "body movement" has been detected during and after each measurement.

Note: It's highly recommended that you measure again if the icon appears.

Hypertension Risk Indication

The National High Blood Pressure Education Program Coordinating Committee has developed a blood pressure standard, classifying blood pressure ranges into 4 stages. This unit is equipped with innovative blood pressure risk indication, which visually indicates the assumed risk level (normal / 1 hypertension / 2 stage 1 hypertension / 3 stage 2 hypertension) of the result after each measurement.

Irregular Heartbeat (IHB) Detection

This unit is equipped with an Irregular Heartbeat (IHB) Detection which allows those who have an irregular heartbeat to obtain accurate measurements alerting the user of the presence of an irregular heart beat during the measurement.

Note: It is strongly recommended that you consult your physician if the IHB icon appears often.

Installing Batteries

1. Press down and lift the battery cover in the direction of the arrow to open the battery compartment.

2. Install or replace 4 "AAA" sized batteries in the battery compartment according to the indications inside the compartment.

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1. Replace the battery cover by clicking in the bottom hooks first, then push in the top end of the battery cover.

2. Replace the batteries in pairs. Remove batteries when unit is not in use for extended periods of time.

3. The ON/OFF/START key is pressed and nothing appears on display. Caution:

1. Batteries are hazardous waste. Do not dispose them together with the household garbage.

2. There are no user serviceable parts inside. Batteries or damage from old batteries are not covered by warranty.

3. Use exclusively brand batteries. Always replace with new batteries together. Use batteries of the same brand and same type.

Applying the Cuff

1. Unwrap the arm cuff, leaving the end of the cuff through the D-ring of the cuff.

2. Put your left arm through the cuff loop. The color strip indication should be positioned closer to you with the tube pointing in the direction of your arm (Fig. ②). Turn your left palm upward and place the edge of the arm cuff at approximately 1.5 to 2.5 cm above the inner side of the elbow joint (Fig. ②). Tighten the cuff by pulling the end of the cuff.

3. Center the tube over the middle of the arm. Press the hook and loop material together securely. Allow room for 2 fingers to fit between the cuff and your arm. Position the artery mark (②) over the main artery (on the inside of your arm) (Fig. ②, ③). Note: Locate the main artery by pressing with 2 fingers approximately 2 cm above the bend of your elbow on the inside of your left arm. Identify where the pulse can be felt the strongest. This is your main artery.

4. Plug in the cuff connecting tube into the unit (Fig. ②).

5. Lay your arm on a table (palm upward) so the cuff is at the same height as your heart. Make sure the tube is not kinked (Fig. ③).

6. This cuff is suitable for your use if the arrow falls within the solid color line as shown on the right (Fig. ②). If the arrow falls outside the solid color line, you will need a cuff with other circumferences. Contact your local dealer for additional size cuffs.

Measurement Procedures

Here are a few helpful tips to help you obtain more accurate readings:

• Blood pressure changes with every heartbeat and is in constant fluctuation throughout the day.

• Blood pressure recording can be affected by the position of the user, his or her physiological condition and other factors. For greatest accuracy, wait one hour after exercising, bathing, eating, drinking beverages with alcohol or caffeine, or smoking to measure blood pressure.

• Before measurement, it's suggested that you sit quietly for at least 5 minutes as measurement taken during a relaxed state will have greater accuracy. You should not be physically tired or exhausted while taking a measurement.

• Do not take measurements if you are under stress or tension.

• During measurement, do not talk or move your arm or hand muscles.

• Take your blood pressure at normal body temperature. If you are feeling cold or hot, wait a while before taking a measurement.

• If the monitor is stored at very low temperature (near freezing), have it placed at a warm location for at least one hour before using it.

• Wait 5 minutes before taking the next measurement.

1. Press the ON/OFF/START key. All displays will appear for approximately one second before returning to "0".

2. The unit will automatically inflate to the appropriate inflation level based on the user's pulse oscillations. Measurement will then begin. It is important to remain still and quiet during measurement. Any significant movement may affect measurement results.

3. When the measurement is completed, systolic, diastolic and pulse will be shown simultaneously and be saved automatically in memory system. Up to 90 memories can be saved.

4. Measurement is now completed. Press the ON/OFF/START key to turn off the power. If no key is pressed, the unit will shut off automatically in 1 minute.

This monitor will re-inflate automatically to approximately 220 mmHg if the system detects that your body needs more pressure to measure your blood pressure.

Note: 1. This monitor automatically switches off approximately 1 minute after last key operation.

2. To interrupt the measurement, simply press the Memory or ON/OFF/START key; the cuff will deflate immediately.

3. During the measurement, do not talk or move your arm or hand muscles.

Recalling Values from Memory

1. To recall stored blood pressure readings from memory, simply press the Memory key, the first reading displayed is the average of the last 3 measurements stored in memory, and then the last set of memorized readings will be displayed.

2. Another press of the Memory key will recall the previous set of readings.

3. All readings stored in memory will be displayed with its sequence number. Every measurement made with an assigned memory sequence number.


Clearing Values from Memory

Press and hold the Memory key for approximately 5 seconds, then the data in the memory zone can be erased automatically.

Troubleshooting

If any abnormality will arise during use, please check the following points.

Symptoms	Check Points	Correction
No display	Have the batteries run down?	Replace them with four new batteries.
No display when the ON/OFF/START key is pressed	Have the batteries' polarities been positioned incorrectly?	Re-insert the batteries in the correct positions.

EE mark shown on display or the blood pressure value is displayed excessively low (high)	Is the cuff placed correctly? Did you talk or move during measurement? Did you vigorously shake the cuff during measurement?	Wrap the cuff properly so that it is positioned correctly. Measure again. Keep arm steady during measurement.
--	--	--

Note: If the unit still does not work, return it to your dealer. Under no circumstance should you disassemble and repair the unit by yourself.

Cautionary Notes

1. The unit contains high-precision assemblies. Therefore, avoid extreme temperatures, humidity, and direct sunlight. Avoid dropping or strongly shocking the main unit, and protect it from dust.

2. Clean the blood pressure monitor body and the cuff carefully with a slightly damp, soft cloth. Do not press. Do not wash the cuff or use chemical cleaner on it. Never use thinner, alcohol or petrol (gasoline) as cleaner.

3. Leaky batteries can damage the unit. Remove the batteries when the unit is not used for a long time.

4. The unit should not be operated by children so to avoid hazardous situations.

5. If the unit is stored near freezing, allow it to acclimate at room temperature before use.

6. This unit is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit.

7. As a common issue for all blood pressure monitors using the oscilometric measurement function, the device may have difficulty in determining the proper blood pressure for users diagnosed with common arrhythmia (atrial or ventricular premature beats or atrial fibrillation), diabetes, poor circulation of blood, kidney problems, or for users suffered from stroke, or for unconscious users.

8. To stop operation at any time, press the ON/OFF/START key, and the air in the cuff will be rapidly exhausted.

9. Once the inflation reaches 300 mmHg, the unit will start deflating rapidly for safety reasons.

10. Please note that this is a home healthcare product only and it is not intended to serve as a substitute for the advice of a physician or medical professional.

11. Do not use this device for diagnosis or treatment of any health problem or disease. Measurement results are for reference only. Consult a healthcare professional for interpretation of pressure measurement results. Contact your physician if you have or suspect any medical problem. Do not change your medications without the advice of your physician or healthcare professional.

12. Electromagnetic interference: The device contains sensitive electronic components. Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave ovens). These may lead to temporary impairment of measurement accuracy.

13. Dispose of device, batteries, components and accessories according to local regulations.

14. This monitor may not meet its performance specification if stored or used outside temperature and humidity ranges specified in Specifications.

Specifications

Measurement Method Oscilometric
Measurement Range Pressure: 30~260 mmHg; Pulse: 40~199 beats/minute

Pressure Sensor Semi-conductor
Accuracy Pressure: ± 3mmHg; Pulse: ± 5% of reading

Inflation Pump Driven
Deflation Automatic Air Release Valve

Memory capacity 90 memories
Auto-shut-off 1 minute after last key operation

Operation Environment 10°C~40°C (50°F~104°F); 40%~85% RH;
Storage and Transportation 10°C~60°C (14°F~140°F); 10%~90% RH;

DC Power Source 24~40 cm (9.4"~15.7")
Dimensions 124 (L) X 85 (W) X 68.6 (H) mm

Weight 345g (G.W.) (w/o Batteries)
Accessories Adult: 24~40 cm (9.4"~15.7")
Limited Users Type BF-Device and cuff are designed to provide special protection against electrical shocks.

IP Classification IP21: Protection against harmful ingress of water and particulate matter

*Specifications are subject to change without notice.

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- EN

1. Arm Cuff
2. LCD Display
3. Air Tube and Connector
4. Memory Key
5. ON/OFF/START Key
6. Battery Cover
- FR

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- EN

1. Movement Mark
2. Hypertension Risk Indicator
3. Weak Battery Mark
4. Irregular Heartbeat Detector (IHB)
5. Memory Mark
6. Memory Average
7. Systolic Pressure
8. Diastolic Pressure
9. Pulse Rate
10. Pulse Mark
- FR

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