Model: K150

www.accumed.ch

English

Introduction

Blood pressure measurements determined with K150 are equivalent to those obtained by a trained observer using cuff/stethoscope auscultation method, within the limits prescribed by the American National Standard, Electronic or Automated Sphygmomanometers. This unit is to be used by adult consumers in a home environment. Do not use this device on infants or neonates. K150 is protected against manufacturing defects by an established International Warranty Program. For warranty information, you can contact your local distributor.

Attention: Consult the accompanying documents. Please read this manual carefully before use. For specific information on your own blood pressure, contact your physician. Please be sure to keep this manual.

How This Unit Works

K150 uses the oscillometric method to detect your blood pressure. Before the cuff starts inflating, the device will establish a baseline cuff pressure, which is equivalent to the air pressure. The measurement of your blood pressure is based on this baseline pressure. After the cuff inflates to block your blood in the artery, the deflation process starts. During the deflation of the wrist cuff, the monitor is detecting the pressure oscillations generated by the beat-to-beat pulsatile. Any muscle movement during this period of time will cause measurement error. After detecting the amplitude and the slope of the pressure oscillations during the deflation process, your K150 will determine for you the systolic and diastolic pressures, and your pulse rate is detected at the same time.

Preliminary Remarks

This Blood Pressure Monitor complies with the European regulations and bears the CE mark "CE 0120". The quality of the device has been verified and conforms to the provisions of the EC council directive 93/42/EEC (Medical Device Directive), Annex I essential requirements and applied harmonized standards.

EN 1060-1; 1995/A2: 2009 Non-invasive sphygmomanometers - Part 1 - General requirements and CE marking

EN 1060-3; 1997/A2: 2009 Non-invasive sphygmomanometers - Part 3 - Supplementary requirements for electro-mechanical blood pressure measuring systems

EN 1060-4; 2004 Non-invasive sphygmomanometers - Part 4- Test Procedures to determine the overall system accuracy of automated non-invasive sphygmomanometers

This blood pressure monitor was designed for long service time. To ensure accurate measurements, this monitor is recommended to be re-calibrated every two years.

Blood Pressure Standard

The National High Blood Pressure Education Program Coordinating Committee has developed a blood pressure standard, classifying blood pressure ranges into 4 stages. (Ref. The Seventh Report of the Joint National Committee on Prevention, Detection, Evaluation, and Treatment of High Blood Pressure - Complete Report JNC-7, 2003). This blood pressure classification are based on historical data and may not be directly applicable to any particular patient. It is important that you consult with your physician regularly. Your physician will tell you your normal blood pressure range as well as the point at which you will be considered at risk. For reliable monitoring and reference of blood pressure, keeping long-term records is recommended. Please download the blood pressure log at our website www.accumed.ch.

Display Explanations

E1 / Measurement Error: Measure again. Wap the cuff correctly and keep wrist steady during measurement. If error keeps occurring, return the device to your local distributor or service center.

E1 / Air Circuit Abnormally: Measure again. If the error keeps recurring, return the device to your local distributor or service center.

E2 / Pressure Exceeding 300 mmHg: Switch the unit off and measure again. If the error keeps occurring, return the device to your local distributor or service center.

E3 / Data Error: Remove and reload the batteries. If the error keeps occurring, return the device to your local distributor or service center.

E4 / Exceeding Measurement Range: Measure again. If the error keeps occurring, return the device to your local distributor or service center.

Hypertension Risk Indication

The National High Blood Pressure Education Program Coordinating Committee has developed a blood pressure standard, classifying blood pressure ranges into 4 stages. This unit is equipped with innovative blood pressure risk indicator, which visually indicates the assumed risk level (normal / prehypertension / stage 1 hypertension / stage 2 hypertension ☹️) of the result after each measurement.

Irregular Heartbeat (IHB) Detection
This unit is equipped with an Irregular Heartbeat (IHB) Detector which allows those who have an irregular heartbeat to obtain accurate measurements alerting the user of the presence of an irregular heart beat during the measurement.

Note: It is strongly recommended that you consult your physician if the IHB icon ☹️ appears often.

Installing Batteries

Open the battery cover in the right direction. Install the battery in the right position. Battery Type: 2 AAA size. Replace the cover and click in the other end to secure the battery cover. Replace the batteries if:

- The weak battery mark appears in the display.
- Nothing appears in the display when the ON/OFF/START key is pressed.
- Use exclusively brand batteries. Always replace with new batteries together. Use batteries of the same brand and same type.
- Apply the Cuff
 - Remove all watches, jewelry, etc. prior to attaching the wrist monitor. Clothing sleeves should be rolled up and the cuff should be wrapped on bare skin for correct measurements.
 - Apply cuff to left wrist with palm facing up as Ⓓ.
 - Make sure the edge of the cuff is about 1 cm from the palm as Ⓓ.
 - In order to ensure accurate measurements, fasten the wrist strap securely around your wrist so there is no extra space between the cuff and the wrist as Ⓖ. If the cuff is not wrapped tight enough, the measurement values may be false.
 - If your physician has diagnosed you with poor circulation on your left arm, carefully place the cuff around your right wrist as shown in Ⓔ.

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
- Electromagnetic interference: The device contains sensitive electronic components. Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave ovens). These may lead to temporary impairment of measurement accuracy.
- Dispose of device, batteries, components and accessories according to local regulations.
- This monitor may not meet a performance specification if stored or used outside temperature and humidity ranges specified in Specifications.
- Batteries can be fatal if swallowed. You should therefore store the batteries and products where they are inaccessible to small children. If a battery has been swallowed, call a doctor immediately.

20. Specifications	
Power source	DC 3V, two AAAA batteries
Measurement method	Oscillometric
Measurement range	Pressure: 40-250 mmHg; Pulse: 40-199 beats/minute
Accuracy	Pressure: ±3 mmHg; Pulse: within ±5% of reading
Pressure sensor	Semiconductor
Inflation	Pump driver
Deflation	Automatic Pressure release valve
Memory capacity	90 memories
Auto-on/Off	Auto-on/Off after 147 key operation
Operation environment	10°C–40°C (50°F–104°F); 40%–85% RH; 700–1060 hPa
Storage and Transportation Environment	–10°C–60°C (14°F–140°F); 10%–90% RH; 700–1060 hPa
Dimensions	82(L) x 69(W) x 66(H) mm
Weight	115g (lwo batteries)
Wrist circumference	132–222mm(5.3"-8.7")
Wrist users	Adults
Adult users	Type: BF; Device and cuff are designed to provide special protection against electrical shocks
IP Classification	IP22, Protection against harmful ingress of water and particulate matter

*Specifications are subject to change without notice.

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