
Digital Blood Pressure Arm Monitor

Model

01-513BL (BLUE); 01-513GR (GREEN);

01-513PK (PINK); 01-513SL (SILVER)


INSTRUCTION MANUAL • ENGLISH & ESPAÑOL

STOP!
PLEASE ENSURE YOU HAVE ALL OF THE
FOLLOWING COMPONENTS BEFORE USING
YOUR DIGITAL BLOOD PRESSURE MONITOR


Digital Monitor


Pre-Formed
Arm Cuff

Storage Case

4 AAA Batteries

Instruction Manual
Quick Start Guide
Blood Pressure Reading
Tracking Log
FAQ Information Sheet

**IF YOU ARE MISSING ANY PARTS, INCLUDING
INSERTS OR INSTRUCTION MANUALS, DO NOT
RETURN TO PLACE OF PURCHASE. CONTACT
CUSTOMER CARE AT 866-326-1313.**

INDEX

Care & Safety Information.....	4-5
Introduction	6
Digital Blood Pressure Monitor Features	7
Getting Started.....	8
Battery Installation	9
Setting Date and Time	10-11
Applying Your Blood Pressure Cuff	12-13
Taking Your Blood Pressure Reading	14-15
Interpreting Your Results.....	16-17
Memory Functions	18-19
Care and Maintenance	20-21
Device & Label Symbols.....	22
Display Symbols.....	23
Troubleshooting and Error Codes	24-25
FCC Statement	26-27
Electromagnetic Compatibility.....	28-31
Product Specifications	32
Warranty	33-34
Instrucciones en Español.....	35-70

Toll-Free Customer Care Help Line: 1-866-326-1313
Monday – Friday 8:30 a.m. – 4:30 p.m. CST

Manufactured for:
Veridian Healthcare
1175 Lakeside Drive
Gurnee, IL 60031
www.veridianhealthcare.com

Made in China
#93-1179 02/14
©2014 Veridian Healthcare, LLC

CARE & SAFETY INFORMATION

Note: Read all instructions carefully before use. The following basic precautions are needed when using an electrical product.

Caution: Failure to read and observe all precautions could result in personal injury or equipment damage.

Improper care or use of your blood pressure monitor may result in injury, damage to the unit or ineffective treatment. Following these instructions will ensure the blood pressure monitor's efficacy and long life.

General Cautions & Warnings

- Measurement results should be discussed with your physician or healthcare professional; never self-diagnose or attempt treatment as this may be dangerous
- Intended for adult use only; this device is not approved for infant or child use
- Keep out of reach of children
- Do not use the unit if it has any damaged parts, if it has been submersed in water or dropped
- If any abnormality occurs, discontinue use until the unit has been examined and repaired
- Only use included attachments and parts; do not use attachments from other brands or models as these may not be properly calibrated for use with this device and may result in measurement error
- Do not disassemble the unit or attempt to repair it; substitution of a component different from that supplied may result in measurement error and will void manufacturer warranty

CARE & SAFETY INFORMATION

- Always follow local regulations for proper disposal of the monitor, cuff and batteries

Operating Cautions & Warnings

- Prolonged over-inflation may cause congestion, swelling or bruising in some people
- If you experience any discomfort or an abnormality occurs, stop using the device immediately; to stop the device during operation, press the START/STOP button and the cuff will automatically deflate
- Operation outside of stated operating temperature may result in measurement error or device malfunction; operation environment temperature is: 41°F – 104°F (5°C – 40°C); Humidity: 80% RH max
- Never use this unit while operating a vehicle or in the bathtub or shower

Storage Cautions & Warnings

- Storage outside of stated storage temperature may result in measurement error or device malfunction; storage environment temperature is: -4°F – 140°F (-20°C – 60°C); Humidity: 10-93%
- Keep the unit out of reach of small children
- Always keep the unit unplugged while not in use; remove the batteries if the unit will not be used for an extended period of time

Cleaning Cautions & Warnings

- Never immerse the unit in water to clean as it may damage the unit
- Follow the 'Cleaning and Maintenance' portion of this manual for instruction on how to clean and care for your monitor

INTRODUCTION

It is recommended that you first seek the advice and recommendation of your physician or healthcare professional when using home diagnostic devices, including blood pressure monitors.


Automatic digital blood pressure monitors use the oscillometric method to electronically measure your blood pressure. The monitor detects your blood's movement through the artery in your arm and converts the movements into a digital reading. The oscillometric method does not require a stethoscope, making the monitor ideal for home use.

Blood pressure readings determined with this device are equivalent to measurements obtained by a trained healthcare professional using the cuff/stethoscope auscultation method, within the limits prescribed by the American National Standard for Electronic or Automated Sphygmomanometers.

Thank you for purchasing a Veridian Healthcare Digital Blood Pressure Monitor. With proper care and use, your monitor will provide you with many years of reliable readings.

**PLEASE READ THIS INSTRUCTION MANUAL
COMPLETELY BEFORE OPERATING THIS UNIT.**

MONITOR FEATURES

DETAILS OF THE DISPLAY AND INTERPRETATION OF THE SYMBOLS ARE INCLUDED THROUGHOUT THIS MANUAL.

GETTING STARTED

Before Taking a Measurement

1. Your readings should only be interpreted by your physician or healthcare professional with access to your individual medical history. Regular use of a home blood pressure monitor will allow you to track and record your readings for discussion with your physician.
2. Conduct your measurement in a quiet place while seated in a relaxed position. Rest for 15 minutes before taking your reading. If an error occurs or you desire to take a second reading, allow 15 minutes between readings for your blood vessels to return to normal.
3. Avoid smoking, eating, taking medication, alcohol, physical activity or any other stressful activity for 30 minutes prior to taking a reading.
4. Always remove any jewelry or constrictive clothing that may interfere with the cuff placement.
5. Keep yourself and the monitor still during measuring; do not talk during the reading.
6. It is recommended that you take your readings at the same time each day to better monitor any indications in your results.
7. Record your daily measurement on the included chart or some other written document to share with your physician.

BATTERY INSTALATION


This unit comes complete with 4 AAA batteries.

Replace the batteries when any of the following occur:

- Low Battery symbol and 'Lo' appears on the display
- The display backlight dims or does not illuminate
- Monitor does not turn on after the START button is pressed

1. The battery compartment is located on the underside of the monitor as shown in the below illustration. Carefully slide the cover away from the housing and the cover will open.
2. Insert or replace 4 AAA batteries into the battery compartment, ensuring to match the indicated polarity symbols. Always use new batteries.
3. Replace the battery cover.
4. Dispose of batteries according to local disposal and recycling regulations.

It is recommended to remove the batteries if the unit will not be used for an extended period of time.

SETTING DATE/TIME

IT IS NECESSARY TO SET THE DATE AND TIME FOR THE UNIT EVERY TIME BATTERIES ARE INITIALLY INSTALLED OR REPLACED.

The date/time feature is on the lower right of the display screen. This unit stores the date and time, including year, allowing for a date/time stamp on all readings in the memory recall function.

1. When the unit is off (a blank display screen), press and hold the SET button for 2-3 seconds to activate the display. The display will illuminate and the year will begin flashing.

2. **MONTH / DATE / YEAR**—press the MEM button to advance the display to the desired year, press SET to confirm your selection.

Repeat for Month and Date-MEM to advance, SET to confirm.

SETTING DATE/TIME

3. **HOUR**—press the MEM button to advance the display to the desired hour, press SET to confirm your selection.
4. **MINUTES**—press the MEM button to advance the display to the desired minutes, press SET to confirm your selection.
5. After setting the final step, the display will briefly show 'done' and then automatically shut off.

The monitor is now ready for use.

APPLYING YOUR ARM CUFF

PROPER CUFF FIT AND PLACEMENT IS CRITICAL IN OBTAINING ACCURATE BLOOD PRESSURE MEASUREMENTS. PLEASE READ THESE INSTRUCTIONS CAREFULLY AND CONSIDER CONSULTING WITH YOUR PHYSICIAN FOR FURTHER DIRECTION.

To determine the proper cuff size, measure the circumference of your left upper arm at the point midway between inside bend of your elbow and your shoulder, Fig 1.

Adult Cuff fits arm circumference:
8-3/4" - 16-1/2" (22-42 cm).


Fig. 1

Reminder, it is important to avoid smoking, eating, taking medication, alcohol consumption or physical activity 30 minutes prior to taking a reading.

If for any reason you are unable to or should not use your left arm, please modify the instructions for cuff application to your right arm. Your physician can help you identify which arm is best for you to take measurements from.

The cuff should not be plugged into the monitor until after the cuff is applied to your arm to avoid pulling or yanking the monitor.

APPLYING YOUR ARM CUFF

1. Remove any constrictive clothing or jewelry that may interfere with cuff placement.
2. Be seated at a table or desk with your feet flat on the floor.
3. Insert your left arm into the cuff cylinder, keeping the air hose on the inside of the arm, Fig 1.
4. Align the hose with your brachial artery, slightly to the right on the inside of the bend of the elbow, Fig 2.
5. The bottom edge of the cuff should be positioned approximately one inch above the elbow joint, Fig 3.
6. The cuff should fit comfortably, yet snugly around your arm. You should be able to insert one finger easily between your arm and the cuff.


Fig. 1


Fig. 2


Fig. 3

If you are not comfortable with applying your cuff, please seek the assistance of another member of your household or work with your physician to practice the cuff application. Incorrectly applied cuffs may result in inaccurate readings.

TAKING YOUR MEASUREMENT

Please read the preceding portions of this manual prior to taking your first reading.

Reminder, it is important to avoid smoking, eating, taking medication, alcohol consumption or physical activity 30 minutes prior to taking a reading.


Fig. 1

1. Rest your elbow on a solid surface with your palm facing upward. Elevate your arm so that the cuff is at the same level as your heart, Fig 1. Be certain that the cuff is elevated to the proper level; you may choose to support your hand with a folded towel or book. Relax your hand.
2. Press the START button to turn the unit on. The unit will run a brief self-test, and the display will briefly show '0'.

3. The cuff will automatically begin inflating. The numbers on the display screen will begin to go up. As the heart rate is being measured, the heart icon will also flash.

TAKING YOUR MEASUREMENT

- When the measurement is complete, your blood pressure and pulse reading results will show on the display screen at the same time.

- The Hypertension Indicator on the left side of the display will indicate Normal, Mild, Moderate or Severe. See page 17 for information regarding interpreting the Hypertension Indicator results.
- If an irregular heartbeat was detected during the reading, the Irregular Heartbeat Detector indicator will appear flashing on the display, Fig 5. See page 16 for more information on the Irregular Heartbeat Detector feature.

- The reading will automatically be stored in memory, up to 60 readings.
- Select START to turn the unit off and conserve energy and battery life. The unit will automatically shut-off after approximately 1 minute.


YOU CAN STOP THE INFLATION OR DEFLATION PROCESS AT ANY TIME BY PRESSING THE STOP BUTTON.

INTERPRETING YOUR RESULTS

Irregular Heartbeat Detection

An irregular heartbeat is defined as a heartbeat rhythm that is more than 25% slower or 25% faster than the average measurement detected while the monitor is conducting the measurement. When this device detects an irregular rhythm two or more times during the reading, the IHB symbol will appear on the display screen.


This monitor will operate if the icon appears, but the results may be compromised, especially if this symbol appears often. Please consult with your physician or trained healthcare professional for further information regarding an irregular heartbeat and if this symbol appears frequently.

INTERPRETING YOUR RESULTS

Hypertension Indicator

The graphic on the left side of the display screen indicates where the measurement falls according to the World Health Organization standards.

	Normal	Mild	Moderate	Severe
SYS	<140	140–159	160–179	≥180
DIA	<90	90–99	100–109	≥110

American Heart Association Classifications

This chart can help you identify where a given blood pressure measurement falls according to American Heart Association standards. This chart should only be considered as a guideline; always consult with your physician to interpret your individual results. Never adjust your medication(s) dosage, or severely alter your diet or exercise routine without consulting with your physician.

Category	Systolic		Diastolic
Normal	< 120	and	<80
Prehypertension	120-139	or	80-89
High Blood Pressure, Stage 1	140-159	or	90-99
High Blood Pressure, Stage 2	≥160	or	<100
Hypertensive Crisis	≥180	or	<110

MEMORY FUNCTIONS

Recalling Measurements in Memory

This monitor stores and recalls up to 60 measurements, plus an average of the last three measurements in memory. These measurements can be shared with your physician or trained healthcare professional for use in evaluating your overall health.

1. When the monitor is off, press and release the MEM button. The average blood pressure and pulse measurement of the most recent three readings will show with AVG icon on the screen.

2. Press the MEM button again to begin scrolling through all of the measurements from most current to oldest; an M icon and the measurement sequence number will appear on the lower right of the display. For example, 8 / 30 indicates you are viewing measurement 8 of 30 recorded)

Each measurement will indicate the date and time taken (alternating), the measurement and pulse reading, Hypertension Indicator results and Irregular Heartbeat if present.

MEMORY FUNCTIONS

NOTE: The most recent record (1) is shown first. When the total measurements in memory exceed 60 the last record (60) is dropped from the list..

Clearing Measurements from Memory

To delete all measurements in memory, follow these steps:

1. When the monitor is off (display blank), press and hold the MEM button 3–4 seconds. The display will flash 'dEL ALL'

2. Press SET to confirm deleting all measurements. The display will flash 'dEL' and 'done'.
3. The memory bank is now clear.

CARE & MAINTENANCE

Proper care and maintenance of your blood pressure monitor will help ensure long life for the product; improper use of the monitor may void the included warranty.

Monitor Care and Maintenance

- Do not use any liquids on the monitor; use a soft, dry cloth to clean the monitor as needed
- Storage outside of stated storage temperature may result in measurement error or device malfunction; storage environment temperature is: -4°F – 140°F (-20°C – 60°C); humidity: 10-93% RH
- Do not store in direct sunlight
- Always keep the unit unplugged while not in use; remove the batteries if the unit will not be used for an extended period of time
- Do not attempt to disassemble your monitor; disassembling your monitor may void the manufacturers warranty
- Do not subject the monitor to strong shocks; take care not to drop the monitor
- With proper care, the monitor will maintain performance for a minimum of 10,000 measurements

CARE & MAINTENANCE

Cuff Care and Maintenance

- Do not attempt to remove the cuff from the monitor housing for any reason
- The cuff should not be submerged or washed in a washing machine, but may be spot cleaned with a mild-detergent as needed; it is suggested that the cuff be cleaned once a month if the device is used daily
- To disinfect the cuff, wipe the inside (skin contact side) of the cuff lightly moistened with 70-90% ethyl rubbing alcohol and allow to air dry; it is suggested to disinfect the cuff frequently when used in clinical settings or if the cuff is contaminated in any way
- Never attempt to iron or dry the cuff; if the cuff becomes wet, allow the cuff to air dry completely before use or storage

DEVICE & LABEL SYMBOLS

These symbols may appear on your device, instructions or packaging and may vary by make and model.

SYMBOL	MEANING
	Read This Manual—All included manuals should be read prior to first use
	Warning—Symbol indicates a warning, prohibition or mandatory action that mitigates a risk that is not necessarily obvious to the device operator
	Type BF Applied Parts—Indicates that a part of this unit comes in contact with the patient in order to carry out it's intended function; in the case of this device the cuff is the Type BF Applied Part
	Environment Protection—dispose of this product properly; consult with your local recycling ordinances for proper recycling and disposal
	Keep Dry—This device should be kept dry; never submerge the unit or cuff. Consult with the Care and Maintenance section of this manual for information on cleaning your monitor
	Polarity Symbol—Polarity of D.C. power connector

DISPLAY SYMBOLS

SYMBOL	MEANING
SYS	Systolic Blood Pressure—The ‘top’ blood pressure result
DIA	Diastolic Blood Pressure—The ‘bottom’ blood pressure result
Pul/min	Pulse Rate—Pulse rate per minute during the measurement
 +Lo	Low Battery—Battery is low and should be replaced
mmHg	Millimeters of Mercury—Unit of blood pressure measurement for this device
	Irregular Heartbeat Icon—An irregular heartbeat was detected during the measurement
M + 	Memory—The display shows the measurements stored in memory
	Average—The display shows the average of the last three measurements in memory
	Movement—There was excessive movement during measurement and results may be affected

TROUBLESHOOTING

PROBLEM	POSSIBLE CAUSE	SOLUTION
Display shows an abnormal result	Cuff position was not correct or was not tight enough	Apply the cuff properly and attempt a new measurement
	Body posture was not correct during measurement	Review the 'Taking a Measurement' portion of this manual and attempt a new measurement
	Speaking, excessive movement, angry, excited or nervous anxiety during measurement	Wait a period of time and attempt a new measurement after calm has returned; do not speak or move during measurement
Display is blank when power is on	Batteries may be expired or installed improperly	Check the polarity of the batteries and reinstall if necessary; replace batteries with new batteries

TROUBLESHOOTING - ERROR CODES

CODE	CAUSE	SOLUTION
E 1	The cuff is not secure/not tight enough.	Re-apply the cuff and attempt new measurement.
E 2	The cuff is wrapped too tightly.	Re-apply the cuff and attempt new measurement.
E 3	The cuff pressure is over-inflated.	Relax for several minutes and attempt new measurement
E 10 E 11	Excessive movement was detected during measurement.	Remain still and attempt new measurement.
E 20 E 21	An error occurred during measurement and the results are inaccurate.	Relax for several minutes and attempt new measurement
Eexx	A calibration error has occurred.	A calibration error has occurred.

If an error code persists or is not remedied, please contact Customer Care at 866-326-1313.

FCC STATEMENT

NOTE:

POTENTIAL FOR RADIO/TELEVISION INTERFERENCE (for U.S.A. only)

This product has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules.

These limits are designed to provide reasonable protection against harmful interference in a residential installation. The product generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If the product does cause harmful interference to radio or television reception, which can be determined by turning the product on and off, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna
- Increase the separation between the product and the receiver.

FCC STATEMENT

- Connect the product into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

POTENTIAL FOR RADIO/TELEVISION INTERFERENCE

(for Canada only)

This digital apparatus does not exceed the Class B limits for radio noise emissions from digital apparatus as set out in the interference-causing equipment standard entitled “Digital Apparatus”, ICES-003 of the Canadian Department of Communications.

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Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

ELECTROMAGNETIC COMPATIBILITY

Table 1

For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacture's declaration - electromagnetic emissions

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

Emmission Test	Compliance	Electromagnetic Environment Guidance
RF emissions CISPR 11	Group 1	This monitor uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	This monitor is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

ELECTROMAGNETIC COMPATIBILITY

Table 2
For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacture's declaration - electromagnetic immunity

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 s	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of this monitor requires continued operation during power mains interruptions, it is recommended that this monitor be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

ELECTROMAGNETIC COMPATIBILITY

Table 3

For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING

Guidance and manufacturer's declaration - electromagnetic immunity

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of this monitor, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance : $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P} \text{ 80 MHz to 800 MHz}$ $d = 2.3 \sqrt{P} \text{ 800 MHz to 2,5 GHz}$ Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,* should be less than the compliance level in each frequency range.** Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

* Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which this monitor is used exceeds the applicable RF compliance level above, this monitor should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating this monitor.

** Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

ELECTROMAGNETIC COMPATIBILITY

Table 4

For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and this monitor

This monitor is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this monitor can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this monitor as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2.3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

PRODUCT SPECIFICATIONS

Name	Automatic Digital Blood Pressure Arm Monitor
Model Number	01-513BL, 01-513GR, 01-513PK, 01-513SL
Display System	Digital LCD display; 1.5" x 5.6"
Measuring Method	Oscillometric Method, Automatic Air Inflation and Measurement
Power Source	4 x 1.5V  SIZE AAA
Measuring Range	Pressure: 0-300 mmHg; . Pulse: 40-199 beats/minute
Accuracy	Pressure: ± 3 mmHg; Pulse: $\pm 5\%$
Deflation	Automatic pressure release valve
Memory	60 measurements
Automatic Power-Off	Approximately 1 minute after last operation
Battery Life	Approximately 1 month with normal usage
Operation Environment	Temperature 41°F – 104°F (5°C – 40°C); Humidity $\leq 80\%$ RH max
Storage Environment	Temperature -4°F – 140°F (-20°C – 60°C); Humidity 10% – 93% RH max
Monitor Dimensions	7" x 3-7/8" x 1-1/2"
Monitor Weight	12.1 oz. (without batteries)
Arm Circumference	8-3/4" - 16-1/2" (22-42 cm)
Accessories	Adult arm cuff, storage case, Instruction Guidebook, four AAA alkaline batteries, Quick Start Guide, Blood Pressure Log

Specifications are subject to change without notice

LIFETIME LIMITED WARRANTY

Congratulations on your purchase of a Veridian Digital Blood Pressure Monitor. Your Digital Blood Pressure Monitor is covered by the following limited warranty commencing upon the date of purchase, and subject to the following terms and conditions:

Veridian Healthcare warrants that its Digital Blood Pressure Monitor will be free from defects in materials and workmanship under normal consumer usage for as long as the original purchaser owns the product.

Periodic maintenance, repair and replacement of parts due to normal wear and tear are excluded from coverage. Defects or damage that result from: (a) improper operation, storage, misuse or abuse, accident or neglect, such as physical damage (cracks, scratches, etc.) to the surface of the product resulting from misuse; (b) contact with liquid, water, rain, extreme humidity or heavy perspiration, sand, dirt or the like, extreme heat, or food; (c) use of the Digital Blood Pressure Monitor for commercial purposes or subjecting the Digital Blood Pressure Monitor to abnormal usage or conditions; or (d) other acts which are not the fault of Veridian, are excluded from coverage. This warranty does not cover batteries or other power sources that may be provided with, or used with the Digital Blood Pressure Monitor.

If the Digital Blood Pressure Monitor fails to conform to this limited warranty, return the Digital Blood Pressure Monitor postage prepaid along with \$5.00 for return shipping, handling and insurance to: Veridian Healthcare, Attn: Repair Department, 1175 Lakeside Drive, Gurnee, IL, 60031. When returning a product, please also include: (i) a copy of your receipt, bill of sale or other comparable proof of purchase; (ii) a written description of the problem; and (iii) your name, address and telephone

LIFETIME LIMITED WARRANTY

number. Carefully package the product to avoid any damage that may occur while in transit; shipping insurance with returned receipt is recommended. At our option, Veridian Healthcare will repair or replace the unit found to be defective in materials or workmanship under normal consumer usage. The purchaser will be notified of any additional repairs required prior to completing the repair, and will be responsible for parts charges, if any, and repair charges not covered by this limited warranty.

EXCEPT AS PROVIDED FOR IN THIS LIMITED WARRANTY, ALL EXPRESS AND IMPLIED WARRANTIES AND CONDITIONS ARE DISCLAIMED, INCLUDING WITHOUT LIMITATION THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. THE REPAIR OR REPLACEMENT AS PROVIDED UNDER THIS LIMITED WARRANTY IS THE EXCLUSIVE REMEDY OF THE CONSUMER, AND IS PROVIDED IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. IN NO EVENT SHALL VERIDIAN HEALTHCARE BE LIABLE, WHETHER IN CONTRACT OR TORT (INCLUDING NEGLIGENCE) FOR DAMAGES IN EXCESS OF THE PURCHASE PRICE OF THE PRODUCT, OR FOR ANY INDIRECT, INCIDENTAL, SPECIAL OR CONSEQUENTIAL DAMAGES OF ANY KIND, OR FOR DAMAGES TO, OR LOSS OF, OTHER PROPERTY OR EQUIPMENT OR PERSONAL INJURIES TO THE FULL EXTENT THESE DAMAGES MAY BE DISCLAIMED BY LAW.

Some states and jurisdictions do not allow the limitation or exclusion of incidental or consequential damages, or limitation on the length of an implied warranty, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that vary from state to state or from one jurisdiction to another.

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