Jumbo Screen Premium Digital Blood Pressure Arm Monitor

Model 01-514


INSTRUCTION MANUAL • ENGLISH & ESPAÑOL


STOP!
PLEASE ENSURE YOU HAVE ALL OF THE FOLLOWING
COMPONENTS BEFORE USING YOUR DIGITAL
BLOOD PRESSURE MONITOR

Digital Monitor


Adult Arm Cuff

4 AA Batteries

Canvas Storage Bag

Instruction Manual
Quick Start Guide
Blood Pressure Reading Tracking Log
FAQ Information Sheet
Replacement Parts Order Form

**IF YOU ARE MISSING ANY PARTS, INCLUDING INSERTS OR
INSTRUCTION MANUALS, DO NOT RETURN TO PLACE OF
PURCHASE. CONTACT CUSTOMER CARE AT 866-326-1313.**

INDEX

Care and Safety Information	4-5
Introduction	6
Digital Blood Pressure Monitor Features	7
Getting Started.....	8
Battery Installation / AC Power	9
Setting Date and Time.....	10
Applying Your Blood Pressure Cuff	11-12
Taking Your Blood Pressure Reading	13-14
Interpreting your Results.....	15-16
Memory Functions	17-18
Care and Maintenance	19-20
Display and Label Symbols	21
Display Symbols.....	22
Error Codes & Troubleshooting	23
FCC Statement.....	24
Electromagnetic Compatibility	25-28
Specifications	29
Warranty	30
Instrucciones en Español.....	33-64

Toll-Free Customer Care Help Line: 1-866-326-1313
Monday – Friday 8:30 a.m. – 4:30 p.m. CST

Manufactured for:
Veridian Healthcare
1175 Lakeside Drive
Gurnee, IL 60031
www.veridianhealthcare.com

Made in China

#93-1058 09/14
©2014 Veridian Healthcare, LLC

CARE & SAFETY INFORMATION

NOTE: Read all instructions and inserts included with this monitor carefully before use. The following basic precautions are needed when using an electrical product.

CAUTION: Failure to read and observe all precautions could result in personal injury or equipment damage.

Improper care or use of your blood pressure monitor may result in injury, damage to the unit or ineffective treatment. Following these instructions will ensure the blood pressure monitor's efficacy and long life.

GENERAL CAUTIONS AND WARNINGS

- Measurement results should be discussed with your physician or healthcare professional; never self-diagnose or attempt treatment as this may be dangerous
- Intended for adult use only; this device is not approved for infant or child use
- Keep out of reach of children
- Do not use the unit if it has any damaged parts, if it has been submersed in water or dropped
- If any abnormality occurs, discontinue use until the unit has been examined and repaired
- Only use included attachments and parts; do not use attachments from other brands or models as these may not be properly calibrated for use with this device and may result in measurement error
- This monitor may not be appropriate for individuals with serious arrhythmias; consult with your physician for proper monitor use and evaluation.
- Check the blood circulation of the user before measurement
- Do not share the cuff with an infected person.
- Do not disassemble the unit or attempt to repair it; substitution of a component different from that supplied may result in measurement error and will void manufacturer warranty
-  Always follow local regulations for proper disposal of the monitor, cuff and batteries

CARE & SAFETY INFORMATION

OPERATING CAUTIONS AND WARNINGS

- Prolonged overinflation may cause congestion, swelling or bruising in some people
- Wait a minimum of 2 minutes between measurements (20-minutes is recommended); excessive frequent measurements may restrict blood circulation
- Do not use this device while connected to any monitoring equipment on the same limb
- Consult with your healthcare professional for guidance for use if you have any of the following: intravascular shunt; mastectomy on the arm corresponding to the measurement side of the body
- If you experience any discomfort or an abnormality occurs, stop using the device immediately; to stop the device during operation, press the ⏻ (POWER) button and the cuff will automatically deflate
- Operation outside of stated operating temperature may result in measurement error or device malfunction; operation environment temperature is: 41°F – 104°F (5°C – 40°C); Humidity: ≤90% RH
- Never use this unit while operating a vehicle or in the bathtub or shower

STORAGE CAUTIONS AND WARNINGS

- Storage outside of stated storage temperature may result in measurement error or device malfunction; storage environment temperature is: -4°F – 131°F (-20°C – 55°C); Humidity: ≤90% RH
- Keep the unit out of reach of small children
- Always keep the unit unplugged while not in use; remove the batteries if the unit will not be used for an extended period of time

CLEANING CAUTIONS AND WARNINGS

- Never immerse the unit in water to clean as it may damage the unit
- Follow the 'Cleaning and Maintenance' portion of this manual for instruction on how to clean and care for your monitor

INTRODUCTION

It is recommended that you first seek the advice and recommendation of your physician or healthcare professional when using home diagnostic devices, including blood pressure monitors.


Automatic digital blood pressure monitors use the oscillometric method to electronically measure your blood pressure. The monitor detects your blood's movement through the artery in your arm and converts the movements into a digital reading. The oscillometric method does not require a stethoscope, making the monitor ideal for home use.

Blood pressure readings determined with this device are equivalent to measurements obtained by a trained healthcare professional using the cuff/stethoscope auscultation method, within the limits prescribed by the American National Standard for Electronic or Automated Sphygmomanometers.

Thank you for purchasing a Digital Blood Pressure Monitor. With proper care and use, your monitor will provide you with many years of reliable readings.

**PLEASE READ THIS INSTRUCTION MANUAL
COMPLETELY BEFORE OPERATING THIS UNIT.**

BLOOD PRESSURE MONITOR FEATURES

DETAILS OF THE DISPLAY AND INTERPRETATION OF THE SYMBOLS ARE INCLUDED THROUGHOUT THIS MANUAL.

GETTING STARTED

Before Taking a Measurement

1. Your readings should only be interpreted by your physician or healthcare professional with access to your individual medical history. Regular use of a home blood pressure monitor will allow you to track and record your readings for discussion with your physician.
2. Conduct your measurement in a quiet place while seated in a relaxed position. Rest for 15 minutes before taking your reading. If an error occurs or you desire to take a second reading, allow 15 minutes between readings for your blood vessels to return to normal.
3. Avoid smoking, eating, taking medication, alcohol, physical activity or any other stressful activity for 30 minutes prior to taking a reading.
4. Always remove any jewelry or constrictive clothing that may interfere with the cuff placement.
5. Keep yourself and the monitor still during measuring; do not talk during the reading.
6. It is recommended that you take your readings at the same time each day to better monitor any indications in your results.
7. Record your daily measurement on the included chart or some other written document to share with your physician.

BATTERY INSTALLATION / AC POWER

This unit comes complete with 4 AA batteries.

It is necessary to replace the batteries when the Low Battery symbol appears on the display or when the display does not turn on after the **POWER** button is pressed. 

Replacing Battery, Fig. 1

1. Press down and lift the cover.
2. Insert or replace 4 AA batteries into the battery compartment, ensuring to match the indicated polarity symbols. Always use new batteries.
3. Replace the battery cover.
4. It is recommended to remove the batteries if the unit will not be used for an extended period of time.


Fig. 1

AC Power Adapter

An AC adapter is available separately for battery-free operation. The AC adapter plugs into the right side of the unit, Fig. 2.

Please remember to practice common sense electrical safety precautions while using the adapter.

To order an AC adapter, please see the enclosed order form.


Fig. 2

AC Adapter Output: 6V DC 600mA

SETTING DATE/TIME

IT IS NECESSARY TO SET THE DATE AND TIME FOR THE MONITOR EVERY TIME BATTERIES ARE INITIALLY INSTALLED OR REPLACED.

DATE AND TIME WILL DEFAULT TO FACTORY SETTINGS, 12:00 A.M. 1M 1D, IF YOU DO NOT CHOOSE TO SET THE CORRECT INFORMATION.

1. When the unit is off, press and release both the  and **M** button at the same time. The month will begin to flash, Fig. 1.
2. Press the **M** button to change until the desired number appears. To scroll through the numbers quicker, hold down the M button as the numbers advance.
3. Press the  **POWER** button to confirm the month entry.
4. Next, the date will blink. Repeat steps 2 and 3 to set the date, then hours, then minutes. This unit does not track the year.
5. After setting the minutes, the unit will automatically exit out of the date/time setting mode.


Fig. 1

When the monitor is off, the date and time will still appear at the top of the screen.

Time is displayed in 24hr Military Time. This chart explains the time conversion.

MILITARY TIME = REGULAR TIME

1:00 = 1am	7:00 = 7am	13:00 = 1pm	19:00 = 7pm
2:00 = 2am	8:00 = 8am	14:00 = 2pm	20:00 = 8pm
3:00 = 3am	9:00 = 9am	15:00 = 3pm	21:00 = 9pm
4:00 = 4am	10:00 = 10am	16:00 = 4pm	22:00 = 10pm
5:00 = 5am	11:00 = 11am	17:00 = 5pm	23:00 = 11pm
6:00 = 6am	12:00 = 12pm	18:00 = 6pm	24:00 = 12am

APPLYING YOUR BLOOD PRESSURE CUFF

PROPER CUFF FIT AND PLACEMENT IS CRITICAL IN OBTAINING ACCURATE BLOOD PRESSURE MEASUREMENTS. PLEASE READ THESE INSTRUCTIONS CAREFULLY AND CONSIDER CONSULTING WITH YOUR PHYSICIAN FOR FURTHER DIRECTION.

Before applying your blood pressure cuff, be sure you have selected the appropriate size cuff.

To determine the proper cuff size, measure the circumference of your left upper arm at the point midway between inside bend of your elbow and your shoulder, Fig 1.

Adult Cuff fits arm circumference 11-3/4" - 16-1/2" (30-42 cm)

Additional cuff sizes are available for order on the enclosed order form.

If for any reason you are unable to or should not use your left arm, please modify the instructions for cuff application to your right arm. Your physician can help you identify which arm is best for you to take measurements from.

1. The cuff should not be plugged into the monitor until after the cuff is applied to your arm.
2. Remove any constrictive clothing or jewelry that may interfere with cuff placement.
3. Be seated at a table or desk with your feet flat on the floor and legs uncrossed.
4. Position the cuff on a solid surface with the tubing facing up and away from you. The metal ring/bar on the cuff should be to the left of the tubing.


Fig. 1

APPLYING YOUR BLOOD PRESSURE CUFF

5. Open the cuff by pulling or rolling the bottom of the cuff to the right, Fig 2. This should open the cuff without fully unrolling it, creating a cylinder. Do not extend the cuff beyond the metal bar.
6. Insert your left arm into the created cuff cylinder. Position the index mark Φ over the main artery on the inside of your arm, Fig 3, 4.
7. The bottom edge of the cuff should be positioned approximately one inch above the elbow joint, Fig 5.
8. The cuff should fit comfortably, yet snugly around your arm. You should be able to insert one finger easily between your arm and the cuff.

The cuff size is suitable for use when the vertical index mark Φ is within the horizontal OK range. A different size cuff is needed if the 'index' mark is outside the OK range.

If you are not comfortable with applying your cuff, please seek the assistance of another member of your household or work with your physician to practice the cuff application. Incorrectly applied cuffs may result in inaccurate readings.

Fig. 2


Fig. 3


Fig. 4


Fig. 5

TAKING YOUR BLOOD PRESSURE READING

Please read the preceding portions of this manual prior to taking your first reading.

Reminder, it is important to avoid smoking, eating, taking medication, alcohol consumption or physical activity 30 minutes prior to taking a reading.


Fig. 1

1. Position the monitor on a flat, stable surface with the digital display in view.
2. Insert the cuff tubing connector into the port on the left side of your monitor, Fig 1.
3. Rest your elbow on a solid surface with your palm facing upward. Elevate your arm so that the cuff is at the same level as your heart, Fig 2. Relax your left hand.
4. Press the **POWER** button to turn the unit on; the display will illuminate and the unit will run a self-test, Fig 3.
5. The display will show the current memory bank in use, U1 or U2, Fig 4. Press the M button to select a different memory bank if desired.


Fig. 2


Fig. 3


Fig. 4

TAKING YOUR BLOOD PRESSURE READING

YOU CAN STOP THE INFLATION OR DEFLATION PROCESS AT ANY TIME BY PRESSING THE START/POWER BUTTON.

6. The unit will briefly display the last reading value for that memory bank. If there is no previous measurement, the unit will display '0' for all of the reading values, Fig 5.
7. The cuff will automatically inflate to approximately 190mmHg or 40 mmHg higher than the last systolic pressure reading stored in memory.
8. When the inflation has reached optimum level, the display will begin to show the decreasing pressure while you feel the pressure of the cuff decrease.
9. When the measurement is complete your blood pressure measurement and pulse will display simultaneously on the screen, Fig 6. Both the systolic and diastolic readings will display at the same time.
10. The Hypertension Indicator will indicate your reading range in the left side of the display, Fig 7. See page 16 for a chart interpreting the Hypertension Indicator results.
11. If an irregular heartbeat was detected during the reading, the Irregular Heartbeat Detector will appear on the display, Fig 7. See page 15 for more information on the Irregular Heartbeat Detector feature.
12. The results will automatically be stored in the memory bank selected before the measurement, up to 60 readings.
13. Select **POWER** to turn the unit off and conserve energy and battery life. The unit will automatically shut-off and revert to time/date mode after approximately 1 minute.
14. Disconnect the cuff and tubing from the monitor prior to storing.

Fig. 5

Fig. 6


Fig. 7

INTERPRETING YOUR RESULTS

Irregular Heartbeat Detection ‹♥›

An irregular heartbeat is defined as a heartbeat rhythm that is more than 25% slower or 25% faster than the average measurement detected while the monitor is conducting the measurement. When this device detects an irregular rhythm two or more times during the reading, the IHB symbol will appear on the display screen.

This monitor will operate if the icon appears, but the results may be compromised, especially if this symbol appears often. Please consult with your physician or trained healthcare professional for further information regarding an irregular heartbeat and if this symbol appears frequently.

INTERPRETING YOUR RESULTS

Hypertension Indicator

This unit features our unique Hypertension Indicator. The World Health Organization has established globally accepted standards for the assessment of high or low blood pressure.

This chart can help you identify where your blood pressure measurements fall according to WHO classification standards. This chart should only be considered as a guideline; always consult with your physician to interpret your individual results. Never adjust your medication(s) dosage, or severely alter your diet or exercise routine without consulting with your physician.

The Hypertension Indicator is located on the left side of the display window

Stage 3 Hypertension

Systolic pressure: ≥ 180 mmHg
or
Diastolic pressure: ≥ 110 mmHg

High Normal

Systolic pressure: ≥ 130 -139 mmHg
or
Diastolic pressure: ≥ 85 -89 mmHg

Stage 2 Hypertension

Systolic pressure: ≥ 160 -179 mmHg
or
Diastolic pressure: ≥ 100 -109 mmHg

Normal

Systolic pressure: ≥ 120 -129 mmHg
or
Diastolic pressure: ≥ 80 -84 mmHg

Stage 1 Hypertension

Systolic pressure: ≥ 140 -159 mmHg
or
Diastolic pressure: ≥ 90 -99 mmHg

Optimal

Systolic pressure: < 120 mmHg
or
Diastolic pressure: < 80 mmHg

MEMORY FUNCTION-RECALL

Recalling Measurements in Memory

This monitor stores and recalls up to 60 measurements per memory bank, 120 total, plus an average of the last three measurements in each bank. These measurements can be shared with your physician or trained healthcare professional.

1. Press and release the **M** (MEMORY) button. The display will illuminate and show the last memory bank used, U1 or U1, and the number of measurements in that bank, Fig 1. Press the **POWER** button if you wish to select a different memory bank.
2. Press the M button again to review the saved measurements. The first reading displayed is the average of the last three measurements in that bank, Fig 2.
3. Press the M button again to scroll through all of the measurements. Each measurement will indicate the date and time taken, the measurement and pulse reading, Hypertension Indicator results and Irregular Heartbeat if present, Fig 3.
4. After you've scrolled through all the readings, the unit will automatically turn off.

When the number of readings exceeds 60, the oldest data will be replaced with the new record.


Fig. 1


Fig. 2

Fig. 3

MEMORY FUNCTION – CLEARING READINGS

Clearing Measurements from Memory

1. Press and release the **M** (MEMORY) button. The unit will display how many measurements are stored within the bank.
2. Select the desired user bank to clear by pressing the **⏻ POWER** button to switch between U1 and U2.
3. Press and release the **M** (MEMORY) button two times to advance to the first measurement in the bank.
4. When the first measurement is on the display (indicated by the 01 at the lower left of the screen), Fig. 1, press and hold the **M** (MEMORY) button until two horizontal lines flash, Fig. 2. This indicates that all measurements have been erased.


Fig. 1


Fig. 2

REMINDER: WHEN USING AN AC ADAPTER, IT IS RECOMMENDED THAT BATTERIES ARE ALSO INSTALLED IN ORDER TO PRESERVE THE DATE AND TIME SETTINGS WHEN THE MONITOR IS NOT PLUGGED IN TO A POWER SOURCE.

CARE & MAINTENANCE

Proper care and maintenance of your blood pressure monitor will help ensure long life for the product; improper use of the monitor may void the included warranty.

Always disconnect the cuff and tubing from the monitor prior to cleaning and storage.

Monitor Care and Maintenance

- Do not use any liquids on the monitor; use a soft, dry cloth to clean the monitor as needed
- Storage outside of stated storage temperature may result in measurement error or device malfunction; storage environment temperature is: -4°F – 131°F (-20°C – 55°C); humidity: ≤90% RH
- Do not store in direct sunlight
- Always keep the unit unplugged while not in use; remove the batteries if the unit will not be used for an extended period of time
- Do not attempt to disassemble your monitor; disassembling your monitor may void the manufacturers warranty
- Do not subject the monitor to strong shocks; take care not to drop the monitor
- With proper care, the monitor will maintain performance for a minimum of 10,000 measurements or three years
- No component of this device can be maintained by the user; the circuit diagrams, component part lists, descriptions, calibration instruction or other information which will assist the user's appropriately qualified technical personnel to repair those parts of equipment which are designated repairable can be supplied
- It is recommended that monitor performance be reviewed every two years or after service or repair; please contact your service provider or call Customer Care

CARE & MAINTENANCE

Cuff Care and Maintenance

- The cuff should not be submerged or washed in a washing machine, but may be spot cleaned with a mild-detergent as needed
- Clean the cuff thoroughly after every 200 uses or more frequently as needed
- To disinfect the cuff, wipe the inside (skin contact side) of the cuff lightly moistened with 70-90% ethyl rubbing alcohol and allow to air dry; it is suggested to disinfect the cuff frequently when used in clinical settings or if the cuff is contaminated in any way
- Never attempt to iron or dry the cuff; if the cuff becomes wet, allow the cuff to air dry completely before use or storage
- With proper care, the cuff will maintain performance for a minimum of 1000 measurements

DEVICE & LABEL SYMBOLS

These symbols may appear on your device, instructions or packaging and may vary by make and model.

Symbol	Meaning
	Read This Manual—All included manuals should be read prior to first use
	Warning—Symbol indicates a warning, prohibition or mandatory action that mitigates a risk that is not necessarily obvious to the device operator
	Type BF Applied Parts—Indicates that a part of this unit comes in contact with the patient in order to carry out it's intended function; in the case of this device the cuff is the Type BF Applied Part
	Environment Protection—dispose of this product properly; consult with your local recycling ordinances for proper recycling and disposal
	Keep Dry—This device should be kept dry; never submerge the unit or cuff. Consult with the Care and Maintenance section of this manual for information on cleaning your monitor
	Polarity Symbol—Polarity of D.C. power connector; note connector may not be included with this device but is available for purchase with the included form

DISPLAY SYMBOLS

Symbol	Meaning
SYS	Systolic Blood Pressure– The 'top' blood pressure result
DIA	Diastolic Blood Pressure– The 'bottom' blood pressure result
mmHg	Milligrams of Mercury– Unit of blood pressure measurement for this device
Pulse/ min	Pulse Rate– Pulse rate per minute during the measurement
	Low Battery– Battery life indicator
	Pulse Rate– Pulse rate per minute
	Ready to Inflate Indicator– The cuff is not properly wrapped or has too much air left in it. The arrow will stop flashing when the unit is ready to use
	Measurement Error– A measurement error has occurred and the measurement should be attempted again.
	Irregular Heartbeat Icon– An irregular heartbeat was detected during the measurement
	Memory– The measurement displayed is from the memory bank
	Average– The measurement displayed is average of the last three measurements taken
	Error– An error has occurred during the measurement
	User– Indicates which user the measurement will be stored for

ERROR CODES & TROUBLESHOOTING

Problem	Possible Cause	Solution
Display shows an abnormal result	Cuff position was not correct or was not tight enough	Apply the cuff properly and attempt a new measurement
	Body posture was not correct during measurement	Review the 'Taking a Measurement' portion of this manual and attempt a new measurement
	Speaking, excessive movement, angry, excited or nervous anxiety during measurement	Wait a period of time and attempt a new measurement after calm has returned; do not speak or move during measurement
Display is blank when power is on	Batteries may be expired or installed improperly	Check the polarity of the batteries and reinstall if necessary; replace batteries with new batteries
Display Shows Er 0	Pressure system is unstable before measurement	Review the 'Taking a Measurement' portion of the manual, wait a minimum of 5 minutes and attempt a new measurement.
Display Shows Er 1	Failed to detect systolic pressure	
Display Shows Er 2	Failed to detect diastolic pressure	
Display Shows Er 3	Pneumatic system blocked or cuff is too tight during inflation	Review the Applying Cuff portion of this manual; re-apply the cuff, wait a minimum of 5 minutes and attempt a new measurement.
Display Shows Er 4	Pneumatic system leakage or cuff is too loose during inflation	
Display Shows Er 5	Cuff pressure above 300 mmHg	Check that the cuff is applied properly, wait a minimum of 5 minutes and attempt a new measurement. If the monitor still does not perform or performs abnormally, contact Customer Care
Display Shows Er 6	Cuff pressure above 15 mmHg for more than 3 minutes	
Display Shows Er 7	EEPROM accessing error	
Display Shows Er 8	Device parameter checking error	
Display Shows Er A	Pressure sensor parameter error	

If your error persists or is not remedied, please contact Customer Care at 866-326-1313

FCC STATEMENT

NOTE:

POTENTIAL FOR RADIO/TELEVISION INTERFERENCE (for U.S.A. only)

This product has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules.

These limits are designed to provide reasonable protection against harmful interference in a residential installation. The product generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If the product does cause harmful interference to radio or television reception, which can be determined by turning the product on and off, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna
- Increase the separation between the product and the receiver.
- Connect the product into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

POTENTIAL FOR RADIO/TELEVISION INTERFERENCE (for Canada only)

This digital apparatus does not exceed the Class B limits for radio noise emissions from digital apparatus as set out in the interference-causing equipment standard entitled "Digital Apparatus", ICES-003 of the Canadian Department of Communications.

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Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

ELECTROMAGNETIC COMPATIBILITY

Table 1
For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacturer's declaration - electromagnetic emissions

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

<i>Emission Test</i>	<i>Compliance</i>	<i>Electromagnetic Environment Guidance</i>
RF emissions CISPR 11	Group 1	This monitor uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. This monitor is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Not Applicable	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not Applicable	

ELECTROMAGNETIC COMPATIBILITY

Table 2
For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacture's declaration - electromagnetic immunity

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

<i>Immunity Test</i>	<i>IEC 60601 Test Level</i>	<i>Compliance Level</i>	<i>Electromagnetic Environment Guidance</i>
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 s	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of this monitor requires continued operation during power mains interruptions, it is recommended that this monitor be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

ELECTROMAGNETIC COMPATIBILITY

Table 3

For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING

Guidance and manufacturer's declaration - electromagnetic immunity

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

<i>Immunity Test</i>	<i>IEC 60601 Test Level</i>	<i>Compliance Level</i>	<i>Electromagnetic Environment Guidance</i>
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of this monitor, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2.3 \sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,* should be less than the compliance level in each frequency range.** Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

* Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which this monitor is used exceeds the applicable RF compliance level above, this monitor should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating this monitor.

** Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

ELECTROMAGNETIC COMPATIBILITY

Table 4

For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and this monitor

This monitor is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this monitor can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this monitor as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter <i>W</i>	Separation distance according to frequency of transmitter <i>m</i>		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2.3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance <i>d</i> in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where <i>P</i> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

PRODUCT SPECIFICATIONS

NAME	Jumbo Screen Premium Digital Blood Pressure Arm Monitor
MODEL NUMBER	01-514
DISPLAY SYSTEM	Digital display/LCD
MEASUREMENT METHOD	Oscillometric
POWER SOURCE	4 x 1.5V  SIZE AA, DC 6V 600mA
BATTERY LIFE	Approximately 270 measurements
MEASURING RANGE	Diastolic: 40-199 mmHg; Systolic: 60-260 mmHg Pulse: 40-180 beats/minute
ACCURACY	Pressure: ± 3 mmHg Pulse: $\pm 5\%$ of reading
INFLATION	Pump Drive
DEFLATION	Automatic Pressure Release Valve
MEMORY	Built-in memory enabling display of up to 60 measurements per bank; 120 total measurements
AUTOMATIC SHUT-OFF	Approximately 1 minute after last button operation
OPERATION ENVIRONMENT	Temperature 41°F – 104°F (5°C – 40°C) Humidity $\leq 90\%$
STORAGE ENVIRONMENT	Temperature -4°F – 131°F (-20°C – 55°C) Humidity $\leq 90\%$
MONITOR DIMENSIONS	4-1/8" x 6-3/4" x 2-3/4"
WEIGHT	10.4 oz. (without batteries)
CUFF ARM CIRCUMFERENCE	Adult size cuff fits arm circumference: 11-3/4" - 16-1/2" (30-42 cm)
ACCESSORIES	Cuff, detailed guidebook, quick start guide, blood pressure log, four "AA" batteries, storage case
OPTIONS	Small adult size cuff fits arm circumference: 8-5/8" - 11-3/4" (22-30 cm), Large adult size cuff fits arm circumference: 16-1/4" - 18-7/8" (42-48 cm)
(AVAILABLE SEPARATELY)	AC Adapter

Specifications are subject to change without notice

LIFETIME LIMITED WARRANTY

Congratulations on your purchase of a Digital Blood Pressure Monitor. Your Digital Blood Pressure Monitor is covered by the following limited warranty commencing upon the date of purchase, and subject to the following terms and conditions:

The warrantor warrants that its Digital Blood Pressure Monitor will be free from defects in materials and workmanship under normal consumer usage for as long as the original purchaser owns the product.

Periodic maintenance, repair and replacement of parts due to normal wear and tear are excluded from coverage. Defects or damage that result from: (a) improper operation, storage, misuse or abuse, accident or neglect, such as physical damage (cracks, scratches, etc.) to the surface of the product resulting from misuse; (b) contact with liquid, water, rain, extreme humidity or heavy perspiration, sand, dirt or the like, extreme heat, or food; (c) use of the Digital Blood Pressure Monitor for commercial purposes or subjecting the Digital Blood Pressure Monitor to abnormal usage or conditions; or (d) other acts which are not the fault of Veridian, are excluded from coverage. This warranty does not cover batteries or other power sources that may be provided with, or used with the Digital Blood Pressure Monitor.

If the Digital Blood Pressure Monitor fails to conform to this limited warranty, return the Digital Blood Pressure Monitor postage prepaid to: **Veridian Healthcare, Attn: Repair Department, 1175 Lakeside Drive, Gurnee, IL, 60031**. When returning a product, please also include: (i) a copy of your receipt, bill of sale or other comparable proof of purchase; (ii) a written description of the problem; and (iii) your name, address and telephone number. Carefully package the product to avoid any damage that may occur while in transit; shipping insurance with returned receipt is recommended. At our option, the warrantor will repair or replace the unit found to be defective in materials or workmanship under normal consumer usage. The purchaser will be notified of any additional repairs required prior to completing the repair, and will be responsible for parts charges, if any, and repair charges not covered by this limited warranty.

EXCEPT AS PROVIDED FOR IN THIS LIMITED WARRANTY, ALL EXPRESS AND IMPLIED WARRANTIES AND CONDITIONS ARE DISCLAIMED, INCLUDING WITHOUT LIMITATION THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. THE REPAIR OR REPLACEMENT AS PROVIDED UNDER THIS LIMITED WARRANTY IS THE EXCLUSIVE REMEDY OF THE CONSUMER, AND IS PROVIDED IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. IN NO EVENT SHALL THE WARRANTOR BE LIABLE, WHETHER IN CONTRACT OR TORT (INCLUDING NEGLIGENCE) FOR DAMAGES IN EXCESS OF THE PURCHASE PRICE OF THE PRODUCT, OR FOR ANY INDIRECT, INCIDENTAL, SPECIAL OR CONSEQUENTIAL DAMAGES OF ANY KIND, OR FOR DAMAGES TO, OR LOSS OF, OTHER PROPERTY OR EQUIPMENT OR PERSONAL INJURIES TO THE FULL EXTENT THESE DAMAGES MAY BE DISCLAIMED BY LAW.

Some states and jurisdictions do not allow the limitation or exclusion of incidental or consequential damages, or limitation on the length of an implied warranty, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that vary from state to state or from one jurisdiction to another.

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