visocor[®] HM40

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Instructions for Use (30-57)

Mode d'emploi (58-85)

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CE 0123

Inhaltsangabe

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Contents

A	Product information/notes on use	32
	Device description/indicators	34
B	Abridged instructions on measurement	36
C	Operation of the device	
	1. Setting up the device	39
	2. Fitting the cuff	39
	3. Positioning the device correctly	40
	4. Inflating the cuff	42
	5. Measuring blood pressure	43
	6. Memory recall	44
	7. Deleting the memory	44
D	Important information	
	1. Basic rules for blood pressure measurement	45
	2. Other points to note	46
	Wrist circumference	46
	Cardiac arrhythmia	46
	Frequency of cardiac pacemakers	46
	Pregnancy	46
	Medicines	46

E	What you should know about blood pressure	
1.	Systolic and diastolic blood pressures	47
2.	Why you have to measure different parameters	47
3.	Why measure regularly?	47
4.	What are normal blood pressure levels?	48
F	Maintenance	
1.	Changing the batteries	48
2.	Cleaning and disinfection	49
3.	Customer service	49
G	Appendix	
1.	What to do in the event of measuring errors and operating errors	50
2.	Specifications/accessories/spare parts	52
3.	Explanation of symbols	54
4.	Metrological inspection	55
5.	Warranty	56

These instructions are designed to help the user operate the digital -visocor HM40 blood pressure monitor safely and efficiently. The device must be used according to the procedures contained in these instructions and must not be applied/used for any other purposes. It is important to read and understand the entire set of instructions. Read through the instructions carefully and familiarise yourself particularly with the section entitled „Basic rules for blood pressure measurement“.

Measuring accuracy can be incorrect if absolute rest is not observed during measurement or if the conditions indicated (temperature or humidity) are outside the limits set out in the specifications. Inaccurate results can occur if the blood pressure monitor is used on children. Consult your doctor if you wish to measure the blood pressure of a child. Under no circumstances may the device be used on a baby.

The standard cuff fits an wrist circumference of 12,5 to 21,5 centimetres. The device covers pressures up to 300 mmHg and a pulse rate of 40 to 160 beats per minute.

Method of measurement

This product uses the oscillometric method for measuring blood pressure and pulse rate. When the Start/Stop button is pressed, the system starts to build up pressure automatically up to ca. 190 mmHg. If cuff pressure is too low, the device automatically increases the pressure in increments of 40 mmHg.

Within a very short space of time the device records the small oscillations within the cuff, which occur owing to expansion and contraction of the arteries in the arm (heart beat). The measuring process stops after establishing the systole, diastole and pulse when air is let out of the cuff.

The amplitude of each pressure wave is measured in millimetres of mercury (mmHg), converted and displayed on the LCD in the form of a digital value. A memory stores the last 30 results for the purposes of comparison.

The visocor HM40 is based on non-invasive measurement of systolic and diastolic blood pressure on the wrist and measurement of pulse rate in adults, i.e. from the age of 15 upwards. The product is recommended to patients with unstable blood pressure or known hypertension, for home use and therapy support.

Device description

A

- ☒ Display for systole, diastole, pulse and control readings
- ☒ Start/Stop button
- ☒ Memory button (memory recall button)
- ☒ Battery compartment
- ☒ Cuff

Kontrollanzeigen/Symbole

No.30	Memory data (No. of measurement)
A 13	Average No. of stored data
	Pulse rate mark
	Battery exchange
SYS	Systole
DIA	Diastole
PUL	Pulse rate
	Inflation / Deflation

visocor HM40:

Automatic digital device with electric pump for measuring blood pressure on the wrist.

Please read the detailed Instructions for Use before measuring for the first time (from page 39 onwards).

1. Insert batteries

(2 x 1.5 V alkaline-manganese batteries)

2. Fit the cuff

Bare your left wrist.

Bare your left wrist, remove any watch or jewellery. Place the device round the wrist. The distance between the top edge of the cuff and the palm should be approx. 1 cm (Fig. 1). Close the cuff. Keep the device at heart level (Fig. 2). Keep your arm still and do not talk.


Fig. 1



Fig. 2

3. Inflate the cuff

Press the Start/Stop button. The device conducts a performance test and all the display elements light up.

Inflation of the cuff is automatic up to approx. 190 mmHg. If cuff pressure is too low, the device automatically increases the pressure in increments of approx. 40 mmHg.

4. Blood pressure measurement

After inflation the air is let out automatically, the falling cuff pressure is displayed and the „☒“ symbol flashes. Automatic deflation of the cuff indicates the end of measurement.

The systolic and diastolic blood pressure readings will appear in the display together with your pulse rate.

5. Switch off the device

Press Start/Stop button. The figures measured are saved automatically.

1. Setting up the device

Place the enclosed batteries in the device if they have not already been inserted (see page 48, Changing batteries).


Important:

Never use the device above a critical part of the body (e.g. wound, intravascular port).

2. Fitting the cuff

Bare your left wrist, remove any watch or jewellery. Place the device round the wrist. The distance between the top edge of the cuff and the palm should be approx. 1 cm.


Left-handed persons can place the cuff round their right arm. The instructions for use apply by analogy. Now pull the free end of the cuff tight and place it over the velcro strip firmly. Owing to its optimised shape the cuff sits tightly in order to prevent any measuring errors.


3. Positioning the device correctly

3 measuring positions are possible:

1. Place the left arm on a dining table and place it at an angle in such a way that the wrist is exactly at heart level.


2. If measurement is not being performed at a table, the following position can be used: press the Start button and place the left arm loosely against the body in such a way that the cuff of the device is at heart level. To assist and relieve the left arm you can use the right arm as a support. Keep both arms still until measurement has been completed.

3. Measuring in the lying position

Lie down on your back. Place the device round your wrist and keep it at heart level, supported by a cushion for example. To relax before measurement, take a deep breath 5 or 6 times. During measurement it is important to keep your arm still and refrain from talking.

If the cuff is lower or higher than heart level, the measured value will tend to be higher or lower as the case may be.


Do not bend your wrist! Correct measurements can only be obtained if the wrist is at heart level. If the wrist is too high or too low, it will lead to false results.

4. Inflating the cuff

Switch on the device by pressing the Start/Stop button. The entire display lights up (segment test). In this way you can check to make sure that all the display elements are working properly. Shortly after that a down arrow appears followed by the display „0“, and then the cuff is automatically inflated to approx. 190 mmHg. If that cuff pressure is not sufficient for measuring blood pressure levels, the device automatically increases cuff pressure in increments of 40 mmHg until sufficient pressure has been reached. The inflation and measuring process can be discontinued at any time by pressing the Start/Stop button. To control the inflation processed manually, keep the Start/Stop button pressed as soon as inflation has started and release it again as soon as cuff pressure has reached 40 mmHg above the maximum systolic pressure to be expected. If cuff pressure has exceeded approx. 190 mmHg,

inflation stops as soon as you let go of the button. The pressure can be increased up to approx. 280 mmHg.

Please note:

The device may only be inflated when it is attached to the wrist.

Operation of the device

5. Measuring blood pressure

As soon as the necessary cuff pressure has been reached, the pump switches off and the air is let out automatically. The falling cuff pressure appears on the display. After a short while the „☹“ symbol flashes at the same rate as your pulse. At the end of measurement air is let out of the cuff automatically. The blood pressure levels measured and the pulse reading appear on the display. The systole and diastole levels are saved automatically when the device is switched off.

Record the figures saved in your blood pressure pass. Over a lengthy period this will provide you with a better picture of your blood pressure situation, which you may wish to discuss with your doctor.

Before commencing a repeat measurement you should wait three to five minutes until the circulation of blood in your arm has returned to normal.

Excessively frequent measurements have an effect on blood flow and should preferably be restricted to the morning and evening measurements.

Important:

Irregular, double or absent pulse signals may be an indication of cardiac arrhythmia.

6. Memory recall

The device has a memory which makes it possible to display the last 30 blood pressure levels and pulse values measured and to indicate the mean (average). If you press the Memory button, the mean value will appear provided 2 measurements have already been saved. The mean is indicated by „A“ and the number of results saved.

Pressing again calls up the last value measured. Pressing again displays the previous measurement as the case may be. In this way you can call up all 30 measurements backwards. If more than 30 measurements have been saved, the oldest measurement (no. 30) is deleted in order to record the latest measurement (no. 1). If the button is pressed again, the mean of all the measurements appears.

7. Deleting the memory

To delete individual measurements, press the Memory button a number of times to call up the measurement required, press again and keep the button pressed for 8 - 10 seconds until the measurement has disappeared. To delete the entire memory, call up the mean, press the Memory button again and keep it pressed until the mean has disappeared (8-10 sec).

If the display flashes, the deleting process has not been completed and must be repeated.

Important information

1. Basic rules for blood pressure measurement
 - Measurement must be performed when you are at rest. Blood pressure is increased by walking and climbing stairs. Wait at least 5-10 minutes until your circulation has calmed down.
 - Check your blood pressure twice a day – in the morning after getting up, and in the evening when you have relaxed after work.
 - Blood pressure must be measured when you are seated. For measurement, sit down at your dining table (never at a low coffee table) in as relaxed and upright a posture as possible, support the arm and back, and do not cross your legs. Keep your wrist at heart level. Do not bend your wrist.
 - The most reliable way of obtaining accurate measurements is to take off your watch or jewellery and perform measurement on your bare wrist.
 - During measurement the arm with the cuff must remain still in order to prevent interfering pulses. Do not talk during measurement.
 - Between successive measurements always take a break of about 3-5 minutes in order to prevent incorrect measurements.

2. Other points to note

- Wrist circumference

The cuff is designed for a wrist circumference of 12,5 – 21,5 cm.

- Cardiac arrhythmia

If you suffer from cardiac arrhythmia, this can lead to incorrect results. (During measurement, make sure the heart symbol on the display flashes at a constant rate). Therefore consult your doctor before measuring your blood pressure.

- Frequency of cardiac pacemakers

The display of pulse rate is not suitable for monitoring the frequency of cardiac pacemakers.

- Pregnancy

Blood pressure levels can be different as a result of pregnancy. Only perform blood pressure self-monitoring after consulting your doctor.

- Medicines

Never use this device without consulting your doctor first if you are undergoing dialysis treatment or taking anti-coagulant, anti-platelet medicines or steroids. The use of a blood pressure monitor under such conditions can cause internal bleeding (possibly causing haematomas).

Self-measurement does not constitute treatment. Never change of your own accord the dose of medicines prescribed by your doctor.

What you should know about blood pressure

1. Systolic and diastolic blood pressures

The cardiovascular system has the important function of supplying all the organs and tissues in the body with sufficient blood and removing metabolic products. The heart contracts and relaxes at a regular rate of about 60-80 times a minute. The pressure of the flowing blood on the arterial walls, which arises when the heart contracts, is termed systolic. The pressure in the ensuing relaxation phase when the heart is refilled with blood is termed diastolic. When you perform daily monitoring you measure both figures.

2. Why you have to measure different parameters

Our blood pressure responds to external and internal influences like a sensitive measuring instrument. It can be affected by even slight changes. That makes it easier to understand

why figures measured at the doctor's or pharmacist's are often higher than those you measure at home in your accustomed environment. Changes in weather, changes in climate, physical work and emotional stress can also have an effect.

3. Why measure blood pressure regularly?

The time of day also has an influence on the level of blood pressure. During the day the figures are usually higher than during periods of rest at night. For this reason, one-off measurements and irregular measurements will tell you little about your actual blood pressure. Reliable assessment is only possible if individual measurements are taken regularly. Discuss the results with your doctor.

4. What are normal blood pressure levels

The World Health Organization (WHO) has compiled the following overview for assessing blood pressure levels:

Normal pressure	Threshold	High pressure
Systolic pressure ≤ 119	Systolic pressure 120–139	Systolic pressure $\geq 140^*$ and/or **
Diastolic pressure ≤ 79	Diastolic pressure 80–89	Diastolic pressure ≥ 90

* mmHg = millimetres of hydragyrum (mercury)

** It is sufficient if one of the figures is too high

1. Changing batteries

- Open the battery compartment
Remove the cover of the battery compartment on the side of the device.
- Insert batteries
Take the old batteries out of the device and insert the new ones. Make sure the polarity is correct (there is a mark inside the battery compartment). The easiest way of inserting the batteries is to push them into the battery compartment while pressing their negative (–) terminal against the spring.
- Close the battery compartment
Close the battery compartment by replacing the cover.
- Please note:
Change the batteries
 - if the battery symbol on the display appears after the display segments have been tested
 - if the display fails to light up despite your having pressed the Start button.

Always change both batteries. Only use branded batteries of the quality specified (see Specifications). If the device is going to be out of use for any lengthy period, please remove the batteries. This prevents leakage. Any battery can leak.

Caution!

Do not use rechargeable batteries.

Batteries and appliances should not be disposed of with domestic waste but should be handed in at appropriate collection and disposal points.

2. Cleaning and disinfection

Only clean the device with a soft, slightly moistened cloth and a small amount of detergent. Stains on the cuff can be removed carefully using a moist cloth and soap and water.

If you have to disinfect the device, we recommend using Descosept AF, Dr. Schumacher GmbH, Melsungen. It can be used for fast wipe-down and spray disinfection.

The cuff must not be washed or dry-cleaned.

3. Customer service

The device may only be repaired by the manufacturer or an authorised service centre.

Please contact:

UEBE Medical GmbH

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97877 Wertheim, Germany

Tel.-Nr.: +49 (0) 9342/924040

Fax-Nr.: +49 (0) 9342/924080

E-mail: info@uebe.com

Internet: www.uebe.com

What to do in the event of measuring errors or operating errors

People suffering from cardiac arrhythmia, vascular constrictions, arteriosclerosis in the extremities and diabetics should consult their doctor before performing blood pressure self-monitoring because in such cases the blood pressure levels may be different.

Error	Cause	Remedy
ERR 1: Measuring error	Arm or body movement during measurement. The cuff is not sitting properly. Talking during measurement.	Repeat measurement after 3-5 minutes. Do not move your arm and do not talk. If the error occurs again, contact customer service.
ERR 2: Inflation error	No pressure is being built up in the system.	Repeat measurement and if the error occurs again, call the service telephone number.
ERR 3: Deflation error	Deflation rate too high or too low.	If the error occurs again, call the service telephone number.
	A fault has been detected.	Repeat measurement after 3-5 minutes. Do not move your arm and do not talk.
ERR 300	Inflation pressure above 300 mmHg	Repeat measurement after 3-5 minutes.
	Battery error	Check batteries
No readings in the display fields.	No batteries inserted. Batteries not inserted properly. Batteries empty.	Check batteries and replace if necessary.
Battery symbol 	Batteries empty.	Check batteries and replace if necessary.

Error	Cause	Remedy
Incorrect measured values	Arm or body movement during measurement.	Repeat measurement after 3-5 minutes. Do not move your arm and do not talk.
	Wrist position too high or too low.	Repeat measurement after 3-5 minutes, with your wrist at heart level.
	Cuff not sitting properly Cuff not fitted properly	Check wrist circumference Repeat measurement after 5 minutes
	Cardiac arrhythmia during measurement	Repeat measurement after 5 minutes
Measurements too high	5-10 minutes' relaxation not observed before measurement	Repeat measurement with at least 5 minutes' relaxation before the next measurement
Cuff is not being inflated	Batteries empty?	Check batteries and replace if necessary.

2. Specifications/accessories/spare parts

- Type of device
 - visocor HM40 – automatic digital device with electric pump for blood pressure monitoring on the wrist.
- The device bears a serial number  which provides clear identification
- Device dimensions:
ca. B = 80 mm x L = 70 mm x H = 30 mm
- Weight: 115 g not including batteries
- Measuring system:
Oscillometric determination of systole, diastole and pulse
- Reference procedure of clinical trial:
Auscultatory measurement
- Display:
LCD display (liquid crystal display)
- Memory:
30 measurements (saved automatically) and average value
- Pressure display range:
0-300 mmHg
- Measuring range:
Systolic: 50 – 250 mmHg
Diastolic: 40 – 180 mmHg
Pulse measurement: 40 – 160 beats/min
- Tolerance:
Blood pressure monitoring:
conforms to EN 1060 Part 3
Pressure measurement: ± 3 mmHg
Pulse measurement: ± 5
- Cuff:
Doubled preformed cuff, Type VHG4
- Deflation valve:
electronic control valve
- Inflation pressure:
ca. 190 mmHg
Additional inflation in increments of 40 mmHg
- Deflation rate:
On average 5,5 mmHg/sec

- Batteries:
 - battery type:
2 x 1.5 V Micro/AAA LR 03 cells
 - Service life:
More than 600 measurements in 2 years
 - Battery indicator: „“ symbol means change the batteries – automatic shutdown:
approx. 3 minutes after the end of measurement
- Operating conditions
Ambient temperature 10-40°C
Rel. humidity up to 85%
- Storage and transport conditions:
Ambient temperature -5 to +50°C
Rel. humidity up to 85%

Subject to technical change without notice.

Accessories/Spare parts

You can obtain the following genuine spare parts and accessories from specialist dealers:

- Preformed cuff “Type VHG 4”
(for arm circumference 12,5 – 21,5 cm)
Art. no. 2204001

Important:

To prevent incorrect measurements, the device may only be used with genuine type VHG 4 cuffs.

Appendix

3. Explanation of symbols

This blood pressure monitor conforms to Council Directive 93/42/EEC of 14 June 1993 on Medical Devices and bears the mark CE 0123. Devices with CE marks are quality-controlled in accordance with this Directive and provide a high level.

Degree of protection against electric shock: Type BF

Observe Instructions for Use

Batteries and appliances must not be disposed of with domestic waste but should be handed in at the appropriate collection and disposal points.

The requirements of the following standards are met:

- DIN EN 1060-1: 1995
Non-invasive sphygmomanometers:
General Requirements
- DIN EN 1060-3: 1997
Non-invasive sphygmomanometers:
Additional Requirements on electro-mechanical
blood-pressure measuring devices
- EN 1060-4:
2004 Non-invasive sphygmomanometers
Test procedure for establishing accuracy
- DIN EN 60601-1: 1990
+ A1: 1993
+ A2: 1995
1. Medical electrical devices
2. General commitment to safety
- DIN EN 60601-1-2: 2001
Medical electric devices:
Electromagnetic compatibility
- This product conforms to Council Directive 93/42/EEC of 14 June 1993 on Medical Devices.

Manufacturer:
UEBE Medical GmbH
Zum Ottersberg 9
97877 Wertheim
Germany

4. Metrological inspection (previously calibration)

Generally, a metrological inspection is to be recommended at intervals of 2 years. Professional users, however, are under obligation in accordance with the „Medical Devices Operator Ordinance“. The inspection can be performed either by UEBE Medical GmbH, an authority responsible for metrology or authorised maintenance service centres.

Competent authorities and authorised maintenance service centres can obtain a „Test Procedure for Metrological Inspections“ from the manufacturer on request.

Important:

No changes may be made to this device without the manufacturer's permission, e.g. the device must not be opened (except to change the batteries).

5. Warranty

The device was manufactured and tested with the utmost care. If it should nevertheless have defects, we give a warranty on the following terms:

1. During the term of the warranty of 2 years as of the date of purchase we agree to remedy such defects, at our discretion and at our expense, by repair (after the product has been returned to our factory) or in the form of a replacement by a non-defective device.
2. The warranty does not cover normal use of wearing parts or any damage which arises as a result of non-observance of the Instructions for Use, improper handling (e.g. unsuitable power sources, breakage, leaking batteries) and/or dismantling of the device by the purchaser. In addition, the warranty does not cover compensation claims filed with us.
3. Warranty claims can only be filed during the term of the warranty by presenting the receipt. In the event of a warranty claim the device must be returned to the following address, together with the receipt and description of complaint: UEBE Medical GmbH, Service-Center, Zum Schlag 18 97877 Wertheim, Germany.
4. The purchaser's contractual defect rights vis-a-vis the seller in accordance with § 437 of the German Civil Code are not restricted by the warranty.

Important:

In the event of a warranty claim please always include the receipt.

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© Copyright UEBE Medical GmbH