visocor[®] HM50

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Table of contents

A	Intended use	34
B	Safety instructions	35
	1. Important instructions for patients	35
	2. Important technical details	36
C	Operating the unit	38
	1. Unit description	38
	2. Displays	39
	3. Important instructions for use	40
	4. Inserting/replacing batteries	41
	5. Activating and setting the date/time	42
	6. Attaching the unit	44
	7. Body posture when taking measurements	45
	8. Measuring blood pressure	46
	9. Irregular pulse waves and cardiac arrhythmia	48
	10. Using the memory	48
D	What you should know about blood pressure	50
	1. Systolic and diastolic blood pressure values	50
	2. Why you measure different values	50
	3. Why you should measure blood pressure regularly	50
	4. What are normal values for blood pressure?	51

Table of contents

E	Technical information	52
	1. Failure and error messages	52
	2. Customer service	55
	3. Technical data	56
	4. Original spare parts and accessories	58
	5. Technical inspection / calibration check	58
F	General provisions	59
	1. Explanation of symbols	59
	2. Disposal	59
G	Maintaining the unit	60
H	Warranty	61

A Intended use

Thank you for choosing the visocor® HM50 wrist blood pressure monitor (also referred to in the following as the unit).

It's determined for measuring blood pressure and pulse at the wrist, those of people 15 years of age and over. This unit is recommended for patients with unstable blood pressure for domestic self administration blood pressure control and treatment support.

visocor® HM50 uses the oscillometric method for measuring blood pressure and pulse rate.

The double-preformed comfort cuff can be easily attached and provides for a perfect fit during measurement.

Using the pressure variations caused by the pulse at falling cuff pressure, the micro-computer calculates the values systole, diastole, and pulse. Afterwards, the cuff is deflated completely.

A measurement value memory stores the most recent 60 measurement results in each case.

These instructions are intended to help you use the unit reliably and effectively. Keep the instructions together with the unit at all times. The unit must be used in accordance with the methods detailed in these instructions and may not be used for other purposes. Read these instructions carefully before using the unit.

Safety instructions **B**

1. Important instructions for patients

- Blood pressure measurements on children require specialist knowledge! Please consult your doctor if you wish to measure the blood pressure of a child. Do not under any circumstances use the unit on a baby/infant.
- The special double-preformed comfort cuff is designed for a wrist circumference of 12.5-21.5 cm..
- The measurement results of automatically measuring blood pressure monitors may be falsified by pregnancy, irregular heartbeat or arteriosclerosis. Measure your own blood pressure in consultation with your doctor.
- Do not under any circumstances place the cuff on or over any critical point, e.g. wound, aneurysm, etc. Risk of injuries! The supply via an intravascular access (infusion) or other medical monitoring devices might be interrupted.
- Do not use the unit without first consulting your doctor if you are undergoing dialysis treatment or taking anticoagulants, platelet aggregation inhibitors or steroids. Internal bleeding may be caused under these conditions.
- If you have had a mastectomy, do not measure blood pressure on the side of the body where the breast has been removed together with the lymph nodes from the armpit.

B Safety instructions

- The unit contains small parts which could be swallowed by children. It should therefore be kept out of the reach of children at all times.
- Measuring your own blood pressure does not constitute treatment. On no account modify of your own accord the dosage of drugs/medication prescribed by your doctor.
- Please refer to the chapter "Important instructions for use" (Page EN-40) before conducting any measurements yourself.

2. Important technical details

- A consistently good power supply to the unit is essential for trouble-free blood pressure measurements.
 - Please use long-life alkaline batteries only (LR03).
 - When replacing batteries, always replace all the batteries at the same time.
 - You will need 2 x 1.5 V AAA/LR03 batteries. Rechargeable batteries have a voltage of only 1.2 V and are thus not suitable.
 - If the unit is not used for an extended period of time, please remove the batteries. Essentially, all batteries can leak.
- If the unit has been stored at temperatures below 0 °C, let it stand at room temperature for at least 1 hour before use.

Safety instructions B

- The unit must only be operated with original parts. The warranty will be invalidated if the unit is damaged by non-approved accessories!
- Using the unit in the immediate vicinity of mobile phones, microwave appliances or other devices with strong electromagnetic fields may result in impaired functioning. Maintain a minimum distance of 3 m from such devices when using this unit.
- The display of the pulse frequency is not suitable for checking the frequency of cardiac pacemakers. Cardiac pacemakers and blood pressure monitors do not influence each other in their mode of operation.
- Never open or modify the unit or the cuff - this is a medical device and may only be opened by authorised specialists. (Exception: replacing batteries). If the unit has been opened, it must be subjected to a metrological inspection by an authorised institution.
- The cuff may be inflated on the wrist only.
- Please comply with the specified operating conditions for measurement. See Technical data, Page EN-56.
- The inflating and measuring procedure can be stopped by pressing the Start/Stop button. In this case, the unit stops the inflation procedure and deflates the cuff.

C Operating the unit

1. Unit description



- 1** Cuff
- 2** Displays
- 3** Start/Stop button
- 4** Memory recall button
- 5** Battery compartment (Removable battery compartment cover)

Operating the unit C

2. Displays

- 1** Storage location or date/time
- 2** SYS = Systole (upper blood pressure value)
- 3** DIA = Diastole (lower blood pressure value)
- 4** PUL 1/min = Pulse calculated pulse frequency per minute
- 5** Pulse signal display or irregular pulse waves (Page EN-48)
- 6** Battery check display
- 7** Deflating cuff
- 8** Inflating cuff
- 9** Memory code

Failure and error messages (Page EN-52)

C Operating the unit

3. Important instructions for use

For a representative blood pressure value, please observe the following points both before and while taking the measurement:

- Refrain from drinking alcoholic or caffeinated beverages and smoking at least one hour before measuring.
- Please rest at least 5 minutes before measuring. Depending on the previous degree of stress/exertion, this can even require up to one hour.
- Expose the wrist, clothing or a watch/jewellery must not impair the blood flow in or from the arm as this influences the blood pressure at the measurement point.
- Body posture must be relaxed:
 - For this purpose, sit at a table (if possible the height of a dining table, not a coffee table!).
 - Rest your back against the backrest of the chair.
 - Support your forearm so that the cuff is at the same height as your heart.
 - Place your feet on the floor and do not cross your legs.
- It is absolutely essential to sit still and remain quiet. Irregular movements, vibrations (shaking), talking and breathing heavily will affect and may even falsify the measurement. Look out for irregular pulse signals in the display; if necessary, repeat the measurement under more favourable conditions.

Operating the unit C

- Unusual measurement results are possible with all automatically measuring blood pressure measurement systems from time to time. Check yourself: did you comply with the instructions for use above? If necessary, repeat the measurement after allowing the blood circulation in the arm to recuperate briefly; rest for about 3-5 minutes for this purpose.

Tip: Remain seated, your monitor will switch off automatically about 3 minutes after a measurement. Afterwards, we recommend that you repeat the measurement.

In order to rule out differences between left and right and to obtain comparable measurement results, it is important always to take blood-pressure measurements on the same arm. Ask your doctor which arm he/she recommends you take your measurements on.


- Blood pressure is not a fixed value. It may fluctuate upwards or downwards by more than 20 mmHg on patients within a few minutes.

4. Inserting/replacing batteries

- Opening the battery compartment:
Remove the cover of the battery compartment on the side of the unit.
- Inserting batteries:
Remove the old batteries from the unit and insert the new ones. Please observe proper polarity (marking in battery compartment).

C Operating the unit

- Closing the battery compartment:
Close the battery compartment by clipping the battery cover back into the unit.

When the display shows the “Flat battery” symbol for the first time, the unit will still have enough charge for about 30 measurements. Please change the batteries during this period.

5. Activating and setting the date/time

The unit has a date/time function which you can activate where necessary. This function is deactivated when the unit is delivered.


Activating the date/time function

To activate the date/time function, press the Start/Stop button and then immediately, while the full display is visible, the Memory button. This function can be deactivated at any time in the same way.



Operating the unit C

Setting the date/time

The four-digit number representing the year appears in the display first. Press the Memory button to set the year. Press the Start/Stop button to confirm the selected year and proceed to setting the date.

Set the month with the Memory button and confirm this with the Start/Stop button. Set the day, hours and minutes in the same way.

Changing the date/time


To change date and time, take out a battery and wait until the display goes blank. Now you can reactivate the date/time function and reset the values.

C Operating the unit

Deactivating the date/time function

To deactivate the date/time function, press the Start/Stop button and then immediately the Memory button.


When the date/time function is activated, the time is also shown in the display when the unit is switched off.

6. Attaching the unit

The special double-preformed comfort cuff is easy to attach and provides for a perfect fit during measurement. The materials used have been tested for skin compatibility in accordance with ISO 10993.

- Expose your wrist, with the palm of your hand facing upwards.
- Slide the unit over your wrist. Make sure the gap between your wrist and the top edge of the cuff is approx. 1 cm.

Operating the unit C

- Now pull the free end of the cuff tight and close the cuff with the Velcro-type fastener.


Thanks to its optimised shape, the cuff fits very tightly to prevent measurement errors.

7. Body posture when taking measurements

Whenever possible, measurements should be taken while sitting. In exceptional cases it is also permitted to take measurements while lying down.

Measurement while sitting


Sit at a dining table, resting your back against the backrest of the chair, and support your forearm so that the cuff is at exactly the same height as your heart. Place your feet on the floor and do not cross your legs.

C Operating the unit

Measurement while lying down

Lie on your back. Attach the unit to your wrist and hold the unit at the same height as your heart (e.g. by resting your arm on a pillow or cushion). Breathe deeply 5 or 6 times to relax before measurement. While measuring keep your arm still and do not talk.


Please note:

- Do not under any circumstances bend your wrist.
- Keep your hand relaxed, do not clench your hand or form a fist.
- If the cuff is lower or higher than heart height, the reading may in turn be higher or lower.

8. Measuring blood pressure


- Switch on the unit by pressing the Start/Stop button. The following display is visible.
- After the unit has completed its calibration against the ambient air pressure, the cuff is automatically inflated to approx. 190 mmHg. If that cuff pressure is not sufficient

Operating the unit C

for measuring the blood pressure level, the device automatically increases the cuff pressure in steps of 40 mmHg until a sufficient cuff pressure has been reached.

- For an inflation pressure exceeding 190 mmHg the inflation process can also be controlled manually: Keep the Start/Stop button pressed immediately after the inflation has started and release it when the cuff pressure has reached 40 mmHg above the expected maximum systolic blood pressure. The inflation stops with releasing the button.
- The actual measurement takes place with the start of deflation. The "♥" symbol in the display flashes until the pulse frequency is displayed.
- The measurement being completed is indicated by a long beep. The determined values for systole, diastole, and pulse are shown on the display.
- The unit switches off automatically after about 3 minutes.

C Operating the unit

9. Irregular pulse waves and cardiac arrhythmia

If the symbol for irregular pulse waves flashes after a measurement, the unit has recorded irregular pulses during the measurement. This can be set off by an irregular heartbeat (arrhythmia), disturbance caused by movement, talking or even breathing in deeply. The symbol is stored with the respective measurement.

If the symbol appears more frequently, this must be discussed with your doctor! Measurement with irregular pulse waves should therefore be viewed as critical and repeated under more favourable conditions.


10. Using the memory

Measured results are automatically stored in the memory. The memory can store up to 60 results and the average value. When more than 60 measured values have been stored, the oldest value (No. 60) is deleted to allow the latest value (No. 1) to be recorded.

Recalling data

Press the Memory button to call up data. The average value of the stored results is displayed with "A".


Operating the unit C

When the Memory button is pressed again, the value measured last (No. 1) appears.

Press the Memory button repeatedly to call up further measured values.

When the date/time function is activated, the measured values are displayed alternately with memory number, date and time.

The unit switches itself off automatically 30 seconds after the memory is called up.

Deleting data

To delete individual values, press the memory recall button several times to obtain the measured value required.

Then press the Memory button again and hold it down for 8-10 seconds until the value first flashes and then disappears.

To delete the entire memory, recall the relevant average value display (A) and press the Memory button again, holding it down until the mean value first flashes and then disappears.

D What you should know about blood pressure

1. Systolic and diastolic blood pressure values

The cardiovascular system has the important function of supplying all organs and tissues in the body with sufficient amounts of blood and of transporting metabolites. To do so, the heart contracts and expands at a regular rate of about 60 to 80 times per minute. The pressure of the flowing blood on the artery walls caused by the heart contracting is termed systolic. The pressure in the ensuing relaxation phase, when the heart refills with blood, is termed diastolic. During daily measurement you determine both values.

2. Why you measure different values

Our blood pressure responds to internal and external influences like a sensitive measuring instrument. It can be affected by even slight changes. This explains why values measured by a doctor or pharmacist are higher than those measured at home in the environment you are used to. Changes in the weather, climate changes, or physical or psychological stress can have effects as well.

3. Why you should measure blood pressure regularly

Even the time of day has an influence on your blood pressure. During the day the values are generally higher than during the periods of rest at night. One-off and irregular measurements therefore say little about your actual blood pressure. A reliable assessment is possible only when measurements are taken regularly. Discuss the measurement values with your doctor.

D What you should know about blood pressure

4. What are normal values for blood pressure?

The World Health Organisation (WHO) established the following limits in mmHg to classify the blood pressure values.

WHO 2003	Systolic pressure = Upper value mmHg**	Diastolic pressure = Lower value mmHg**
High pressure	from 140*	from 90*
Limit	120 to 139	80 to 89
Normal value	below 120	below 80

* It is sufficient for one of the values to be raised.

** millimetres of mercury

E Technical information

1. Failure and error messages

Failure encountered	Possible cause	Corrective action
Display Err - 300	Excessive cuff pressure. Cuff inflated to maximum due to movement of the arm or body.	Repeat measurement. Do not move arm. Do not talk.
Display Err - 1	Measurement error. Measurement could not be carried out due to movement of the arm or body.	Repeat measurement. Do not move arm. Do not talk.
Display Err - 2	Inflation error, pressure built up too quickly or too slowly. Cuff is not correctly fitted.	- Repeat measurement - Do not move - Do not talk - Check wrist circumference Call the Service hotline if this occurs repeatedly.
Display Err - 3	Deflation rate too high or too low. Movement or talking during measurement.	- Repeat measurement - Do not move - Do not talk Call the Service hotline if this occurs repeatedly.
	Irregular heartbeat, disruptive movements, shaking, wobbling (objects), breathing in deeply, etc.	Repeat measurement after 3-5 minutes rest. Possible factors (depending on the severity of the arrhythmia) affecting the measurement results must be discussed with your doctor.

Technical information E

Failure encountered	Possible cause	Corrective action
Measured values are too high or too low	Too little rest before the measurement. Movement or talking during measurement.	Repeat measurement after a break of approx. 3-5 mins. Observe instructions for use P. EN-40.
	Has the proper cuff been selected?	Check wrist circumference.
	Incorrect measurement position.	Repeat measurement, ensuring correct body posture. (Page EN-45)
Consecutive measurements produce different measured values (readings).	Blood pressure is not a fixed quantity. It may vary in the upward or downward direction within a few minutes by more than 20 mmHg.	-----
Unusual measured values	Movement or talking during measurement, resting time not observed, feet possibly crossed, smoking or drinking coffee.	Check conditions and repeat measurement. Observe instructions for use Page EN-40.

E Technical information

Failure encountered	Possible cause	Corrective action
Display 0 0	Start/Stop button accidentally pressed while battery was being changed.	Switch unit off and then on again with Start/Stop button, if necessary.
Measurement was interrupted and cuff first deflated and then inflated	Unit detects a fault or that diastolic value is too low.	Unit shows no fault. Repeat measurement.
	Movement during measurement.	If necessary, interrupt measurement, measure again after a break of 5 minutes.
No display after unit is switched on	Batteries inserted incorrectly.	Check position of batteries.
	Batteries flat.	Change batteries.
	Battery contacts dirty.	Clean battery contacts.
No display after Memory button is pressed	No measured values stored.	Take a new measurement.
Measurement interrupted	Batteries flat.	Change batteries.
Failure encountered	Possible cause	Corrective action

Technical information E

Failure encountered	Possible cause	Corrective action
Measured values extremely high or low	Incorrect measurement position.	Repeat measurement, ensuring correct body posture. Do not move your arm and do not talk.
Date/time not shown in display	Function accidentally deactivated. Function was not reactivated after battery replacement.	Activate date and time as described on Page EN-42.

2. Customer service

The unit may only be repaired by the manufacturer or by an expressly authorised organisation. Please contact:

UEBE Medical GmbH
 Zum Ottersberg 9
 97877 Wertheim, Germany
 Tel.: +49 (0) 9342/924040
 Fax: +49 (0) 9342/924080
 E-mail: info@uebe.com
 Internet: www.uebe.com

E Technical information

3. Technical data

Unit type:	Digital automatic unit with electric pump for measuring blood pressure at the wrist
Dimensions:	approx. W = 70 mm x H = 70 mm x D = 80 mm with cuff
Weight:	121 g with cuff, not including batteries
Display:	LCD display (liquid crystal display) for measured values and check displays
Memory:	60 measured values (stored automatically) and average value (A)
Measurement procedure:	Oscillometric determination of systole, diastole and pulse
Reference procedure of clinical testing:	Auscultatory measurement
Inflating pressure:	Initially up to approx. 190 mmHg, then in stages of 40 mmHg
Pressure display range:	0-300 mmHg
Measurement range	Systolic: 50-250 mmHg Diastolic: 40-150 mmHg Pulse: 40-160 pulses/min

Technical information E

Error limits:	Blood pressure measurement: corresponds to EN 1060 Part 3 Pressure measurement: ± 3 mmHg Pulse measurement: ± 5 %
Serial number:	In Battery compartment is a serial number SN which provides clear identification.
Power supply:	Battery type: 2 x 1.5 V, AAA, LR 03, durability: over 700 measurements in 2 years
Cuff:	Cuff type H5, wrist circumference 12.5 cm - 21.5 cm
Operating conditions:	Ambient temperature 10 to 40 °C Relative humidity up to 85%, non-condensing Air pressure 700 to 1040 hPa
Storage and transport conditions:	Ambient temperature -20 to +50 °C, relative humidity up to 85 %, non-condensing
Air deflation valve:	Electronically controlled
Automatic switch-off	approx. 3 minutes after end of measurement
Protection class IP:	IP 20: Protected against solid foreign particles with a diameter of more than 12.5 mm, no protection against water.
Expected service life:	5 years
Classification:	Internally powered ME Equipment

E Technical information

4. Original spare parts and accessories

The following original spare parts and accessories are available from specialist dealers:

- Cuff type H5 (for wrist circumference 12.5-21.5 cm)
Part no. 2205001, PZN 09101582

Subject to technical modifications.

Manufacturer: UEBE Medical GmbH
Zum Ottersberg 9
97877 Wertheim
Germany

5. Technical inspection / calibration check

Basically, a metrological inspection is recommended at intervals of 2 years. However, professional users in Germany are obligated to comply with the aforementioned according to "Regulation for Operators of Medical Devices". This can be implemented either by UEBE Medical GmbH, an authority responsible for metrology, or authorised maintenance services. For this, please observe your national provisions. Upon request, responsible authorities or authorised maintenance services receive a "Test instruction for metrological inspection" from the manufacturer. Important: No modifications, e.g. opening the unit (except to replace the batteries) may be made to this unit without the manufacturer's prior permission.

F Explanation of symbols

1. Explanation of symbols

This product complies with the Council Directive 93/42/EC from 5 September 2007 regarding medical devices and bears the mark CE 0123 (TÜV SÜD Product Service GmbH).

Degree of protection against electric shock: TYPE BF

Please observe instructions for use

Storage and transportation conditions
Ambient temperature -20 to +50 °C

Relative humidity up to 85 %, non-condensing

Keep dry

IP20

Protected against solid foreign particles with a diameter of more than 12.5 mm, no protection against water.

Manufacturer

2. Disposal

Batteries and technical appliances must not be disposed of with domestic waste, but should be handed in at the appropriate collection and disposal points.

G Maintaining the unit

- The unit contains sensitive parts and must be protected against strong variations in temperature, air humidity, dust and direct sunlight.
- Make sure when storing the unit that no heavy objects are left lying on the unit.
- The unit is neither impact-resistant nor shock-proof. We recommend that you have the intactness and accuracy of the display checked after heavier falls or impacts.
- The unit is not waterproof.
- If the unit becomes dirty during use, use a soft, dry cloth to clean the unit. Do not use benzene, thinners or other strong solvents.
- Do not scrub or machine-wash the cuff. If the cuff gets dirty during use, use a synthetic cleaner and softly rub the surface.
- If the unit has been opened, it must be subjected to a metrological inspection conducted by an authorised organisation.

Warranty H

The device has been manufactured and tested with great care. However, in the unlikely event of a defect being detected after delivery, we provide warranty in accordance with the following terms and conditions:

During the warranty period of 2 years from the date of purchase we reserve the right either to repair any such defect at our expense or to supply a perfect replacement unit. The cost of returning the unit to our factory shall be borne by the sender. UEBE shall refuse to accept return deliveries that have not been paid for by the sender.

Excluded from the warranty are parts subject to normal wear and tear as well as damage caused by non-compliance with the instructions for use, improper handling (e.g. unsuitable power sources, breakages, leaking batteries) and/or disassembly of the unit by the purchaser. Furthermore, no claims for damages against us are substantiated by the warranty.

Warranty claims can only be advanced in the warranty period and by presenting proof of purchase. In the event of a warranty claim, the unit must be sent to the following address together with the proof of purchase and a description of the complaint:

UEBE Medical GmbH, Service-Center
Zum Ottersberg 9, 97877 Wertheim, Germany

The cost of returning the unit to our factory shall be borne by the sender. UEBE shall refuse to accept return deliveries that have not been paid for by the sender.

In the case of defectiveness of the goods, the contractual rights of the purchaser to claim against the seller in accordance with § 437 German Civil Code are not limited by the warranty.

Please note: In the event of a warranty claim it is essential to attach the proof of purchase.

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© Copyright 2014 UEBE Medical GmbH