
# Welch Allyn OSZ 3 Easy

---

## Directions for Use

Wrist Watch Blood Pressure Monitor  
Handgelenk-Blutdruckmessgerät  
Tensiomètre bracelet-montre  
Misuratore di pressione da polso  
Tensiómetro Reloj de Pulsera  
Monitor de Pressão Arterial de Pulso

**WelchAllyn®**

Advancing Frontline Care™

Copyright 2007 Welch Allyn. All rights are reserved. No one is permitted to reproduce or duplicate, in any form, this manual or any part thereof without permission from Welch Allyn.

Welch Allyn assumes no responsibility for any injury, or for any illegal or improper use of the product, that may result from failure to use this product in accordance with the instructions, cautions, warnings, or indications for use published in this manual.

Welch Allyn® is a registered trademark of Welch Allyn.


Caution: Changes or modifications not expressly approved by Welch Allyn could void the purchaser's authority to operate the equipment.

Mat. Number: OSZ3 Easy 55-A3-127, Ver: A  
Model 7052-40

Welch Allyn GmbH & Co. KG  
Zollerstrasse 2-4  
D-72417 Jungingen, Germany  
Tel. (+49) 7477-9271-0  
Fax (+49) 7477-9271-90  
e-mail: [info@welchallyn.de](mailto:info@welchallyn.de)

The CE mark on this product indicates that it is in conformity with the provisions of Council Directive 93/42/EEC.

European Regulatory Manager, Medical Division, Welch Allyn Ltd.,  
Navan Business Park, Dublin Road, Navan, County Meath,  
Republic of Ireland  
Tel.: 353 46 90 67700  
Fax: 353 46 90 67755

Printed in USA

**WelchAllyn®**  
Advancing Frontline Care™

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# 1 Wrist Watch Blood Pressure Monitor

---

## 1. Introduction

### 1.1. Features of the device

The blood-pressure monitor is a fully automatic, digital blood-pressure measuring device for use on the wrist, which enables very fast and reliable measurement of the systolic and diastolic blood-pressure as well as the pulse frequency by way of the oscillometric method of measuring.

The device offers a very high and clinically tested measurement accuracy and has been designed to provide a maximum of user-friendliness.

Before using, please read through this instruction manual carefully and then keep it in a safe place. For further questions on the subject of blood-pressure and its measurement, please contact your doctor.

**Attention!**

### 1.2. Important information about self-measurement

- Do not forget: self-measurement means **control**, not diagnosis or treatment: Unusual values must always be discussed with your doctor. **Under no circumstances should you alter the dosages of any drugs prescribed by your doctor.**
- The pulse display is **not** suitable for checking the frequency of heart pacemakers!

- In cases of cardiac irregularity (Arrhythmia), measurements made with this instrument should only be evaluated after consultation with the doctor.

**Electromagnetic interference**

The device contains sensitive electronic components (Microcomputer). Therefore, avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave cookers). These can lead to temporary impairment of the measuring accuracy.

## 2. Important information on the subject of blood-pressure and its measurement

### 2.1 How does high/low blood-pressure arise?

The level of blood-pressure is determined in a part of the brain, the circulatory centre, and adapted to the respective situation by way of feedback via the nervous system.

To adjust the blood-pressure, the strength and frequency of the heart (Pulse), as well as the width of circulatory blood vessels is altered. The latter is effected by way of fine muscles in the blood vessel walls.

The level of arterial blood-pressure changes periodically during the heart activity: During the «blood ejection» (Systole) the value is maximal (systolic blood-pressure value), at the end of the heart's «rest period» (Diastole) it is minimal (diastolic blood-pressure value).

The blood-pressure values must lie within certain normal ranges in order to prevent particular diseases.

### 2.2 Which values are normal?

Blood pressure is too high if at rest, the diastolic pressure is above 90 mmHg and/or the systolic blood-pressure is over 160 mmHg. In this case, please consult your doctor immediately.

Long-term values at this level endanger your health due to the associated advancing damage to the blood vessels in your body.

Should the systolic blood-pressure values lie between 140 mmHg and 160 mmHg and/or the diastolic blood-pressure values lie between 90 mmHg and 95 mmHg, please consult your doctor. Furthermore, regular self-checks will be necessary.

With blood-pressure values that are too low, i.e. systolic values under 105 mmHg and/or diastolic values under 60 mmHg, please consult your doctor.

Even with normal blood-pressure values, a regular self-check with your blood-pressure monitor is recommended. In this way you can detect possible changes in your values early and react appropriately.

If you are undergoing medical treatment to control your blood pressure, please keep a record of the level of your blood pressure by carrying out regular self-measurements at specific times of the day. Show these values to your doctor. **Never use the results of your measurements to alter independently the drug doses prescribed by your doctor.**

Table for classifying blood-pressure values (units mmHg) according to World Health Organization:

| Range                              | Systolic Blood-pressure | Diastolic Blood-pressure | Measures                        |
|------------------------------------|-------------------------|--------------------------|---------------------------------|
| Normal                             | 105-129                 | 60-84                    | Self-check                      |
| High normal                        | 130-139                 | 85-89                    | Self-check                      |
| Hypertension Stage 1 (Mild)        | 140-159                 | 90-99                    | Consult your doctor             |
| Hypertension Stage 2 (Moderate)    | 160-179                 | 100-109                  | Consult your doctor immediately |
| Hypertension Stage 3 (Severe)      | 180-209                 | 110-119                  | Consult your doctor immediately |
| Hypertension Stage 4 (Very Severe) | >210                    | >120                     | Consult your doctor immediately |

**Further information**

- If your values are mostly standard under resting conditions but exceptionally high under conditions of physical or psychological stress, it is possible that you are suffering from «labile hypertension». Please consult your doctor if you suspect that this might be the case.
- Correctly measured **diastolic blood-pressure values above 120 mmHg require immediate medical treatment.**

## 2.3. What can be done if regular increased/low values are obtained?

1. Please consult your doctor.
2. Increased blood-pressure values (various forms of hypertension) are associated long- and medium term with considerable risks to health. This concerns the arterial blood vessels of your body, which are endangered due to constriction caused by deposits in the vessel walls (Arteriosclerosis). A deficient supply of blood to important organs (heart, brain, muscles) can be the result. Furthermore, with long-term continuously increased blood-pressure values, the heart will become structurally damaged.
3. There are many different causes of the appearance of high blood pressure. We differentiate between the common primary (essential) hypertension, and secondary hypertension. The latter group can be ascribed to specific organic malfunctions. Please consult your doctor for information about the possible origins of your own increased blood pressure values.
4. There are measures which you can take, not only for reducing a medically established high blood pressure, but also for prevention. These measures are part of your general way of life:

**Eating habits**

- Strive for a normal weight corresponding to your age. Reduce overweight!
- Avoid excessive consumption of common salt.
- Avoid fatty foods.

**Previous illnesses**

- Follow consistently any medical instructions for treating previous illness such as:
  - Diabetes (Diabetes Miletus)
  - Fat metabolism disorder
  - Gout

**Habits**


- Give up smoking completely
- Drink only moderate amounts of alcohol
- Restrict your caffeine consumption (Coffee)

**Physical constitution:**

- After a preliminary medical examination, do regular sport.
- Choose sports which require stamina and avoid those which require strength.
- Avoid reaching the limit of your performance.
- With previous illnesses and/or an age of over 40 years, please consult your doctor before beginning your sporting activities. He will advise you regarding the type and extent of types of sport that are possible for you.

### 3. The various components of the blood-pressure monitor

The illustration shows the blood-pressure monitor, consisting of:


Wrist cuff Type WC1 13.5--19.5 cm, for wrist circumference of 13.5 to 19.5 cm

## 4. Putting the blood-pressure monitor into operation

### 4.1. Inserting the batteries

After you have unpacked your device, first insert the batteries. The battery compartment is located on the left lateral side of the device (see illustration).

1. Remove cover as illustrated.
2. Insert the batteries (2 x size AAA 1.5V), observing the indicated polarity.
3. If the battery warning appears in the display, the batteries are empty and must be replaced by new ones.

**Attention!**

- After the battery warning appears, the device is blocked until the batteries have been replaced.
- Please use «AAA» Long-Life or Alkaline 1.5V Batteries. The use of 1.2V Accumulators is not recommended.
- If the blood-pressure monitor is left unused for long periods, please remove the batteries from the device.

**Functional check:** Press and hold the O/I button to test all the display elements. When functioning correctly all segments appear for about 2 seconds (If the O/I button is pressed and held for up to 3 seconds, the last measured value will be displayed).

## 5. Carrying out a measurement

### 5.1. Before the measurement:

- Avoid eating, smoking, as well as all forms of exertion directly before the measurement. All these factors influence the measurement result. Try and find time to relax by sitting in an armchair in a quiet atmosphere for about ten minutes before the measurement.
- Always measure on the same wrist (normally left).
- Attempt to carry out the measurements regularly at the same time of day, since the blood-pressure changes during the course of the day.

### 5.2. Common sources of error.

**Note:**

Comparable blood-pressure measurements always require **the same conditions!** These are always **quiet conditions**.


- All efforts by the patient to support the arm can increase the blood-pressure. Make sure you are in a comfortable, relaxed position and do not activate any of the muscles in the

measurement arm during the measurement. Use a cushion for support if necessary.

- If the wrist artery lies considerably lower (higher) than the heart, an erroneously higher (lower) blood-pressure will be measured! (Each 15 cm difference in height results in a measurement error of 10 mmHg!)
- A loose cuff causes false measurement values.
- With repeated measurements, blood accumulates in the respective arm, which can lead to false results. Correctly executed blood-pressure measurements should therefore first be repeated after a 5 minute pause or after the arm has been held up in order to allow the accumulated blood to flow away (after at least 3 minutes).

### 5.3. Fitting the cuff

1. Remove all objects and jewellery (e.g. wristwatch) from the wrist in question.
2. Draw the cuff over the wrist.
3. The distance between the cuff and the hand should be approximately 10 mm.
4. Secure the cuff with the Velcro fastener, so that it lies comfortably and not too tight whereby no space should remain between the cuff and the wrist.
5. Lay the arm on a table, with the palm upwards. Support the arm with a rest (cushion), so that the cuff rests at about the same height as the heart. Take care, that the cuff lies free. Remain so for 2 minutes sitting quietly, before beginning with the measurement.


## 5.4. Measuring procedure

After the cuff has been appropriately positioned, the measurement can begin:

1. Press the O/I-button, the pump begins to inflate the cuff. In the display, the increasing cuff-pressure is continually displayed.

2. After reaching the inflation pressure, the pump stops and the pressure slowly falls away. The cuff-pressure (large characters) is displayed during the measurement. When the device has detected the pulse, the heart symbol in the display begins to blink and a beep tone is audible for every pulse beat.

3. When the measurement has been concluded, a long beep tone sounds. The measured systolic and diastolic blood-pressure values as well as the pulse frequency are now displayed.

Example (Fig.): Systole 118, Diastole 73, Pulse 75

The measurement results are displayed, until you switch the device off. If no button is pressed for 5 minutes, the device switches automatically off, to save the batteries.


## 5.5. Discontinuing a measurement

If it is necessary to interrupt a blood pressure measurement for any reason (e.g. the patient feels unwell), the O/I button can be pressed at any time. The device then immediately lowers the cuff-pressure automatically.

## 5.6. Memory - storage and recall of the measurements

The blood-pressure monitor automatically stores the last measurement value. By pressing the O/I button (when the power is off) for up to 3 seconds, the last measurement can be displayed.

### Further information

Measurements should not occur soon after each other, since the results will be falsified. Wait for several minutes in a relaxed position, sitting or lying, before you repeat a measurement.

## 6. Error messages/malfunctions

If an error occurs during a measurement, the measurement is discontinued and a corresponding error code is displayed (Example: Error No. 2).

| Error Number | Possible Cause(s)                                                                                                             |
|--------------|-------------------------------------------------------------------------------------------------------------------------------|
| Error 1      | No pulse has been detected.                                                                                                   |
| Error 2      | Unnatural pressure impulses influence the measurement result.<br>Reason: The arm was moved during the Measurement (Artifact). |

| Error Number | Possible Cause(s)                                                           |
|--------------|-----------------------------------------------------------------------------|
| Error 3      | The inflation of the cuff takes too long. The cuff is not correctly seated. |

### Other possible malfunctions and their elimination

If problems occur when using the device, the following points should be checked and if necessary, the corresponding measures are to be taken:

| Malfunction                                                                                                            | Remedy                                                                                                                                                                                                       |
|------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The display remains empty when the instrument is switched on although the batteries are in place.                      | <ol style="list-style-type: none"> <li>1. Check batteries for correct polarity and , if necessary, insert correctly.</li> <li>2. If the display is unusual, re-insert batteries or exchange them.</li> </ol> |
| The device frequently fails to measure the blood pressure values, or the values measured are too low (too high).       | <ol style="list-style-type: none"> <li>1. Check the positioning of the cuff.</li> <li>2. Measure the blood-pressure again in peace and quiet under observance of the details made in point 5.</li> </ol>     |
| Every measurement results in different values, although the device functions normally and normal values are displayed. | Please pay attention to the following information as well as the points listed under «Common Sources of Error» and repeat the measurement.                                                                   |
| The blood pressure values measured differ from those measured by the doctor.                                           | Record the daily development of the values and consult your doctor about them.                                                                                                                               |

### Further Information

The level of blood-pressure is subject to fluctuations even with healthy people. **Comparable measurements always require the same conditions (Quiet conditions)!**

If, in spite of observing all these factors, the fluctuations are larger than 15 mmHg, and/or you hear irregular pulse tones on several occasions, please consult your doctor.

This device has been subjected to strict clinical tests. The conformity of this product with the requirements of the CE-

guideline 93/42 EWG, council from June 14, 1993 concerning Medical products and Medical-Product-Legislation from August 2, 1994, was established through a conformity check process.

You must consult your dealer if there are technical problems with the blood-pressure instrument. **Never attempt to repair the instrument yourself!**

Any unauthorized opening of the instrument invalidates all guarantee claims!

## 7. Care and maintenance, recalibration

- a. Do not expose the device to extreme temperatures, humidity, dust, or direct sunlight.

- b. The cuff contains a sensitive air-tight bubble. Handle this carefully and avoid all types of straining through twisting or buckling.

- c. Clean the device with a soft, dry cloth. Do not use petrol, thinners, or similar solvent. Spots on the cuff can be removed carefully with a damp cloth and soapsuds. The cuff must not be washed!

- d. Do not drop the instrument or treat it roughly in any way. Avoid strong vibrations.

- e. Never open the device! Otherwise the manufacturer calibration becomes invalid!

### Periodical recalibration

Sensitive measuring devices must from time to time be checked for accuracy. We recommend

a periodic inspection of the static pressure display **every 2 years**.  
Your dealer will provide more extensive information about this.

## 8. Guarantee

The blood-pressure monitor is **guaranteed for 2 years** from date of purchase. This guarantee includes the instrument and the cuff. The guarantee does not apply to damage caused by improper handling, accidents, not following the operating instructions, or alterations made to the instrument by third parties. The guarantee is only valid upon presentation of the guarantee card filled out by the dealer.

## 9. Reference to standards

|                                       |                                                                                                                                                            |
|---------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Device standard:</b>               | Device corresponds to the requirements of the European standard for non-invasive blood- pressure monitor<br>EN1060-1 / 12:95<br>EN 1060-3 / 09:97          |
| <b>Electromagnetic compatibility:</b> | Device fulfils the stipulations of the European standard<br>EN 60601-1-2                                                                                   |
| <b>Clinical testing:</b>              | The clinical performance test was carried out in Germany according to the DIN 58130/1997 procedure N6 (sequential)<br>ANSI/AAMI SPI O. NIBP - requirements |

The stipulations of the EU-Guidelines 93/42/EEG for Medical Products Class IIa have been fulfilled.

## 10. Technical specifications:

|                                     |                                                                |
|-------------------------------------|----------------------------------------------------------------|
| <b>Weight:</b>                      | 150 g (with batteries + cuff)                                  |
| <b>Size:</b>                        | 84 x 82 x 70 mm (including cuff)                               |
| <b>Storage temperature:</b>         | -5 to +50° C                                                   |
| <b>Humidity:</b>                    | 15 to 85% relative humidity maximum                            |
| <b>Operation temperature:</b>       | 10 to 40° C                                                    |
| <b>Display:</b>                     | LCD-Display (Liquid Crystal Display)                           |
| <b>Measuring method:</b>            | oscillometric                                                  |
| <b>Pressure sensor:</b>             | capacitive                                                     |
| <b>Measuring range:</b>             |                                                                |
| <b>SYS/DIA:</b>                     | 30 to 280 mmHg                                                 |
| <b>Pulse:</b>                       | 40 to 200 per minute                                           |
| <b>Cuff pressure display range:</b> | 0-299 mmHg                                                     |
| <b>Memory:</b>                      | storing the last measurement automatically                     |
| <b>Measuring resolution:</b>        | 1 mmHg                                                         |
| <b>Accuracy:</b>                    | Pressure within $\pm 3$ mmHg<br>Pulse $\pm 5\%$ of the reading |
| <b>Power source:</b>                | 2 dry cells (Batteries) UM-4, size AAA 1.5V                    |
| <b>Accessories:</b>                 | storage case                                                   |

Technical and color alterations reserved!

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Printed in U.S.A.