
Welch Allyn OSZ5

Directions for Use

Automatic Blood Pressure Monitor
Automatisches Blutdruckmessgerät
Tensiomètre automatique
Misura pressione automatico
Tensiómetro automático
Monitor de Pressão Arterial Automático

WelchAllyn®

Advancing Frontline Care™

Copyright 2007 Welch Allyn. All rights are reserved. No one is permitted to reproduce or duplicate, in any form, this manual or any part thereof without permission from Welch Allyn.

Welch Allyn assumes no responsibility for any injury, or for any illegal or improper use of the product, that may result from failure to use this product in accordance with the instructions, cautions, warnings, or indications for use published in this manual.

Welch Allyn® is a registered trademark of Welch Allyn.


Caution: Changes or modifications not expressly approved by Welch Allyn could void the purchaser's authority to operate the equipment.

Mat. Number: OSZ5 55-XX-XXX, Ver: A

Welch Allyn GmbH & Co. KG
Zollerstrasse 2-4
D-72417 Jungingen, Germany
Tel. (+49) 7477-9271-0
Fax (+49) 7477-9271-90
e-mail: info@welchallyn.de

The CE mark on this product indicates that it is in conformity with the provisions of Council Directive 93/42/EEC.

European Regulatory Manager, Medical Division, Welch Allyn Ltd.,
Navan Business Park, Dublin Road, Navan, County Meath,
Republic of Ireland
Tel.: 353 46 90 67700
Fax: 353 46 90 67755

Printed in USA

WelchAllyn®
Advancing Frontline Care™

CONTENTS

Semi-Automatic Blood Pressure Monitor	1
Introduction.	1
Features of the device	1
Important information about self-measurement	1
Important information on the subject of blood pressure and its measurement	2
How does high/low blood pressure arise?	2
Which values are normal?	2
What can be done if regular high/low values are obtained?	4
The various components of the blood pressure monitor	6
Operating the blood pressure monitor	7
Inserting the batteries	7
Tube connection.	8
Carrying out a measurement	8
Before the measurement:	8
Common sources of error	9
Fitting the cuff	10
Measuring procedure	11
Inflating further.	13
Discontinuing a measurement	13
Memory - displaying the last measurement.	14
Error messages/malfunctions.	14
Care and maintenance, recalibration.	16
Guarantee	17
Reference to standards	18
Technical specifications.	18
 Halbautomatisches Blutdruckmessgerät.	 21
Einführung	21
Besonderheiten des Gerätes:.	21
Wichtige Hinweise für die Selbstmessung	21

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Introduction

Features of the device

The blood pressure monitor is a fully automatic digital blood pressure measuring device for use on the upper arm. It enables fast and reliable measurement of the systolic and diastolic blood pressure, as well as the pulse rate, by way of the oscillometric method. The device has been clinically verified to provide very accurate measurements and have been designed to be user friendly.

Operational status and cuff pressure are visible for the entire duration of the measurement on the large LCD display.

Before using the device, read through this instruction manual carefully, then keep it in a safe place. For further information on the subject of blood pressure and its measurement, contact your doctor.

Attention!

Important information about self-measurement

- Do not forget that self-measurement means «control», not diagnosis or treatment. Unusual values must always be discussed with your doctor. Under no circumstances should you alter the dosages of any drugs prescribed by your doctor.
- The pulse display is **not** suitable for checking the pulse rate of heart pacemakers.

- In cases of cardiac irregularity (arrhythmia), measurements made with this instrument should be evaluated only after consultation with the doctor.

Electromagnetic interference:

The device contains sensitive electronic components (a microcomputer). Therefore, avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g., mobile telephones, microwave ovens). These can temporarily impair the measuring accuracy.

Important information on the subject of blood pressure and its measurement

How does high/low blood pressure arise?

The level of blood pressure is determined in a part of the brain, the «circulatory center», and adapted to the respective situation by way of feedback via the nervous system. Blood pressure changes when the strength and pulse rate of the heart and the width of circulatory blood vessels are altered. Blood vessel width is affected by the fine muscles in the blood vessel walls.

The level of arterial blood pressure changes periodically with heart activity. During the blood ejection (systole), the value is maximal (systolic blood pressure value); at the end of the heart's rest period (diastole), it is minimal (diastolic blood pressure value). The blood pressure values must lie within certain normal ranges to prevent particular diseases.

Which values are normal?

While there is no «normal» blood pressure or one combination of systolic or diastolic pressures which is right for everyone, it is generally accepted by physicians that 120/80 is a healthy blood pressure. During the day, the blood pressure fluctuates due to

changes in activity and stress levels. These fluctuations are normal.

Blood pressure is too high if, at rest, the diastolic pressure is higher than 90 mmHg and/or the systolic blood pressure is higher than 140 mmHg. Blood pressure is too low if, at rest, the diastolic pressure is lower than 60 mmHg and/or the systolic blood pressure is lower than 105 mmHg. Long-term values at these levels may endanger your health because of the continuing damage to the blood vessels in your body. However, only a doctor can correctly interpret your blood pressure.

Even with normal blood pressure values, a regular self-check with your blood pressure monitor is recommended. By doing this you can detect possible changes in your values early and react appropriately.

If you are undergoing medical treatment to control your blood pressure, keep a record of the level of your blood pressure by carrying out regular self-measurements at specific times of the day. Show these values to your doctor. **Never use the results of your measurements to alter the drug doses yourself that were prescribed by your doctor.**

Table for classifying blood pressure values (units mm Hg):

Range	Systolic Blood-pressure	Diastolic Blood-pressure	Measures
Normal	105-129	60-84	Self-check
High normal	130-139	85-89	Self-check
Hypertension Stage 1 (Mild)	140-159	90-99	Consult your doctor
Hypertension Stage 2 (Moderate)	160-179	100-109	Consult your doctor immediately
Hypertension Stage 3 (Severe)	180-209	110-119	Consult your doctor immediately
Hypertension Stage 4 (Very Severe)	>210	>120	Consult your doctor immediately

Further information

- If your values are mostly standard under resting conditions but exceptionally high under conditions of physical or psychological stress, it is possible that you are suffering from «labile hypertension». Consult your doctor if you suspect this might be the case.

What can be done if regular high/low values are obtained?

1. Consult your doctor.

Hypertension is one of the major modifiable risk factors for stroke, coronary heart disease, congestive heart failure, renal failure, and peripheral vascular disease.¹ An estimated 50 million people have a blood pressure level consistently at or above 140/90 mm Hg, and therefore, are at increased risk for developing complications.² The risk from hypertension is present whether only systolic pressure, only diastolic pressure, or both are elevated.³

There are many different causes of high blood pressure. Consult your doctor for information about the possible origins of your own increased blood pressure values.

2. You can take measures to prevent medically established high blood pressure as well as to reduce it. These measures are part of your general way of life.

Eating habits

- Strive for a normal weight for your age.
- Avoid excessive consumption of common salt.
- Avoid fatty foods.

Previous illnesses

- Consistently follow any medical instructions for treating such previous illness as:
 - Diabetes (diabetes mellitus).

- Fat metabolism disorder.
- Gout.

Habits

- Give up smoking completely.
- Drink only moderate amounts of alcohol.
- Restrict your caffeine consumption.

Physical constitution

- After a preliminary medical examination, participate in regular physical activity.
- Choose physical activities that require stamina and avoid those that require strength.
- Avoid overexerting yourself.
- Consult your doctor before beginning a program of physical activity if you've had previous illnesses or are older than 40. Your doctor will advise you on the types of activities and extent of exercise that are appropriate for you.


¹ An epidemiological approach to describing risk associated with blood pressure levels: final report of the Working Group on Risk and High Blood Pressure. *Hypertension*. 1985;7:641-651.

² The fifth report of the Joint National Committee on Detection, Evaluation, and Treatment of High Blood Pressure (JNCV). *Arch Intern Med*. 1993;153:154-183.

³ Stamler J. Neaton JD. Wentworth DN. Blood pressure (systolic and diastolic) and risk of fatal coronary heart disease. *Hypertension*. 1989; 13(suppl):1-2-1-12.


The various components of the blood pressure monitor

Measuring unit

Upper arm cuff:


Type AC1- M for arm circumference 22-30 cm Available as special accessory;
Type AC1 -L for arm circumference 30-38 cm

Operating the blood pressure monitor

Inserting the batteries


After you unpack your monitor, insert the batteries. The battery compartment is located on the back of the monitor (see illustration).

1. Open cover as illustrated.
2. Insert the batteries (4 AA, 1.5V), observing the correct polarity.

Attention!

- If the battery warning symbol appears (as shown in the diagram), the monitor is unusable until the batteries have been replaced.
- Use AA long-life or alkaline 1.5V batteries. The use of rechargeable batteries is not recommended.
- If the blood pressure monitor is left unused for long periods, remove the batteries from the device.


Functional check: Hold down the Start/Stop button to test all the display elements. When functioning correctly, all segments will appear.

Using an AC adaptor (special accessory)

It is possible to operate this blood pressure monitor using an AC adaptor (output 6V DC/600mA).

- Push the power plug into the socket at the back of the monitor.
- Plug the voltage adaptor into a 110V /230V power socket.
- Test that power is available by pressing the Start/Stop button.


Note:

- No power is taken from the batteries while the adaptor is connected to the instrument.
- If the voltage is interrupted during the measurement (e.g., by removing the adaptor), the instrument must be reset by removing the plug from its socket and replacing it.
- Consult a dealer if you have questions relating to the AC adaptor.

Tube connection

Insert the cuff tube into the opening provided on the left side of the monitor, as shown in the diagram.

Carrying out a measurement

Before the measurement:

- Avoid eating, smoking, and all forms of exertion directly before performing the measurement. All these factors influence the measurement results. Relax by sitting in an armchair in a quiet atmosphere for about 10 minutes before the measurement.
- Remove any garment that fits closely to your upper arm.
- Always measure on the same arm (preferably left).
- Attempt to carry out the measurements regularly at the same time of day, since blood pressure changes during the course of the day.

Common sources of error

Note:

Comparable blood pressure measurements always require the same conditions. These should be quiet conditions.



- All efforts by the patient to support the arm can increase the blood pressure. Make sure you are in a comfortable, relaxed position and do not activate any of the muscles in the arm

being measured during the procedure. Use a cushion for support if necessary.


- If the arm artery lies considerably lower or higher than the heart, an erroneous blood pressure will be measured. (Each 15 cm difference in height results in a measurement error of 10 mmHg.)
- Cuffs that are too narrow or too short result in false measurement values. Selecting the correct cuff is very important. The cuff size is dependent upon the circumference of the arm (measured midway between the armpit and elbow). The permissible range is printed on the cuff. If this is not suitable for your use, contact your dealer.
Note: Only use clinically approved original cuffs.
- A loose cuff or a sideways protruding air pocket causes false measurement values.
- With repeated measurements, blood accumulates in the tested arm, which can lead to false results. Correctly executed blood pressure measurements, therefore, should be repeated at least 30 seconds after previous attempts.

Fitting the cuff


1. Pass the end of the cuff (with sewn-in rubber stopper) through the metal stirrup so a loop is formed. The Velcro closure must be facing outward. (Ignore this step if the cuff has already been prepared.)
2. Push the cuff over the left upper arm so the tube points in the direction of the lower arm.


3. Lay the cuff on the arm as illustrated. It is important that the tube coming out of the cuff lies over the brachial artery, located just to the inside of the bicep. Make certain that the lower edge of the cuff lies approximately 2 to 3 cm above the elbow and that the rubber tube leaves the cuff on the inner side of the arm. Tighten the free end of the cuff and close the cuff with the Velcro closure.

4. There must be no free space between the arm and the cuff. Clothing must not restrict the arm. Any piece of clothing that does must be taken off.

5. Secure the cuff with the Velcro closure in such a way that it lies comfortably and is not too tight. Lay the arm on the table (palm upward) so the cuff is at the same height as the heart. Make sure the tube is not kinked.

6. Remain seated quietly for two minutes before you begin the measurement.

Note: If it is not possible to fit the cuff to the left arm, it can be placed on the right one. However, all measurements should be made using the same arm.

Measuring procedure


After the cuff has been positioned appropriately, the measurement can begin.

1. Set the pressure pre-selection.
2. Press the Start/Stop button. The pump begins to inflate the cuff. The rising pressure in the cuff is shown on the display.
3. After the suitable inflation pressure has been reached, the pump stops and the pressure gradually falls. The cuff pressures are displayed.
4. When the instrument detects a pulse, the heart symbol in the display starts to flash and a beep sounds for every heart beat.
5. A longer beep is sounded when the measurement has been completed. The systolic, diastolic blood pressures, and pulse rate now appear in the display.
6. The results of the measurement remain on the display until you turn off the instrument. If no button is pressed for a period of five minutes, the instrument turns itself off to preserve the batteries.

Discontinuing a measurement

If it is necessary to interrupt a blood pressure measurement for any reason (e.g., the patient feels ill), the Start/Stop button can be pressed at any time. The device immediately lowers the cuff pressure.


Further information

Measurements should not occur soon after each other, otherwise the results will be false. Therefore, wait at least 30 seconds in a relaxed position, sitting or lying, before you repeat a measurement.

Memory - displaying the last measurement

The blood pressure monitor automatically stores the last measurement value. Press the Start/Stop button for at least three seconds, the last measurement will be displayed.

Error messages/malfunctions

If an error occurs during a measurement, the measurement will be discontinued and a corresponding error code will be displayed.

Error No.	Possible cause(s)
ER 1	The systolic pressure was determined but afterwards the cuff pressure fell below 20 mmHg. This situation can occur, for example, if the tube becomes unplugged after the systolic blood pressure has been measured. No pulse could be detected.
ER 2	Unnatural pressure impulses impair the result of the measurement Cause: The arm was moved during the measurement (artefact).
ER 3	Inflation of the cuff takes too long. The cuff is not correctly seated, or the hose connection is not sealed.

Other possible malfunctions and their remedies

If problems occur when using the device, the following points should be checked and if necessary, the corresponding measures are to be taken:

Malfunction	Remedy
The display remains blank when the instrument is switched on, although the batteries are in place.	1. Check whether the batteries are installed with the correct polarity, then correct if necessary.
	2. If the display remains blank, remove the batteries, then replace them, or exchange them for new ones.

Malfunction	Remedy
The device frequently fails to measure the blood pressure values, or the values measured are too low or too high.	1. Fit the cuff correctly on the arm.
	2. Before starting the measurement, make sure the cuff is not fitted too tightly and that there is no rolled-up sleeve exerting pressure on the arm above the measuring position. Take off articles of clothing if necessary.
	3. Measure the blood pressure again in complete silence.
Every measurement results in different values, although the device functions normally and normal values are displayed.	Read the following information and the points listed under «Common sources of error». Repeat the measurement.
The blood pressure values measured differ from those measured by the doctor.	Record the daily development of the values and consult your doctor about them.
After the instrument has inflated the cuff the pressure falls very slowly or not at all (no reasonable measurement possible).	Missing plastic ring (distance piece) at tube connection. This may only be removed for servicing reasons. Push this ring over the tube connector so it fits properly and repeat the measurement.

Further information

The level of blood pressure is subject to fluctuations, even with healthy people. **It is important to understand that comparable measurements always require the same conditions.** If, in spite of observing all these factors, the fluctuations are larger than 15 mmHg and/or you hear irregular pulse tones on several occasions, consult your doctor.

This device has been subjected to strict clinical tests. The conformity of this product with the requirements of the CE-guideline 93/42 EWG, council from June 14, 1993 concerning Medical Products and Medical Product-Legislation from August 2, 1994, was established through a conformity check process.

Consult your dealer if there are technical problems with the blood pressure monitor. **Never attempt to repair the instrument**

yourself. Any unauthorized opening of the monitor invalidates all guarantee claims.

Care and maintenance, recalibration

- a. Do not expose the device to extreme temperatures, humidity, dust, or direct sunlight.
- b. The cuff contains a sensitive air-tight bladder. Handle this carefully and avoid all types of straining through twisting or buckling.
- c. Handle the tube carefully. Do not pull on it, allow it to be buckled, or laid over sharp edges.
- d. Clean the device with a soft, dry cloth. Do not use gasoline, thinners, or other petroleum products. Spots on the cuff can be removed carefully with a damp cloth and soap. The cuff must not be machine washed.
- e. Do not drop the instrument or treat it roughly in any way. Avoid strong vibrations.

Never open the device, otherwise the manufacturer's guarantee becomes invalid.

Periodic recalibration

Sensitive measuring devices must be checked for accuracy from time to time. Therefore, we recommend a periodic inspection of the static pressure display every two years. Your dealer can provide more information about this.

Guarantee

The blood pressure monitor is guaranteed for **two years** from date of purchase. This guarantee includes the instrument and cuff. The guarantee does not apply to damage caused by improper handling, accidents, not following the operating

instructions, or alterations made to the instrument by third parties.

Reference to standards

Device standard:	Device corresponds to the requirements of the European standard for non-invasive blood pressure monitor EN1060-1/12:95 EN 1 060-3/09:97 DIN 58130, NIBP - clinical investigation ANSI/AAMI SP10, NIBP - requirements
Electromagnetic compatibility:	Device fulfils the stipulations of the European standard EN 60601-1-2
Clinical testing:	The clinical performance test was carried out in Germany according to the DIN 58130/1997 procedure N6 (sequential).

The stipulations of the EU-Guidelines 93/42/EEG for Medical Products Class IIa have been fulfilled.

Technical specifications

Size:	125 x 155 x 80 mm
Weight:	465 g (with batteries)
Storage temperature:	-5 to +50°C
Humidity:	15 to 85% relative humidity maximum
Operation temperature:	10 to 40°C
Display:	LCD-Display (Liquid Crystal Display)
Measuring method:	Oscillometric
Pressure sensor:	Capacitive
Measuring range:	SYS/DIA: 30 to 280 mmHg Pulse: 40 to 200 per minute
Cuff pressure display range:	0-299 mmHg
Memory:	Storing the last measurement automatically
Measuring resolution:	1 mmHg
Accuracy:	Pressure within ± 3 mmHg Pulse $\pm 5\%$ of the reading
Power source:	a) 4 dry cells (Batteries); size AA, 1.5V b) AC adapter 6V DC 600mA (voltage 4.5V DC to 6V DC)
Accessories:	Medium Adult Cuff: AC1-M for arm circumference 22 to 30 cm Available by separate order: Large Adult Cuff: AC1-L for arm circumference 30 to 38 cm AC Voltage Adapter, 6V DC, 600mA

Technical and color alterations reserved!

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Printed in U.S.A.