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Gamma
TECHNICS FOR HEALTH

Automatic Blood Pressure Monitor Gamma™ Model M1-3

INSTRUCTION BOOK

Thank you for choosing our Automatic Blood Pressure Monitor Gamma™ designed on the basis of the most advanced technologies. The device you are holding now is your reliable assistant. You can measure your blood pressure with this monitor as accurately and precisely as an experienced doctor or a skilled nurse would do it. We are sure that having appraised worthily the high quality and reliability of this device you will become a regular user of the products of Gamma Trade Mark.

Before starting to use this device please study the instruction book carefully. The instruction book offers all information you need to correctly measure your blood pressure and pulse.

Regular blood pressure control allows for reliable and easy health monitoring and for effective combating possible complications caused by blood pressure abnormalities.

ATTENTION

This blood pressure monitor is designed to carry out self-control over blood pressure but NOT to make self-diagnosis of hypertension/hypotension. Please DO NOT diagnose by yourselves basing on the measurement results obtained with the blood pressure monitor. Please DO NOT execute self-treatment of high/low blood pressure and DO NOT change the methods prescribed without consulting your doctor.

CONTENT

1. Features of Automatic Blood Pressure Monitor Gamma™ Model M1-3

- 1.1. Components of your blood pressure monitor kit
- 1.2. Description of the blood pressure monitor

2. Getting ready for measurement

- 2.1. Inserting the batteries
 - 2.1.1. Using an A/C power adapter
- 2.2. Cuff connection
- 2.3. Setting the time and date

3. Measurement Procedure

- 3.1. Fitting the cuff
- 3.2. Measuring procedure
- 3.3. Reading measurement results
- 3.4. Discontinuing a measurement
- 3.5. Irregular Heartbeat Detector

4. Memory function

5. Error Messages/Troubleshooting

6. Care and maintenance

7. Technical specifications

8. Warranty

9. Important information on blood pressure and its measurement

- 9.1. What is blood pressure?
- 9.2. Which values are normal?
- 9.3. What can be done to lower/maintain the level of blood pressure?


1. Features of Automatic Blood Pressure Monitor Model M1-3

1.1. Components of your blood pressure monitor kit


The blood pressure monitor kit includes Automatic Blood Pressure Monitor Model M1-3, the cuff, 4 batteries (AA), bag for storing of devices, Instruction Book, Warranty Card, kit box.

1.2. Description of the blood pressure monitor

Measuring unit Model M1-3

SYMBOLS ON DISPLAY


2. Getting ready for measurement

2.1. Insert or replace batteries

1. Remove the battery cover.
2. Insert new batteries into the battery compartment as shown, taking care that the polarities(+) and (-)are correct.
3. Close the battery cover. Use only Alkaline batteries (LR6, AA).

CAUTION

- Insert the batteries as shown in the battery compartment. If not, the device will not work.
- When  (LOW BATTERY mark) blinks in the display, replace all batteries with new ones. Do not mix old and new batteries. It may cause the device to malfunction.  (LOW BATTERY mark) does not appear when the batteries run out.

4X1.5V Batteries
(LR6 or AA)

- After the battery warning  appears, the device is blocked until the batteries have been replaced.
- Battery life varies with the ambient temperature and may be shorten at low temperatures.
- If the blood pressure monitor is not used for more then a week, remove the batteries from the device.
- The battery may leak and cause a malfunction of the device.
- Do not use rechargeable batteries.

Disposal and empty battery to the authorized collection party subject to the regulation of each individual territory

Use the specified batteries only. The batteries provided with the device are for testing monitor performance and may have a shorter life.


2.1.1. Using an A/C power adapter

It is possible to operate Automatic Blood Pressure Monitor Model M1-3 with an AC/DC adaptor.

1. Push the plug into the socket at the bottom part of the device. No power is taken from the batteries while the AC/DC adaptor is connected to the device.

2. Plug the AC adaptor into a 220/110 V power socket. The device is ready for measurement.

3. Test that power is available by pressing the START button.


Note:

Please use only the AC adaptor supplied in the set. Any other adaptor being used, the service centre is not considered responsible for the warranty servicing/repair of the monitor. No power is taken from the batteries while the AC/DC adaptor is connected to the device. Please consult if you have questions relating to the AC/DC adaptor.

2.2. Cuff connection


Insert the cuff tube into the opening provided on the side of the device as shown in the picture.

2.3. Setting the time and date

1. Under Power Off. Press "SET" button, automatically power off in 3 seconds.

2. Press and hold "SET" button until the year number displays and flashes on LCD to enter setting mode. And press "121" button to adjust the year. Confirm your selection by pressing «SET» button.

3. Following the same steps to chose year, month, day, hour, minutes, unit conversion (mmHg/kPa), press “” to change the value.

year

month

day

hour

minute

unit

ATTENTION

At installation of unit conversion be attentive. Choose units of measure mmHg.

3. Measurement Procedure


ATTENTION

- Find time (5-10 min.) to relax by sitting in a quiet atmosphere for some time before measurement.
- Efforts by the patient to support the arm can increase the blood pressure. Make sure you are in a comfortable, relaxed position and do not activate any muscles in the arm during measurement.
- Always measure on the same arm (normally left).
- Remove any garment that fits closely to your upper arm. Do not roll the sleeve since it can squeeze your hand and this can lead to false results.
- Use only clinically approved original cuff.
- If you want to follow the results of your blood pressure measurements, always perform measurements at the same time of day, since blood pressure changes during the course of day.
- Measurements should be done after a 5 minute rest to ensure accuracy.

3.1. Fitting the cuff


1. Pass the end of the cuff through the flat metal ring so that a loop is formed, and the velcro faces outwards.

2. Place the cuff over the left upper arm


so that tube is closer to your lower arm.


3. Lay the cuff on the arm so that the lower edge of the cuff lies approximately 2 to 3 cm above the elbow.


4. Tighten the cuff by pulling the end and close the cuff by affixing the velcro. There should be little free space between the arm and the cuff. Cuff that does not fit properly results in false measurement values. Clothing must not restrict the arm. Any piece of clothing which does must be removed.

5. Lay your arm on a table so the cuff is at the same height as your heart. If the upper arm artery lies considerably lower or higher than the heart, a false higher or lower blood pressure will be measured! A variation of 15 cm between cuff and heart level can result in a reading error of + or - 10mm Hg.

Make sure the tube is not kinked.

6. Remain seated quietly for two minutes before you begin the measurement.


3.2. Measuring procedure

1. Press  or  button, the display shows "0", Inflation icon  flashes on display when it begins to inflate.

During the inflation, please do not move


2. Automatic deflation after measurement, and display the blood pressure, heart rate, and blood pressure indicator

- Mistake, display "E", please refer the instruction about troubleshooting.

- When cuff remove during the measurement, please fix the cuff and test again.
- Stop the measurement in emergency, please press  or  power. Automatic deflation at the same time.

3. Power off. Press  or  ,to cut off the power.



3.3. Reading measurement results

When the measurement has been concluded, the measured blood pressure and pulse values are now displayed.

The appearance of this symbol «♥» signifies that an irregular heartbeat was detected. This indicator is only a caution. It is important that you be relaxed, remain still and do not talk during measurements.

On the graphic indicator the element showing level of arterial pressure on a scale of World Health Organization is displayed


Note: We recommend contacting your physician if you see this indicator frequently.

To prolong the batteries' life the device switches off automatically if no button is pressed for 1 minute. Otherwise you can switch it off by pressing the  or  button.

3.4. Discontinuing a measurement

If it is necessary to interrupt a blood pressure measurement for any reason (e.g. the patient feels unwell), the  or  button can be pressed at any time. The device then immediately lowers the cuff pressure automatically.

3.5. Irregular Heartbeat Detector (IHD-technology)

This function indicates allorhythmic heartbeating. If the symbol of IHD «♥» appears on the display that means that certain abnormality in heart beat frequency was detected during the measurement. In this case, the result may deviate from your normal blood pressure – repeat the measurement. In most cases, this is no cause for concern. However, if the symbol «♥» appears on a regular basis (e.g. 3-5 times a week with measurements taken daily) we advise you to tell your doctor. Please show your doctor the following explanation:

Information for the doctor on frequent appearance of the Irregular Heartbeat Detector

This instrument is an oscillometric blood pressure monitor that also analyses pulse frequency during measurement. The instrument is clinically tested.

The symbol «♥» is displayed after the measurement, if pulse irregularities occur during measurement. If the symbol «♥» appears more frequently (e.g. several times per week on measurements performed daily) we recommend the patient to seek medical advice.

The instrument does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

4. Memory function

1. Under power off. Hold  or  for 3 seconds, get into relative memory mode, it reads the average of latest 3 measurements, press "SET" to find more memory.

- After measurement, press "SET" to read the average of latest 3 measurements.

2. Press  or  to cut off the power.


- If forget to turn it off. Automatically power off in one minute.
- More than 60 measurements, it will deletes the earliest one automatically
- The stored memory will not lose after replacing battery.

3. Delete memory (all the memory will be deleted). Showing memory, hold "SET" for 5 seconds, display "No" means all the memories have been deleted.

5. Error Messages/Troubleshooting

If an error occurs during a measurement, the measurement is discontinued and a corresponding error code is displayed (picture: error no. 2)

- Please consult your doctor, if this or any other problem occurs repeatedly

Other possible errors and their solutions

If problems occur when using the device, the following points should be checked:

ERROR DISPLAY	POSSIBLE CAUSE	HOW TO CORRECT
Nothing is displayed When you push the POWER button or battery icon  flash	No battery installation Battery worn out The polarities of batteries placed wrongly	Insert batteries Replace new batteries Insert new batteries in accordance to this instruction
E1 :can't normally Increase pressure	Check your cuff if any air leakage	Replace cuff with new one
E2 E4 : have shaking while measurement	Hand or body shaking while measurement	keeping static and correct gesture to measure again
E3 inflate pressure too high		Re-measurement or send back to Service Center for re-calibrate pressure
Battery icon  on	Battery low power	Replace battery and measure again

The systolic pressure Value or diastolic Pressure value too high	1. The arm cuff was held lower than your heart 2. The arm cuff was not attached properly 3. You moved your body or spoke during measurement	keeping correct position and gesture to measure again
The systolic pressure Value or diastolic Pressure value too low	1. The arm cuff was held higher than your heart 2. you moved your body or Spoke during measurement	keeping correct position and gesture to measure again

Note

Blood pressure is subject to fluctuations even in healthy people. Please remember that **comparable blood pressure measurements always require the same time and the same conditions! These are normally quiet conditions.**




ATTENTION

Should any technical malfunctions arise in the blood pressure monitor, please contact the point of sale where you bought the device or your nearest Service Center listed in the warranty card.

Never attempt to repair the instrument yourself!

Any unauthorized opening of the instrument invalidates all warranty claims.

6. Care and maintenance

<p>Do not expose the device to either extreme temperatures, humidity, dust or direct sunlight. Store the device at least 2 m from sources of electromagnetic fields.</p>	
<p>Handle the cuff carefully and avoid all types of stress through twisting or buckling in order not to damage the sensitive air-tight bubble.</p>	
<p>Clean the device with a soft, dry cloth. Do not use gas, thinners or similar solvents. Spots on the cuff can be removed carefully with a damp cloth and soapsuds. Do not submerge the cuff in water!</p>	
<p>Handle the tube carefully. Keep the hose away from sharp edges.</p>	
<p>Do not drop and hit the monitor or treat it roughly in any way.</p>	
<p>Never open the monitor! Otherwise the manufacturer's calibration will be spoiled.</p>	

7. Technical specifications*

Measuring Method:	Oscillometric Measurement
Indication:	Digital LCD display
Measuring Range:	Pressure: (30-280) mmHg Pulse: (40-199) Beat/min
Accuracy:	Static Pressure: ± 3 mmHg Pulse: $\pm 5\%$
Memory:	60/60 Memories
Power supply:	4x1.5V Batteries(LR6 or AA)
Operating condition:	+5°C~+40°C. 30%RH~80%RH
Storage condition:	-20r~+55t:. 10%RH~93%RH
Dimensions:	Approx 156(W)X105(H)X65(D) mm
Weight:	Approx,450g, excluding batteries
Classification:	Type BF
Upper arm circumference:	(22-32) cm

*Specifications may be changed without notice in the event of improvement being made!

8. Warranty

Your Automatic Blood Pressure Monitor Model M1-3 is warranted for 3 years from date of purchase. The warranty doesn't cover the cuff.

Service life of the device is 10 years after date of purchase with following all operating instructions.

The warranty does not apply to damage caused by improper handling, accidents, not following the operating instructions or self-maintained alterations made to the device.


The warrantee is only valid upon presentation of the warrantee card which was correctly filled in and sealed.

9. Important information on blood pressure and its measurement

9.1. What is blood pressure?


Blood pressure is the pressure that your blood extends to the vascular walls. Blood pressure is necessary to provide for constant blood flow inside the body. Thanks to it the cells get oxygen that provides for their normal functioning. The heart performs the function of a “pump”, sending blood to the blood vessels. Each heart beat creates a certain level of the blood pressure.

There are 2 kinds of blood pressure: a systolic (upper) one, which corresponds to the heartbeat pushing blood into the arteries; and a diastolic (lower) one, which means the blood pressure between two heartbeats.

Blood pressure is subject to fluctuations during the day even in healthy people. The fluctuations are influenced by a number of factors – time of day, person’s condition, physical or mental activity, environment, etc.

Day-Night Fluctuations of the Blood Pressure

An increase of blood pressure increases the burden onto the heart, affects blood vessels making their walls thick and less elastic. One of the features of the hypertension is its ability to remain unnoticed for the patient at its early stages. That's why the self-control of the blood pressure is so important. With the illness progressing, headaches and regular dizziness appear, the sight declines, the functioning of vitals (encephalon, heart, kidneys, blood vessels) breaks down. Without special treatment the complications of hypertension might be kidney damages, breast-pang, paralytic stroke, aphasia, dementia, heart attack and stroke.

9.2. Which values are normal?

The world standard as for the norms of the blood pressure is the Classification* of the World Health Organization (WHO):

Blood Pressure Category	Systolic Blood Pressure (mmHg)	Systolic Blood Pressure (mmHg)	Recommendation
blood pressure too low	< 100	< 60	Consult your doctor
blood pressure optimum	100 ... 120	60 ... 80	Self-check

blood pressure normal	120 ... 130	80 ... 85	Self-check
blood pressure slightly high	130 ... 140	85 ... 90	Consult your doctor
blood pressure too high	140 ... 160	90 ... 100	Seek medical advice
blood pressure far too high	160 ... 180	100 ... 110	Seek medical advice
blood pressure dangerously high	> 180	> 110	Urgently seek medical advice!

* printed with curtailments


- The diagnosis of hypertension requires from the patient to combine medical treatment prescribed by the doctor and mode of life correction.
- People with normal pressure and high normal pressure are recommended to carry out self-control of their tension in order to timely take measures to decrease the blood pressure level down to the optimal one without using any medications.
- For people more than 50 years old high level of systolic blood pressure (higher than 140 mmHg) is more crucial than diastolic pressure.
- Even with blood pressure being normal, people run the bigger risk of hypertension development with advancing age.

ATTENTION

If you have normal results of blood pressure measured under calm conditions but your results are excessively high when measured under the conditions of physical or mental exhaustion, this might be a sign of so called brittle (that is unstable) hypertension. If you suspect that, please consult your doctor.

When measured correctly, if diastolic blood pressure is more than 120 mmHg, it is necessary to call the doctor immediately.

9.3. What can be done to lower/maintain the level of blood pressure?

Recommendation		How much the systolic blood pressure can be lowered
- Try to maintain the weight which is considered standard for your age. Decrease overweight if any!		5 – 20 mmHg with every 10 kilos of weight reduction
- Avoid salt overuse (it's recommended to use no more than 6 g per day)!		2 – 8 mmHg
- Avoid eating fat food! Introduce to your diet some food that is rich in dietary fiber (fruits, vegetables)		8 – 14 mmHg
- Refuse from smoking outright! - Use only limited amount of alcohol (beer – 350 g, or wine – 150 g, or vodka – 50 g)! - Restrain the usage of caffeine (coffee)!		2 – 4 mmHg

<p>- After medical examination please regularly go in for sports, with the following restrictions:</p> <ul style="list-style-type: none"> • go in for endurance exercises not strengthening exercises. • do not train till complete exhaustion. • if you have any chronic diseases and/or you are older than 40 please consult your doctor before you start going in for sports. <p>- Spend at least half an hour out in the open air every day.</p>		<p>4 – 9 mmHg</p>
---	---	-------------------

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