
MISURATORE DI PRESSIONE AUTOMATICO DA BRACCIO / AUTOMATIC ARM BLOOD PRESSURE MONITOR / MESUREUR DE PRESSION AUTOMATIQUE PAR LE BRAS / AUTOMATISCHES ANDGELENK-BLUTDRUCKMESSGERÄT / MEDIDOR DE TENSIÓN AUTOMÁTICO DE BRAZO / MEDIDOR DE TENSÃO AUTOMÁTICO DE BRAÇO / АВТОМАТИЧЕН АПАРАТ ЗА ИЗМЕРВАНЕ НА КРЪВНО НАЛЯГАНЕ, ЗА РЪКА / TENSIOMETRU AUTOMAT PENTRU BRAȚ / AUTOMATA KAR-VÉRYNOMÁSMÉRŐ

Descrizione dell'apparecchio e simboli sul display / Description of the device and symbols on the display / Description de l'appareil et symboles affichés sur l'écran / Gerätebeschreibung und Symbole in der Anzeige / Descripción del aparato y símbolos de la pantalla / Descrição do aparelho e símbolos no visor / ОПИСАНИЕ НА АПАРАТА И НА СИМВОЛИТЕ ВЪРХУ ДИСПЛЕЯ / Descrierea aparatului și simbolurile de pe ecran / A készülék leírása és a kijelzőn megjelenő szimbólumok

SYS

Pressione sistolica / Systolic blood pressure / Pression systolique / Systolischer Blutdruck / Presión sistólica / Tensão sistólica / Систолно налягане / Tensiune sistolică / Szisztolés vérnyomás

DIA

Pressione diastolica / Diastolic blood pressure / Pression diastolique / Diastolischer Blutdruck / Presión diastólica / Tensão diastólica / Диастолно налягане / Tensiune diastolică / Diasztolés vérnyomás

Pul/min

Pulsazioni al minuto / Pulse- Beats per minute / Pulsations par minute / Puls- Herzschläge pro Minute / Pulsaciones por minuto / Pulsações por minuto / Сърдечни удари в минута / Bătăi pe minut / Percenkénti pulzusszám

Sgonfiaggio: Il bracciale sta terminando a sgonfiarsi / Deflation: The arm cuff is deflating / Dégonflement: Le bracelet termine de se dégonfler / Entleeren: Die Manschette hat sich entleert / Desinflado: El brazalete está terminando de desinflarse / Desinsuflamento: A braçadeira está a acabar de desinsuflar-se / Спадане: Маншетът приключва с изпускането на въздуха / Dezumflare: Manșeta finalizează dezumflarea / A karöltő majdnem leengedett

RM 8:59

Date (ora-minuti) / Date (hour-minutes) / Date (heure-minutes) / Datum (Stunde-Minuten) / Fecha (hora-minutos) / Data (hora-minutos) / Дата (час-минути) / Dată (oră-minute) / Dátum (óra-perc)

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AUTOMATIC ARM BLOOD PRESSURE MONITOR

The Blood Pressure Monitor JOYCARE detects the Systolic, Diastolic, and the Pulse rate using the oscillometric method. All values are displayed on the large LCD screen. Please read carefully this entire instruction manual before use and keep it. For further questions about blood pressure and its measurement, please contact your doctor.

IMPORTANT INFORMATION BEFORE MEASURING

Blood pressure changes during the day and is also influenced by many factors such as smoking, drinking alcohol, taking drugs and physical activity. The blood pressure measurement should be performed by a doctor or a healthcare assistant who knows your medical history thoroughly. The regular use of the JOYCARE'S blood pressure meter will allow you to keep your doctor abreast of the continuing trend of your blood pressure, and therefore provide a useful tool for your health diagnosis. Ensure you wear the cuff correctly, paying particular attention to its position relative to your heart, since a position above or below it could affect the measurement. Do not move the device, body or arm and do not talk during the measurement, otherwise the readings will not be reliable. Perform the measurement while in a comfortable position, sitting, resting quietly and at ease. Never wrap the cuff around the sleeve of your jacket or shirt, otherwise the readings are not reliable. Please allow at least 3-5min between measurements to restore normal circulation. Measure the pressure always in the same arm, preferably the left one, as the pressure between arms may differ. **Certain factors may affect the measurement, cause variations, and unreliable readings, including bathing, talking during the measurement, alcohol, exercise or agitated movement, large meals, sudden temperature changes, nervous tension, stress, etc...**

CAUTION

This device is intended for adult use only.

This device is intended for no-invasive measuring and monitoring of arterial blood pressure.

It is not intended for use on extremities other than the arm or for functions other than obtaining a blood pressure measurement.

Do not confuse self-monitoring with self-diagnosis. This unit allows you to monitor your blood pressure. Do not begin or end medical treatment based solely physician for treatment advice.

If you are taking medication, consult your physician to determine the most appropriate time to measure your blood pressure. Never change a prescribed medication without consulting your Physician.

When the device was used to measure patients who have common arrhythmias such as atrial or ventricular premature beats or atrial fibrillation, the best result may occur deviation. Please consult your physician about the result.

If the cuff pressure exceeds 40 kPa (300 mmHg), the unit will automatically deflate. Should the cuff not deflate when pressures exceeds 40 kPa (300 mmHg), detach the cuff from the arm and press the START/STOP button to stop inflation.

The equipment is not AP/AGP equipment and not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide. The operator shall not touch output of batteries /adapter and the patient simultaneously.

To avoid measurement errors, please avoid the condition of strong electromagnetic field radiated interference signal or electrical fast transient/burst signal. The user must check that the equipment functions safely and see that it is in proper working condition before being used.

This device is contraindicated for any female who may be suspected of, or is pregnant. Besides provided inaccurate readings, the affects of this device on the fetus are unknown.

Manufacturer will make available on request circuit diagrams, component parts list etc.

This unit is not suitable for continuous monitoring during medical emergencies or operations.

Otherwise, the patient's arm and fingers will become anaesthetic, swollen and even purple due to a lack of blood.

Please use the device under the environment which was provided in the user manual. Otherwise, the performance and lifetime of the device will be impacted and reduced.

During using, the patient will contact with the cuff. The materials of the cuff have been tested and found to comply with requirements of ISO 10993-5:2009 and ISO 10993-10:2010. It will not cause any potential sensitization or irritation reaction.

Please use ACCESSORIES and detachable partes specified/ authorised by MANUFACTURE.

Otherwise, it may cause damage to the unit or danger to the user/patients.

The device doesn't need to be calibrated in two years of reliable service.

Please dispose of ACCESSORIES, detachable parts, and the ME EQUIPMENT according to the local guidelines.

If you have any problems with this device, such as setting up, maintaining or using, please contact with Authorized Service Center.

Don't open or repair the device by yourself.

Please report to an authorized Service Center if any unexpected operation or events occur.

Please use the soft cloth to clean the whole unit. Don't use any abrasive or volatile cleaners.

The device must be installed and put in service in compliance with the EMC information contained in this manual.

The portable and mobile radio communication devices may affect the functioning of the device

Do not use near (e.g. 1 m) a DEVICE for short wave or microwave therapy;

Do not connect the patient simultaneously to a high frequency electrosurgery device;

For all technical information please contact the Authorized Service Center JOYCARE

CAUTION: Do not forget: while self-measurement represents control, it does not constitute a proper diagnosis or treatment. Any unusual values should always be discussed with your doctor. **Under no circumstances should you change the dosage of any medication prescribed by your doctor. This blood pressure monitor is not suitable for people with diagnosed cardiac arrhythmia.** Should the user should suffer from heart problems, kidney problems, diabetes and circulatory problems, it shall be necessary to seek the doctor's advice before using the device. This device is not a toy. Keep out of reach of children. Do not leave it unattended and make sure that children do not play with it. This device shall not be used should you have any arm wounds. This device should only be used for measuring blood pressure in adults. The device contains sensitive electronic components (Microcomputer). Therefore, avoid strong electrical or electromagnetic fields in the immediate vicinity of the device (e.g. cell phones, microwave ovens). These can cause a temporary loss of measurement accuracy. To avoid possible damage to the device, always use and store it within the temperature and humidity ranges contained in the technical specifications.

POWER SPECIFICATIONS:

4 AA batteries included

DC 6 V power adaptor (not included)

Network adaptor connector

INSERTING AND REMOVING THE BATTERIES

Remove the cover on the battery compartment at the rear of the device and insert 4 new AA batteries, respecting the polarity indicated, then replace the cover. Replace the batteries when the battery symbol  appears on the display, and when the display starts to wane or doesn't power on. If the product is not used for an extended period of time, it is advisable to remove the batteries, since a leak from the batteries can cause damage to the device.

SETTING THE DATE AND TIME:

Before using your blood pressure monitor, you must set the date and time so that readings can be memorized with their date and time. Once the batteries have been inserted, press and hold the "SET" key for approximately 3 seconds; the year will begin to flash on the display; now press the "MEM" key to set the year, followed by "SET" to confirm; the month will start to flash, press the "MEM" key to set the month, then the "SET" key to confirm. Repeat the same sequence to set the day, time and minutes. Once all the settings have been set, the display shows , and the device will power itself off automatically after a few seconds.

RECOMMENDATIONS FOR TAKING READINGS

5 minutes before taking a blood pressure reading, sit down, relax and remain calm. Wait at least 3 minutes between one measurement and the next. Keep your body and arm cuff still and do not speak during the reading. Avoid touching the device. Always take readings on the same arm. Uncover your arm so that it is exposed in the area in which the reading is being taken, and do not tighten the arm cuff excessively, allowing a space equivalent to the width of a finger. Take your blood pressure reading while seated in an upright position, resting your lower arm onto a table or other surface. If necessary, place something under your arm to adjust the height (e.g. a book). Relax your arm with the palm of your hand facing upwards. The arm cuff must be positioned at the same level as your heart. Because blood pressure varies throughout the day, it is advisable to measure your blood pressure always at the same time of day, so as to be able to make a reliable comparison.

TAKING BLOOD PRESSURE READINGS

Connect the arm cuff's air tube to the main unit. Fit your arm in through the arm cuff and pull the cuff so that it is positioned around 1 or 2 cm above the fold of your elbow; the air tube must be positioned in line with your little finger. Press the START/STOP button. All of the digits on the display will appear simultaneously, then "0" will appear on the display, together with the arrow

symbol . Immediately after, the arm cuff will begin to inflate and the device will begin the blood pressure measurement. Avoid speaking or moving while taking readings. After taking a reading, the display will show the systolic and diastolic pressure, as well as the heart rate, storing them in memory automatically. The date and time also appear on the display, along with the symbol  next to a different colour depending on the World Health Organization classification (see paragraph on "Blood pressure classification standards"). The results of the reading will appear on the measurement bar, according to the WHO classification:

GREEN: blood pressure from normal to slightly high

YELLOW: moderately high blood pressure

ORANGE: very high blood pressure

RED: severely high blood pressure

In addition, in the case of an irregular heart rate, the symbol shown here will appear on the display. In these conditions, the device will take a reading anyhow, but the data obtained may not be exact; it is advisable in this case to contact your physician . This device is not suited to patients with a proven cardiac arrhythmia condition.

Press the START/STOP button to switch off the device, otherwise it will power off automatically after a few seconds.

SAVING, RECALLING AND CANCELLING DATA


This device is capable of memorizing data for up to 60 readings. After each reading, the device stores in memory all data from the blood pressure measurement (Systolic pressure, Diastolic pressure, heart rate, date, time and WHO classification). When the device is off, press the "MEM" key to recall memorized readings. The data from the last measurement appears under "MEM 1". To view all other memorized readings simply press the SET" or "MEM" keys. The device can memorize 60 readings in total, after which it will progressively delete the oldest readings. (starting from no."60"). To delete/cancel the data stored in the memory, simply press and hold the "MEM" key for a few seconds, when the device is off. The heading "DEL ALL" will appear on the display; press the "SET" key to delete all memorized data, and "dEL dOnE " will appear on the display. If no data is stored in the memory, "----" will appear on the display.

BLOOD PRESSURE CLASSIFICATION STANDARDS

The World Health Organization (WHO) has set the standards for blood pressure classification, which you may find in the table below.

Note: This device is suitable for measuring blood pressure exclusively in adults. Risk factors such as diabetes, obesity, smoking etc. are not considered in these tables. This table is only indicative, these values indicate exclusively the different levels of blood pressure, and under **no circumstances you should change the dosage of any medication prescribed by your doctor. For a final evaluation of your blood pressure, it is essential to seek always your doctor's advice.**

<i>Blood pressure classification</i>	<i>Pressione sistolica mmHG</i>	<i>Pressione diastolica mmHG</i>	<i>Colour</i>
Excellent	< 120	<80	GREEN
Normal	120-129	80-84	GREEN
Moderately high	130-139	85-89	GREEN
Hypertension degree 1	140-159	90-99	YELLOW
Hypertension degree 2	160-179	100-109	ORANGE
Hypertension degree 3	≥180	≥110	RED

MAINTENANCE AND CARE OF THE DEVICE

Clean the device using a cloth dampened with water or a delicate detergent, and dry using a supple cloth. Never use aggressive detergents or sprays to clean the device. Never immerse the arm cuff or the device in water. In the eventuality the device is not used for a long period of time, it is advisable to remove the batteries, since leaked battery fluid can damage the appliance. Do not expose the device to direct sunlight. Never store the device in a humid or dusty environment.

PREVENTING MALFUNCTIONS

Do not twist the arm cuff in any way. Do not inflate the arm cuff if it hasn't been positioned around the arm beforehand. Never attempt to disassemble or modify any component/s on the device, including the arm cuff, in an attempt to repair it. In the eventuality of a malfunction, contact an authorized service centre. Avoid dropping the device, and do not expose it to any sort of shock. Always store the device away from any heat source, the sun's rays, excessive humidity, sharp objects, etc.

Problem	Situation	Verify	Intervention
The product does not work	The display is dimmed or not lit.	The batteries are exhausted	Install new batteries
		The batteries were not installed correctly	Insert the batteries as required
The batteries are exhausted	The symbol  +Lo appears on the display	The batteries are exhausted	Install new batteries
Error Message	The display shows E1	The cuff was not properly connected	Fasten the cuff correctly and repeat the measurement
	The display shows E2	The cuff is too tight	Fasten the cuff correctly and repeat the measurement
	The display shows E3	The cuff pressure is too high	Relax for a few minutes and repeat the measurement
	The display shows E10 or E11	The monitor has detected a movement during the measurement	Shifting your position during the measurement may adversely affect the measurement. Relax for a few minutes and repeat the measurement
	The display shows Eexx	Incorrect measurement	relax for a few minutes and repeat the measurement
		Calibration error	Repeat the measurement and contact an authorized service centre should the problem persist.

Should the problems persist, please contact an authorized service centre.
Should the readings be excessively high or low, please seek your doctor's advice.

TECHNICAL SPECIFICATIONS

Model	JC-610 (TMB 995)
Description	Arm blood pressure monitor
Method of measurement	Oscillometric
Mode of operation:	Continuous operation
Protection against ingress of water:	IP21
Degree of protection:	Type BF applied part
Measurement range	Rated cuff pressure: 0mmHg-300mmHg (0kPa-40kPa)
Measurement pressure:	40mmHg-230mmHg (5.3kPa-30.7kPa)
Pulse value:	(40-199) beat/minute
Accuracy	Pressure: 5°C-40°C within ± 3 mmHg (0.4kPa)- Pulse value: $\pm 5\%$
Power	appliance powered by internal electric source, 4 AA batteries (included) DC 6V 1A power pack (not included) from +5°C to +40°C , humidity $\leq 85\%$ from -20°C to +60°C , humidity 10% RH~ 93%RH
Operating conditions	LCD V.A. 80X60mm
Storage and transportation conditions	60 measurements with date and time
Display	approx 102x143x73 mm
Memory	about 500gr, batteries included
Product dimensions	approx. 22~42cm
Product weight	automatic after around 1 minute unused
Circumference of arm cuff	arm cuff with tube, instruction manual, 4 "AA" batteries(test batteries are included).
Power-off	The supplied batteries are for trial use only and they may have a shorter life.
Accessories	

Expected service life of device: about 9 years if calculated by 3 times measurement per day.

Symbols:  Type BF applied part  See the documentation attached  Date of manufacture

Warning: The Blood Pressure Monitor JOYCARE was designed to measure blood pressure exclusively in adults.

JC-610 (TMB 995) has been manufactured in compliance with the European directives 93/42/ECC 0197

MADE
IN CHINA

Distributed under the brand JOYCARE, by JOYCARE S.R.L., sole partner - Via Massimo D'Antona, 28 - 60033 Chiaravalle (AN) - Italy

GUANGDONG TRANSTEK MEDICAL ELECTRONICS CO.,LTD.

Zone A, No.105 ,Dongli Road, Torch Development District,Zhongshan,528437,Guangdong,China

MDSS - Medical Device Safety Service GmbH, Schiffgraben 41, 30175, Hannover - Germany.

DISPOSAL The device (including its removable parts and accessories) must not be disposed of together with municipal waste at the end of its life, but in compliance with European Directive. Since it must be handled separately from household waste, either carry it to a separately-collected waste disposal centre for electrical and electronic appliances or give it back to the retailer on purchasing a new device with the same purpose. Any infringement will be severely prosecuted. **The batteries used in this device must be disposed of in the special bins at the end of their life.**

All the instructions and drawings are based on the latest information available at the time of going to print and may be subject to change.

COMPLIED EUROPEAN STANDARDS LIST	
Risk management	ISO/EN 14971:2012 Medical devices — Application of risk management to medical devices
Labeling	ISO/EN 15223-1:2012 Medical devices. Symbols to be used with medical device labels, labelling and information to be supplied. General requirements
User manual	EN 1041: 2008 Medical equipment manufacturers to provide information
General Requirements for Safety	EN 60601-1: 2006 Medical electrical equipment – Part 1: General requirements for basic safety and essential Performance IEC/EN 60601-1-11: 2010 Medical electrical equipment – Part 1-11: General requirements for basic safety and essential performance - Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment
Electromagnetic compatibility	IEC/EN 60601-1-2:2007 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests
Performance requirements	EN 1060-1:1995+A2:2009 Non-invasive blood pressure Part 1: General requirements EN 1060-3:1997+A2:2009 Non-invasive blood pressure . Part 3: Supplementary requirements for electromechanical blood pressure measuring system
Clinical investigation	EN 1060-4: 2004 Automatic Blood Pressure Monitor overall system Interventional accuracy of the testing process
Usability	IEC/EN 60601-1-6: 2010 Medical electrical equipment -- Part 1-6: General requirements for basic safety and essential performance - Collateral Standard: Usability IEC/EN 62366: 2007 Medical devices - Application of usability engineering to medical devices
Software life-cycle processes	IEC/EN 62304:2006+AC: 2008 Medical device software - Software life cycle processes

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GUANGDONG TRANSTEK MEDICAL ELECTRONICS CO.,LTD.

Zone A, No.105 , Dongli Road, Torch Development District,Zhongshan,528437,Guangdong,China

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**EMC Guidance - Table 1 Guidance and MANUFACTURER's declaration
ELECTROMAGNETIC EMISSIONS- for all ME EQUIPMENT and ME SYSTEMS**

Guidance and manufacturer's declaration - electromagnetic emissions

*The device is intended for use in the electromagnetic environment specified below.
The customer or the user of the device should assure that it is used in such an environment.*

<i>Emissions test</i>	<i>Compliance</i>	<i>Electromagnetic environment - guidance</i>
<i>RF emissions CISPR 11</i>	<i>Group 1</i>	<i>The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.</i>
<i>RF emissions CISPR 11</i>	<i>Class B</i>	
<i>Harmonic emissions IEC 61000-3-2</i>	<i>Not applicable</i>	
<i>Voltage fluctuations/ flicker emissions IEC 61000-3-3</i>	<i>Not applicable</i>	

**Table 2 Guidance and MANUFACTURER's declaration – electromagnetic IMMUNITY
for all ME EQUIPMENT and ME SYSTEMS**

Guidance and manufacturer's declaration – electromagnetic immunity

*The device is intended for use in the electromagnetic environment specified below.
The customer or the user of the device should assure that it is used in such an environment*

<i>IMMUNITY test</i>	<i>IEC 60601 test level</i>	<i>Compliance level</i>	<i>Electromagnetic environment - guidance</i>
<i>Electrostatic discharge (ESD) IEC 61000-4-2</i>	<i>±6 kV contact ±8 kV air</i>	<i>±6 kV contact ±8 kV air</i>	<i>Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.</i>
<i>Electrical fast transient/burst IEC 61000-4-4</i>	<i>±2 kV for power supply lines ±1 kV for input/output lines</i>	<i>±2 kV for power supply lines</i>	<i>Mains power quality should be that of a typical commercial or hospital environment.</i>
<i>Surge IEC 61000-4-5</i>	<i>±1 kV line(s) to line(s) ±2 kV line(s) to earth</i>	<i>±1 kV line(s) to line(s)</i>	<i>Mains power quality should be that of a typical commercial or hospital environment.</i>
<i>Voltage dips, short interruptions and voltage variations on power supply i nput lines IEC 61000-4-11</i>	<i><5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 s</i>	<i><5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 s</i>	<i>Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains inter- ruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.</i>
<i>Power frequency (50/60Hz) magnetic field IEC 61000- 4-8</i>	<i>3A/m</i>	<i>3A/m</i>	<i>Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.</i>

NOTE UT is the a.c. mains voltage prior to application of the test level.

Table 4 Guidance and MANUFACTURER's declaration – electromagnetic IMMUNITY for ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING

Guidance and manufacturer's declaration – electromagnetic immunity

*The device is intended for use in the electromagnetic environment specified below.
The customer or the user of the device should assure that it is used in such an environment.*

IMMUNITY test	IEC 60601 TEST LEVEL	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	3 Vrms 3 V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance</p> $d = 1.167\sqrt{P}$ <p>$d = 1.167\sqrt{P}$ 80 MHz to 800 MHz</p> <p>$d = 2.333\sqrt{P}$ 800 MHz to 2,5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range b. Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

Table 6 Recommended separation distances between portable and mobile RF communications equipment and the ME EQUIPMENT or ME SYSTEM – for ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and the device.

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 1.167\sqrt{P}$	80 MHz to 800 MHz $d = 1.167\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.333\sqrt{P}$
0.01	0.117	0.117	0.233
0.1	0.369	0.369	0.738
1	1.167	1.167	2.333
10	3.690	3.690	7.378
100	11.67	11.67	23.33

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80MHz and 800MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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