LD7

Little Doctor®

Digital Blood Pressure Monitor

Instruction Manual

ENG

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
POL

TABLE OF CONTENTS

ENG

PARTS AND COMPONENTS	3
GENERAL INFORMATION	4
RECOMMENDATIONS ON CORRECT MEASUREMENTS	5
POWER SUPPLY	6
CORRECT POSITION DURING MEASUREMENT	8
SETTING OF DATE AND TIME, TEMPERATURE	8
CUFF PREPARATION	9
MEASUREMENT PROCEDURE	11
MEMORY FUNCTION	12
INFORMATION ABOUT ERRORS	13
CARE, STORAGE, REPAIR AND DISPOSAL	14
TROUBLESHOOTING TIPS	15
WARRANTY	16
TECHNICAL SPECIFICATIONS	16
CERTIFICATION AND STATE REGISTRATION	18
INFORMATION IN POLISH	20

1. Electronic Unit
2. LCD
3. Support
4. Arm cuff jack
5. Cuff CUFF-LDA
6. Button PRESET
7. MEM1 (memory 1)
8. MEM2 (memory 2).
9. START

10. Electrical power supply LD-N057
11. Electrical power supply jack
12. CR2020 battery (backup battery)
13. Batteries (main batteries)
14. Instruction Manual
15. Warranty card

This Instruction Manual is designed to assist the user with safe and effective operation of the automatic digital Device for measurement of blood pressure and heartbeat rate LD, (modification LD7) (hereinafter – the “Device”). Use this Device according to the rules described in this Manual. Operate the Device only as intended. Do not use the Device for any other purposes. Read and understand the whole Instruction Manual, in particular “Recommendations on Correct Measurement”.

INDICATIONS FOR USE

Use this Device to measure your systolic and diastolic blood pressure and heartbeat rate in patients aged from 15. This Device is recommended for use by persons with unstable blood pressure or known arterial hypertension at home as an addition to medical surveillance.

OPERATION PRINCIPLE

This Device uses the oscillometric method of blood pressure and pulse rate measurement. Wrap the cuff around your upper arm and it starts to be inflated automatically. The sensitive element of the Device feels the weak pressure oscillations in the cuff generated by widening and contraction of the brachial artery in response to every heartbeat. Pumping in is ceased when cuff is adequately pumped in to determine diastolic and systolic pressure (the amplitude of the pressure waves is measured, converted into millimeters Hg and shown on the display as figures). This device has 2 memories 60 cells each used for storage of measurement results and a clock with calendar and alarm features. An individual digital thermometer is integrated in the device for measurement of room temperature. Remember that the Device will not maintain the mentioned accuracy of a measurement if it is used or stored at a temperature or humidity other than those specified in Technical Specifications of this Manual. We are warning about possibility of mistakes in blood pressure measurement with this Device in persons with pronounced cardiac arrhythmia. Consult the doctor concerning blood pressure measurement of your child.

APPLIED NEW LD TECHNOLOGIES

Fuzzy Algorithm – algorithm for processing the measurement values with regard to peculiarities of the man’s heartbeat, thus, ensuring high measurement accuracy.

Scale WHO – classification of measurement results according to recommendation of World Health Organization (WHO).


Indication of arrhythmia – special symbol «♥♥» on Device display informs about availability of irregular pulse; in this case.

WARNING! This Device may be used only with cuff:

- Cuff-LDA, size 25-36 cm (delivered in a set with the Device);

RECOMMENDATIONS ON CORRECT MEASUREMENTS

1. For correct measurement you should know that THE BLOOD PRESSURE IS SUBJECT TO SHARP VARIATIONS EVEN WITHIN THE SHORT TIME INTERVALS. The blood pressure depends on many factors. It is usually lower in summer and higher in winter. The blood pressure varies together with the atmospheric pressure, depends on physical loads, emotional excitement, stresses and dietary regime. Drugs, drinking alcohol and smoking produce significant effect. Even the very procedure of blood pressure measurement in a polyclinic sends the blood pressure high in many people, thus, the blood pressure measured at home often differs from the values received in a polyclinic. As the blood pressure tends to rise at low temperatures, make measurements at an indoor temperature (approximately 20° C). If this Device stayed under a low temperature, keep it for at least 1 hour at an indoor temperature before use, otherwise the measurement result may be incorrect. During a day the difference in readings for healthy people may be 30-50 mmHg of systolic pressure and to 10 mmHg of diastolic pressure. The dependence of the blood pressure on various factors is individual for each person. Accordingly, it is recommended to keep a special book with blood pressure records. ONLY A CERTIFIED DOCTOR USING YOUR RECORDS IS CAPABLE TO ANALYZE THE TENDENCY OF YOUR BLOOD PRESSURE VARIATIONS.

2. At cardiovascular and some other diseases requiring blood pressure monitoring make measurements in the hours fixed by your attending doctor. REMEMBER THAT THE DIAGNOSTIC AND ANY TREATMENT OF HYPERTENSION MAY BE CONDUCTED ONLY BY A CERTIFIED DOCTOR ON THE BASIS OF BLOOD PRESSURE VALUES OBTAINED BY THIS DOCTOR. TAKING OF DRUGS AND THEIR DOSES SHOULD BE PRESCRIBED ONLY BY YOUR ATTENDING DOCTOR.

3. At such disorders as deep vascular sclerosis, weak pulse wave and also in patients with the prominent distortions of cardiac rhythm it may be difficult to measure the blood pressure accurately. IN SUCH CASES CONSULT A CERTIFIED DOCTOR ABOUT APPLICATION OF THE ELECTRONIC DEVICE.

4. KEEP QUIET DURING A MEASUREMENT TO OBTAIN THE ACCURATE VALUES OF YOUR BLOOD PRESSURE WITH THE ELECTRONIC DEVICE. Measure your blood pressure in the calm and comfortable conditions at the indoor temperature. No eating an hour before measurement; no smoking, taking tonic agents, alcohol 1.5-2 hours before measurement.
5. The accuracy of blood pressure measurement depends on whether the cuff matches the size of your arm. THE CUFF SHOULD NOT BE TOO SMALL OR TOO LARGE.
6. Wait 3 minutes between measurements for the blood to restore its circulation. However, the persons with prominent atherosclerosis due to considerable loss of vascular elasticity may need to increase the wait time between measurements (10-15 minutes). This also refers to the patients suffering for long from diabetes. For more accurate determination of blood pressure it is recommended to make a series of 3 consecutive measurements and to use the average value.

POWER SUPPLY OF THE DEVICE

This device is powered from two types of sources:

- 1. MAIN BATTERIES:** four AA batteries (1.5V) or ~220-240V supply network (hereinafter referred to as "power supply"). The main batteries support tonometer operation (cuff inflation, measurement and indication of results) with such functions as current time and date displayed.
- 2. BACKUP BATTERY:** one CR2020 (3V) battery that supplies power to an internal clock (time/date) when main batteries fail operating (e.g. in case of replacement of discharged batteries for new ones).

TONOMETER, CLOCK AND ALARM FUNCTIONS ARE ACTIVATED WHEN THE DEVICE IS POWERED FROM THE MAIN BATTERIES ONLY. WHEN THE DEVICE IS POWERED FROM A BACKUP BATTERY, CLOCK AND CALENDAR DATA ARE NOT DISPLAYED. BEFORE THE FIRST USE REMOVE A PLASTIC GASKET THAT IS FITTED BETWEEN CONTACTS AND BACKUP BATTERY.

For replacement of MAIN batteries:

1. Remove a support sliding it down off the electronic unit (Fig. 1). Open the battery compartment cover and insert 4 AA batteries as outlined in the drawing inside the compartment (Fig. 2).
2. Make sure that the battery polarities match those indicated inside the compartment. Do not apply excess force when removing a support and battery compartment cover. Close the battery compartment cover and fix the cuff support.

For replacement of BACKUP battery:

1. Remove a support and open the battery compartment cover. Remove the MAIN batteries, if they are installed. Open a round cover under the MAIN batteries. Replace a backup battery. Make sure that contact "+" faces up. Do not use much force to remove the cover of the battery compartment.
 2. Close the BACKUP battery compartment cover. Insert MAIN batteries, as shown above. Close the battery compartment cover and install the support.
- Replace MAIN batteries when the Low Battery Indicator "☐" appears on the screen or when there is no any indication on the screen. The Low Battery Indicator does not show the discharge level.
 - The BACKUP battery capacity is enough for powering the internal clock and calendar for approximately one year.
 - The MAIN batteries supplied with the Device are intended for check of the Device performance at sale and their service life may be shorter than of the recommended batteries.
 - Replace all four MAIN batteries at the same time. Do not use the waste batteries.
 - If the Device is unused for a long time, remove all batteries.
 - Do not leave the waste batteries in the Device.
 - Recharged AA batteries may be used.


Fig.1


Fig.2

USE OF THE DEVICE WITH THE POWER SOURCE

The manufacturer recommends application of the stabilized power source LD-N057 (it is attached to modification LD7). The jack for connection of the power source is on the right side of the Device.

To use the device with electric power source (EPS), connect plug connector of EPS to the Device and insert plug of EPS into mains socket, press «START».

Measurement over, switch off the Device having depressed button кнопки «START», remove plug of EPS from mains socket and disconnect plug connector from the Device. To avoid resetting of date and time, do not remove the batteries when using the device with power cells.

CORRECT POSITION DURING MEASUREMENT

1. Sit at a table so that during blood pressure measurement your hand rests on its surface. Be sure that the cuff is placed approximately at the level of your heart and that your arm lies freely on the table and does not move (fig.3).


Fig.3

2. You can measure the blood pressure lying on the back (fig. 4). Look at the ceiling, keep quiet and do not move during measurement. Be sure that the cuff is placed approximately at the level of your heart.


Fig.4

SETTING OF DATE, TIME AND TEMPERATURE

On replacing a backup battery, current date/time data will be reset. Enter correct date/time data to make it possible to save measurement results along with fixed date/time.

Attention: Note: Blood pressure and pulse may be measured without setting up date/time data.

This device may be used as a clock with an alarm, calendar and room thermometer functions (Fig. 5). Set up current date/time data and temperature measurement unit (i.e. Celsius or Fahrenheit):

1. For switching to setup mode press button "PRESET".
2. Selecting a year
Blinking year figures will be displayed. Press button MEM2 to increment the value of a year and button MEM1 to decrement the value of a year and set up a current year. Press button "START" to confirm the entry and go to the next step.
3. Selecting a month
Press button MEM2 to increment the value of a month and button MEM1 to decrement the value of a month and set up a current month.
Press button "START" to confirm the entry and go to the next step.

Fig.5

4. Section of a date

Press button MEM2 to increment the value of a date and button MEM1 to decrement the value of a date and set up a current date. Press button "START" to confirm the entry and go to the next step.

5. Selecting current time

A 24 hour format is used in the clock. Press buttons MEM1 and MEM2 to enter an hour value. Press button "START" to confirm the entry. Press button MEM2 to increment and button MEM1 to decrement the values of a minute and set up a current minute.

6. Press button "PRESET" to go to the alarm settings. Press button "START" to save the settings and to leave the preset mode.

7. Selecting an alarm

This device has three individual alarms. An alarm number (01, 02 or 03) is displayed near to the alarm symbol «» Press button "PRESET" to change the alarm value 01 to 02 and the value 02 to 03. Use buttons MEM1 and MEM2 also to adjust the values selected (i.e. HOURS-MINUTES). For adjustment of minute values press button "START" to change HOURS to MINUTES.

8. If the values of alarm actuation time are not set up, no figures will be displayed but the symbol of four dashes « - - : - - » will blink.

For switching off any preset alarm, select the required alarm and press button MEM1 holding it for 5 sec.

9. Setting the thermometer

On setting the alarm, press the preset button to select a thermometer measurement unit. Press button MEM1 to select a °C or °F value.

CUFF PREPARATION

1. Insert the cuff end for about 5 cm into a metal ring as shown in the figure (fig. 6).


Fig.6

2. Apply the cuff to your left upper arm so that the air tube is directed to your palm (fig. 7). If the measurement on your left arm is difficult, you may use your right arm. In this case remember that the readings may differ by 5-10 mmHg and even more.


Fig.7

3. Wrap the cuff around your upper arm so that the bottom of the cuff is approximately 2-3 cm above your elbow. The sign "ARTERY" should be over the arm artery (fig. 8).


Fig.8

4. Fix the cuff so that it fits tightly to the arm, but see that it is not overtight. Too tight or too free placement of the cuff may give inaccurate readings.


Fig.9

5. On the fixed cuff the sign «INDEX» should point to the area «NORMAL (25-36 cm)». It means that the cuff is chosen correctly and fits the size of your upper arm. If the sign points to the area marked «◀|||» the cuff is too small and the readings will be higher. If the sign points to the area marked «|||▶» the cuff is too large and the readings will be lower (fig. 10).


Fig.10

6. If the arm has a conic form, the cuff should be put on with a spiral movement as shown in the figure (fig. 11).


Fig.11

7. If the rolled-up sleeve squeezes the arm interfering with free blood flow the Device may give inaccurate figures not corresponding to your actual blood pressure (fig. 12).


Fig.12

MEASUREMENT PROCEDURE

1. Insert the Air Tube Plug into the Cuff Jack. Make 3-5 deep inhales and exhales before taking a measurement and relax. Do not move, do not speak and do not toughen your arm.

2. Press «START».

3. Press button MEM1 or MEM2 to set up the value of inflation pressure (189, 210, 240, 270, 300 mm Hg). The selected values are saved on default for further measurement. To finish setting press button "START". Should no button be pressed for the period of 3 sec., the device will automatically start measuring with reference to the current value of inflation pressure. Do not move, do not speak and do not toughen your arm during measurement.


Fig.13

4. All symbols will appear on the display screen for a short time, two short sound signals will be given and the Device will inflate automatically the air into the cuff.

5. On achieving a selected value, the cuff inflation pressure will gradually go down. The figures on the screen will count back. The pulse symbol « ♥ » will start flickering (fig. 14).

AS THE BLOOD PRESSURE AND PULSE ARE MEASURED DURING AIR DEFLATION FROM THE ARM CUFF KEEP QUIET AND DO NOT MOVE YOUR ARM AND DO NOT TOUGHEN YOUR ARM MUSCLES.

6. When the measurement is complete the sound signal is given, the arm cuff completely deflates and your measurement results flash on the screen (fig. 15).

Fig.14

If irregular pulse rhythm is detected during measurement, symbol of arrhythmia « ♥ » will appear upon measurement end. Appearance of arrhythmia indicator may also be caused by body movement during measurement. During periodical appearance of this indication apply to You attending doctor.

Apart from numerical value of pressure, result is also displayed on scale WHO (fig. 16). Scale WHO – three-color scale of classification of received value of arterial pressure, according to recommendation of World Health Organization. The scale is available from the left.

7. Symbols « 👤 » is blinking on display, reminding that to retain results, it is necessary to choose memory 1 or 2, having depressed MEM1 or MEM2. Press MEM1 (or MEM2) to retain results.

Fig.15

8. Press «START» to switch over to mode of date and time indication.

TO OBTAIN THE ACCURATE RESULT MAKE INTERVAL BETWEEN MEASUREMENTS TO RESTORE THE BLOOD CIRCULATION. WAIT FOR AT LEAST 3 MINUTES BEFORE MAKING A NEW MEASUREMENT.

If the Unit is ON and is unused for 3 minutes it will be switched off automatically.

AUTOMATIC RE-INFLATION

When during the first blood pressure measurement the cuff inflation to a level of 190 mmHg is not sufficient or you move your arm the Device stops measurement and re-inflates the cuff to the higher level. The Device has 4 fixed levels of the arm cuff inflation: 190, 230, 270 and 300 mmHg.

The automatic re-inflation is repeated until the measurement is completed successfully. This is not a defect.

FORCED DEFLATION FROM A CUFF

For rapid air release from of the arm cuff during arm cuff inflation or during a measurement (slow deflation) press the «START» button. The device will quickly release all air from cuff and will switch over to mode of date and time indication.


Fig.16

MEMORY FUNCTION

1. Result of each measurement (pressure, pulse, time and date of measurement) may be kept in the device memory. For this purpose, after measurement, within not more than 3 minutes, memory M1 or M2 shall be selected for memory storage MEM1 or MEM2.

IF THE NOTICE ON ERROR APPEARS THE MEASUREMENT RESULT WILL NOT BE STORED.

2. Up to 60 measurement results and mean value of last 3 measurements may be kept in each device memory. When the number of measurements exceeds 60, the oldest record is deleted to save the most recent values.
3. Press the MEM1 or MEM2 to see the figures stored in the memory. At the first depression of button MEM1 (or MEM2) mean value of 3 last measurements, kept in memory M1 (or M2), will appear on the screen, marked by symbol «A» (fig. 17). At repeated depression of button MEM1 (or MEM2) indicator of selected memory MEM1 «M1» (or MEM2 «M2») and number of memory cell will appear on the screen, and in 1 second its content is displayed (fig. 18).


Fig.17

When content of memory cell is displayed, date and time of measurement are displayed alternately in the display top line. Each depression of button MEM1 (or MEM2) causes shifting to the next memory cell.

If no measurement results are stored in the memory and MEM1 (or MEM2), buttons are pressed for the first time, current date and time will be displayed.


Fig.18

MEMORY CLEARING

To delete all values stored in the memory, press the MEM1 (or MEM2) Button and hold it down for more than 3 seconds. Symbols "Clr" will appear on the screen and all values will be deleted from the memory.

Measurement results (pressure, pulse, time and date of measurement) kept in the device memory even without power elements.

INFORMATION ABOUT ERRORS

<i>Indication</i>	<i>Likely cause</i>	<i>Methods of correction</i>
	<p>The arm cuff is applied incorrectly or the air tube plus is inserted not tightly enough.</p> <p>Measurements cannot be made due to hand movement or talking during measurements.</p>	<p>Be sure that the arm cuff is applied correctly and the plug is inserted tightly. Repeat the whole measurement procedure.</p> <p>Repeat the measurement following strictly the recommendations of this Manual.</p>
	Batteries are discharged.	Replace the batteries for new ones.

CARE, STORAGE, REPAIR AND DISPOSAL

1. Keep this Device from exposure to higher humidity, direct sunlight, shocks, vibration. **THIS DEVICE IS NOT WATERTIGHT!**
2. Do not keep and use this Device near heating installations and open fire.
3. If the Device was stored at a temperature below the freezing point, keep it at least for 1 hour in some warm place before use.
4. Remove the batteries if the Device will be unused for a long time. Battery leaking may damage the Device. **KEEP BATTERIES OUT OF REACH OF CHILDREN!**
5. Keep the Device clean and protect it from dust. Use the dry soft cloth to clean the Device.
6. Keep the Device and its components away from water, solvents, spirit and petrol.
7. Protect the arm cuff from contacting on sharp things; do not stretch or fold tightly the arm cuff.
8. Do not subject the Device to strong shocks, such as dropping on the floor.
9. This Device does not contain special controls to adjust the measurement accuracy. It is prohibited to open individually the electronic block. Repair the Device only in authorized organizations.
10. Do not leave the device unattended when plugged in the power supply line.
11. On expiration of the warranted service life apply from time to time to authorized repair organizations to check the technical condition of the Device.
12. Dispose of the Device and its components according to the application local regulations. No special requirements to disposal of this Device are defined by the manufacturer.
13. The arm cuff may withstand multiple sanitary treatments. The internal tissue surface of the arm cuff (contacting on arm) may be cleaned with cotton ball soaked in 3%-solution of hydrogen peroxide. At long use the partial color fading of the tissue coating of the arm cuff is possible. Washing and ironing of the arm cuff are not allowed.

TROUBLESHOOTING TIPS

PROBLEM	LIKELY CAUSE	METHOD OF CORRECTION
No indication on the display.	<p>Discharge of batteries.</p> <p>Polarity of batteries is not observed.</p> <p>Battery terminals are contaminated.</p> <p>Power source not plugged in an electrical outlet.</p>	<p>Replace all batteries for new ones.</p> <p>Install batteries correctly.</p> <p>Clean the terminals with dry cloth.</p> <p>Plug the power source into an electrical outlet.</p>
Inflation is stopped and resumed.	<p>Automatic re-inflation to obtain correct measurements.</p> <p>Perhaps you talk or move your arm during the measurement?</p>	<p>See MEASUREMENT PROCEDURE</p> <p>Calm down and repeat the measurement.</p>
Every time the blood pressure is different. Measurements are too low/high.	<p>Check that the arm cuff is level with your heart?</p> <p>Check that the arm cuff is applied correctly?</p> <p>Perhaps your arm muscles are tough?</p> <p>Perhaps you talk or move your arm during the measurement.</p>	<p>Take the correct position for measurement.</p> <p>Take the correct position for measurement.</p> <p>Calm down, apply the arm cuff correctly.</p> <p>Keep silence and quiet during measurement.</p>
Impossibility to make a large of number of measurements.	Application of poor batteries.	Use only alkali batteries of well-known manufacturers.
Measurements of the pulse rate are too high/low.	<p>Perhaps you talk or move your arm during the measurement.</p> <p>Perhaps you make measurement directly after physical load?</p>	<p>Keep silence and quiet during measurement.</p> <p>Repeat the measurement at least in 5 min.</p>

Independent switching to date, time and thermometer setting.	The device is automatically switched to a date/time display mode.	This is not a failure. The device is automatically switched to a date/time display mode in 3 minutes after using.
--	---	---

If regardless of the above recommendation you are unable to get correct measurement results, stop the use of this Device and apply to a maintenance organization (addresses and telephones of authorized organizations may be found in the warranty card). Do not try to adjust the internal mechanism by yourself.

WARRANTY

1. The following LD product is covered by warranty for the period specified in the warranty card.
2. The warranty liabilities are contained in the warranty card given at the sale of this unit to a purchaser.
3. The addresses of organizations for warranty maintenance are given in the warranty card.

TECHNICAL SPECIFICATIONS

Measurement method	oscillometric with Fuzzy Algorithm
Display	LCD, 4-lines
Pressure indication range in an arm cuff, mmHg	from 0 to 300
Measurement range: pressure in an arm cuff, mmHg pulse rate, 1/min	from 40 to 260 from 40 to 160
Temperature display range, °C	from 10 to 40
Range of admissible absolute error at measurement of air pressure in an arm cuff, mmHg	±3
Range of admissible absolute error at measurement of pulse rate measurement, %	±5
Inflation	automatically (air pump)

Pressure release	automatically
Date and time	yes
Memory	2x60 recent measurements + average value of the last three measurements
Power source, V	6
Type of power supply	4 AA (LR6) size batteries or power source or adapter no less than 600 mA
Max power intake, W	3,6

ADAPTER LD-N057

Output voltage, V	$6 \pm 5\%$
Max load current, mA	no more than 600
Input voltage, V/Hz	200-240/50
Dimensions, mm	80 x 79 x 50
Weight, g	293
Length of connecting wire, mm	1800 ± 100
Plug:	
Polarity of terminals	«-» internal
Internal diameter, mm	2.1 ± 0.1
External diameter, mm	5.5 ± 0.1
Length of plug contact, mm	10 ± 0.3
Operation conditions:	
Temperature, °C	from 10 to 40
Relative humidity, % Rh	85 and lower
Storage and transportation conditions:	
Temperature, °C	from -20 to 50
Relative humidity, % Rh	85 and lower
Cuff size:	larger for adults (upper arm circumference 25-36 cm)
Dimensions:	
Size (electronic block), mm	120 x 68 x 166
Weight (without package, case, batteries and adapter), g	611
Completeness	Electronic block, Support, Cuff Cuff-LDA (with air tube and air plug), Electrical power supply LD-N057, 4 Power elements, Power element CR2020, Instruction Manual, Warranty card, Package

Service life of the unit (without cuff), years	7
Service life of Electrical power supply, years	1
Service life of cuff, years	3
Year of manufacture	Year the manufacture is given in the bottom of the Device body in a serial number after symbols "AA"

Technical characteristics may be changed without preliminary notification to improve the operation and quality of the product.

CERTIFICATION AND STATE REGISTRATION

This unit manufacturing is certified according to international standard ISO 13485:2003. Units LD7 comply with the requirements of European Directive MDD 93/42/EEC, international standards, EN980, EN1041, EN1060-1, EN1060-3, EN10601-1-2, ISO 14971, EMC (IEC 60601-1-2:2001/A1:2004, CISPR 11:2003/A2:2006, GROUP 1, CLASS B, IEC 61000-3-2:2005, IEC 61000-3-3:1994/A2:2005).

Power source LD-N057 complies with international standard EN 55022 Class A, protection level against electric shock: Class II, Type B.

✉ **Complaints and requests should be addressed to:**

Little Doctor Europe Sp. z o.o.

57G Zawila Street, 30-390, Krakow, Poland

Service phone: +48 12 2684748, 2684749

Manufactured under control:

Little Doctor International (S) Pte. Ltd., 35 Selegie Road # 09-02 Parkline Shopping Centre, Singapore 188307

Manufacturer:

Little Doctor Electronic (Nantong) Co. Ltd., No.8, Tongxing Road Economic & Technical Development Area, Nantong 226010, Jiangsu, PEOPLE'S REPUBLIC OF CHINA

Distributor in Europe:

Little Doctor Europe Sp. z o.o.

57G Zawila Street, 30-390, Kraków, Poland

Sales Office phone: +48 12 2684746, 12 2684747, fax: +48 12 268 47 53

E-mail: biuro@littledoctor.pl

www.LittleDoctor.pl

Authorized Representative in the EU:
Shanghai International Trading Corp.GmbH
Eiffestrasse 80, 20537 Hamburg, Germany.

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ul. Zawila 57G, 30-390, Kraków, Polska

Serwis tel.: +48 12 2684748, 2684749.

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Little Doctor International (S) Pte. Ltd., 35 Selegie Road #09-02 Parkline Shopping Centre, Singapore 188307

Producent:

Little Doctor Electronic (Nantong) Co. Ltd., No.8, Tongxing Road Economic & Technical Development Area, 226010 Nantong, Jiangsu, PEOPLE'S REPUBLIC OF CHINA

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Little Doctor Europe Sp. z o.o., ul. Zawila 57G, 30-390 Kraków Polska

Biuro handlowe

tel.: +48 12 2684746, 12 2684747, fax: +48 12 268 47 53.

E-mail: biuro@littledoctor.pl

www.LittleDoctor.pl

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
Shanghai International Trading Corp.GmbH

Eiffestrasse 80, 20537 Hamburg, Germany.

www.LittleDoctor.sg

CE 0123

LITTLE DOCTOR INTERNATIONAL (S) PTE. LTD.

Yishun Central P.O. Box 9293 Singapore 917699,

Fax: 65-62342197, E-mail: ld@singaporemail.com

Shanghai International Trading Corp.GmbH
Eiffestrasse 80, 20537 Hamburg, Germany. Tel: 0049-40-2513175

® Registered trade marks of Little Doctor International (S) Pte. Ltd.

© Little Doctor International (S) Pte. Ltd., 2014

E349/1405/1