



Wrist Watch Blood Pressure Monitor



Characteristic:

W.H.O.
Indication

High
Normal High
Normal

iPD



Colored
Backlight



Multi-
-User



Specifications

KP-7030

Measuring Range	Pressure: 20~300 mmHg Pulse Rate: 40~200 pulse/min
Accuracy	Pressure: ± 3 mmHg Pulse Rate: ± 5 % of reading
Measuring Method	Oscillometric method
Inflation Method	Electrical rolling pump
Deflation Method	Mechanical release valve
Rapid Exhaust	Electrical solenoid valve
Display	Digital liquid crystal display
Memory	2x50, 3x33, 4x25, 5x20 sets
Operation Condition	5 ~ 40°C (41 ~ 104°F), 15 ~ 90 % RH
Storage Condition	-20 ~ 55°C (-4 ~ 131°F), < 93 % RH
Power Source	AAA alkaline battery x 2
Battery Life	Around 250 measurements with alkaline batteries
Power Saving	Auto-off after 150 sec non-operation
Size	75 x 75 x 32 mm
Weight	Approx. 134 g including batteries
Cuff Range	13.5 ~ 19.5 cm or 13.5 ~ 22 cm

Specifications are subject to change for improvement without notice.

K-jump Health Co., Ltd.

No.56, Wu Kung 5th Rd., New Taipei Industrial Park,
New Taipei City 24890, Taiwan


Tel: 886-2-22991378 Fax: 886-2-22991386

E-mail: service@kjump.com

<http://www.kjump.com>



CE 0197



HealthCare at Home
We Care for Your Health