

# K-jump®

## Multi-user Wrist Watch Blood Pressure Monitor



KP-7050

### Specifications

|                        |                                                                                 |
|------------------------|---------------------------------------------------------------------------------|
| Model                  | KP-7050                                                                         |
| Measuring Range        | 20 ~ 300 mmHg<br>40 ~ 200 pulse/minute                                          |
| Resolution             | 1 mmHg                                                                          |
| Accuracy               | ±3 mmHg or 2% (pressure)<br>±5% of readings (pulse rate)                        |
| Measuring Method       | Oscillometric method                                                            |
| Operation Key          | 4Keys: Power/Start, M, U, realtime IPD                                          |
| Display                | Liquid crystal display                                                          |
| Inflation              | Micro rolling pump                                                              |
| Exhaust                | Electrical Solenoid Valve                                                       |
| Memory                 | 2x50, 3x33, 4x25, 5x20, Max: 9 x 255 sets                                       |
| Time                   | Displaying date, time and year                                                  |
| Operation Condition    | Temperature: 10°C ~ 40°C (50°F ~ 104°F)<br>Relative Humidity: less than 85%     |
| Storage Condition      | Temperature: -10°C ~ +50°C (+14°F ~ +122°F)<br>Relative Humidity: less than 85% |
| Power Source           | 3Vdc; "AAA" type battery x 2                                                    |
| Battery Life           | Around 250 measurements(alkaline batteries)                                     |
| Low Battery Indication | Replace the batteries if the low battery indicator appears                      |
| Auto Off               | Automatically after 150 seconds without operation                               |
| Size                   | 83 x 69 mm x 34 mm                                                              |
| Weight                 | Approx. 142g (including batteries)                                              |
| Standard               | Complies with EN 1060-1, -3 and ANSI/AAMI SP10:2002                             |

Specifications are subject to change for improvement without notice.

### K-jump Health Co., Ltd.

No.56, Wu Kung 5th Rd., New Taipei Industrial Park,

New Taipei City 24890, Taiwan


Tel: 886-2-22991378 Fax: 886-2-22991386

E-mail: service@kjump.com

<http://www.kjump.com>



CE 0197



HealthCare at Home  
We care for Your Health