

Automatic Upper Arm Blood Pressure Monitor

KP- 7590

Folding stand



Specifications

Measuring Range	Pressure: 20~300 mmHg Pulse Rate: 40~200 pulse/min
Accuracy	Pressure: ± 3 mmHg Pulse Rate: ± 5 % of reading
Measuring Method	Oscillometric method
Inflation Method	Electrical rolling pump
Deflation Method	Mechanical release valve
Rapid Exhaust	Electrical solenoid valve
Display	Digital liquid crystal display
Memory	2x20
Operation Condition	5 ~ 40°C (41 ~ 104°F), 15 ~ 90 % RH
Storage Condition	-20 ~ 55°C (-4 ~ 131°F), < 93 % RH
IP Classification	IP20 (Protection against touching foreign bodies <12.5mm and no water-proof test)
Power Source	AA alkaline battery x 4 or AC adaptor with DC 6V/0.8A output(Optional)
Battery Life	Around 250 measurements with alkaline batteries
Power Saving	Auto-off after 150 sec non-operation
Size	110.6 x 187.3 x 56.3 mm
Weight	Approx. 400 g including batteries
Cuff Range	20 ~ 35 cm , 22 ~ 32 cm , 32 ~ 42 cm

Specifications are subject to change for improvement without notice.

K-jump Health Co., Ltd.

No.56,Wu Kung 5th Rd.,New Taipei Industrial Park,

New Taipei City 24890,Taiwan

Tel: 886-2-22991378

Fax:886-2-22991386

E-mail: service@kjump.com

http://www.kjump.com



CE 0197

Characteristic:

W.H.O.
Indication



iPD



Colored
Backlight



Multi-
User



Jumbo
DISPLAY



Low Battery
Indication



Optional:

Multi OS
Supported

Android
iOS

MOMI



Data
Measurement



Graphic
Analysis



Near Field
Communication

NFC

BlueTooth



1. Read memory back to portable device
2. Write in Settings to monitor
3. Real time measurement

The Bluetooth® word mark and logo are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by K-JUMP HEALTH CO., LTD. is under license. Other trademarks and trade names are those of their respective owners.

HealthCare at Home

We Care for Your Health