

K-jump[®]

Automatic Upper Arm Blood Pressure Monitor



KP-7720

Specifications

Measuring Range	Pressure: 20~300 mmHg Pulse Rate: 40~200 pulse/min
Accuracy	Pressure: ± 3 mmHg Pulse Rate: ± 5 % of reading
Measuring Method	Oscillometric method
Inflation Method	Electrical rolling pump
Deflation Method	Mechanical release valve
Rapid Exhaust	Electrical solenoid valve
Display	Digital liquid crystal display and portable device
Memory	Portable Device: Based on memory size Blood pressure monitor: 42, 85, 98 sets
Operation Condition	5 ~ 40 °C(41 ~ 104°F), 15 ~ 90 % RH
Storage Condition	-20 ~ 55 °C(-4 ~ 131°F) < 93 % RH
Power Source	AA alkaline battery x 4
Battery Life	Around 250 measurements with alkaline batteries
Power Saving	Auto-off after 150 sec non-operation
Size	120 x 120 x 75 mm
Weight	Approx. 335 g including batteries
Cuff Range	22 ~ 32 cm or 32 ~ 42 cm

Specifications are subject to change for improvement without notice.

K-jump Health Co., Ltd.

No.56, Wu Kung 5th Rd., New Taipei Industrial Park,
New Taipei City 24890, Taiwan

Tel: 886-2-22991378 Fax: 886-2-22991386

E-mail: service@kjump.com

<http://www.kjump.com>



CE 0197



Characteristic:

W.H.O.
Indication

High
Normal High
Normal

Portable
Device



Fuzzy
Logic



iPD



Graphic
Analysis



Multi-
-User



Data
Measurement



Talking



Optional:

Colored
Backlight

LCD

1. Read memory back to portable device
2. Write in Settings to monitor



HealthCare at Home

We Care for Your Health